
PULSE OXIMETRY

PRINCIPLES of OPERATION

PULSE OXIMETRY – SUBSYSTEM

The monitor is based on the principles of Spectrophotometry and Plethysmography. It includes an electro-optical sensor and a microprocessor-based monitor. The sensor has two low-voltage light-emitting diodes (LED's) as light sources, and one photodiode as a photo-detector. One LED emits red light (nominal 660nm) and the other emits infrared (nominal 920nm). When the light from the LED's passes through the sensor site, part of it is absorbed. The photo-detector measures the light that passes through, which indicates red and infrared absorption.

With each heartbeat, a pulse of oxygenated arterial blood flows to the sensor site. Oxygenated haemoglobin differs from deoxygenated haemoglobin in its relative red and infrared absorption, and the monitor measures red and infrared absorption to determine the percentage of functional haemoglobin that is saturated with oxygen.

Light absorption that is measured when pulsatile blood is not present, reflects absorption by tissue and nonpulsatile blood absorption that does not change substantially during the pulse. This is analogous to the reference measurement of a Spectrophotometer. Absorption is also measured when pulsatile, arterial blood is in the tissue. The monitor then corrects this measurement for absorption when the pulsatile blood is not present. The ratio of the corrected absorption at each wavelength determines arterial oxygen saturation (SpO₂).

AUTOMATIC CALIBRATION

The oximetry subsystem incorporates automatic calibration mechanisms. It is automatically calibrated each time it is turned on, at periodic intervals thereafter, and when ever a new sensor is connected. Also the intensity of the sensors LED's is adjusted automatically to compensate for differences in tissue thickness.

Each sensor is calibrated when manufactured: the effective mean wavelength of the red LED is determined and encoded into a calibration resistor in the sensor plug. The instruments software reads this calibration resistor to determine the appropriate calibration coefficients for the measurements obtained by that sensor.

FUNCTIONAL versus FRACTIONAL SATURATION


Because the monitor measures functional SpO₂, it may produce measurements that differ from those of instruments that measure fractional SpO₂. Functional SpO₂ is oxygenated haemoglobin expressed as a percentage of the haemoglobin that is capable of transporting oxygen. Because the monitor uses two wavelengths, it measures oxygenated and deoxygenated haemoglobin, yielding functional SpO₂. It does not detect dysfunctional haemoglobin, such as carboxyhaemoglobin or methemoglobin.

In contrast, some laboratory instruments such as the Instrumentation Laboratory 282 Co-Oximeter, report fractional SpO2-oxygenated haemoglobin expressed as a percentage of all measured haemoglobin, whether or not that haemoglobin is available for oxygen transport. Measured dysfunctional haemoglobins are included. Consequently, to compare this monitor's measurements directly with those of another instrument, that other instrument must measure functional SpO2. If it measures fractional SpO2, those measurements can be converted using the following equation:

$$\text{Functional Saturation} = \text{Fractional Saturation} \times \frac{100}{100 - (\% \text{carboxyhaemoglobin} + \% \text{methemoglobin})}$$

MEASURED versus CALCULATED SATURATION

When SpO2 is calculated from a blood gas measurement of the partial pressure of arterial oxygen (PaO2), the calculated value may differ from the monitors SpO2 measurement. This is because the calculated SpO2 may not have been corrected for the effects of variables that shift the relationship between PaO2 and SpO2 (see figure below): temperature, pH, the partial pressure of carbon dioxide (PaCO2), and the concentrations of 2,3-DPG and foetal haemoglobin.



Oxyhaemoglobin Dissociation Curve

SAFETY INFORMATION**WARNINGS**

DANGER!! Explosion Hazard. Do not use in the presence of flammable anaesthetics.

Carefully read any operator manuals, accessory directions for use, all precautionary information, and specifications before application and use of the equipment.

Monitors are intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.

Do not silence any audible alarms if patient safety could be compromised.

Tissue damage can be caused by incorrect application, or use, of a sensor (e.g. wrapping the sensor too tightly, applying supplemental tape, or failing to inspect the sensor site periodically).

Loss of pulse signal can occur if:

- The sensor is too tight;
- There is excessive illumination: e.g. a surgical or bilirubin lamp or sunlight;
- The sensor is placed on an extremity with a blood pressure cuff, arterial catheter, or intravascular line;
- The patient is in shock, has hypotension, severe vasoconstriction or anaemia, hypothermia, arterial occlusion proximal to the sensor, or cardiac arrest.

Inaccurate measurements may be caused by:

- Incorrect application or use of a sensor;
- Significant levels of dysfunctional haemoglobins;
- Significant levels of indocyanine green, methylene blue or other such intravascular dyes;
- Exposure to excessive illumination, such as surgical lamps, bilirubin lamps, fluorescent lights, infrared heating lamps or direct sunlight;
- Excessive patient movement;
- Venous pulsations;
- Electro surgical interference;
- Placement of the sensor on an extremity that has a blood pressure cuff, arterial catheter or intravascular line.

Do not use oximeters / sensors during Magnetic Resonance Imaging (MRI) scanning. Conducted current could potentially cause burns. Also the pulse oximeters may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.

Do not use a damaged sensor, or one with exposed electrical contacts. Do not immerse sensors in liquids. Do not sterilise by irradiation, steam or ethylene oxide.