DL500FS: Foetal Heart Simulator evaluation: SW: 25-5-01.

<u>Performance.</u> DL500FS trialed against HP3010A & Sonicaid ultrasound monitors:

Hp3010A

- (\checkmark) able to pick up and display exact simulated heart rate.
- (\checkmark) relatively realistic heart beat sound.
- (x) max. amplitude setting on simulator swamps monitor.
- (x) oscilloscope display shows twin peaked trace for every heart beat.
- (x) oscilloscope display very clean not enough background noise.
- (x) heart beat sound possibly too clear.
- (x) 30 bpm setting outside performance of monitor.

Sonicaid

- (\checkmark) able to pick up and display exact simulated heart rate.
- (\checkmark) relatively realistic heart beat sound.
- (x) heart beat sound possibly too clear.
- (x) 30 bpm setting outside performance of monitor.

User appeal.

- (\checkmark) Charge / fast charge / low battery indicator.
- (x) Large amount of vibration / audible noise from unit particularly at high amplitude setting
- (?) Black case rather than white less prone to looking dirty.
- (?) Tactile facia panel incorporating on / off button rather than switch.
- (?) Amplitude pot on facia possibly cosmetically recessed slightly.
- (?) Too large ?
- (?) Substantial enough?

Engineering questions.

- (?) Battery duration?
- (?) Charge time?
- (?) Case to be gel proof?
- (?) Unit defaults to 120 bpm at switch on.
- (?) Battery to be sealed in or purposely removable.
- (?) EMC shielded case required.
- (?) Can case be modified to allow full circle of rubber pad for circular ultrasound transducer
- (?) plastic of case immediately under rubber pad requires more amplitude from internal transducer and hence more vibration and audible noise.
- (?) How easy to modify unit to accommodate next generation model giving abnormal conditions.
- (?) Microprocessor based device
- who will be allowed to alter software?

- documentation record of software alterations?

- software based on which language?
- software based on which language
- software disclosed to Viamed?
- (?) Unit draws 42 200mA from battery dependent on amplitude pot setting.