

Training Feedback Form

Training Course Completed: TOF 3D Introductory Training			
Date: 20/02/2023	Time/Length:	Trainer: Catrin Hollings	
Content	Yes	No	Unsure
Was the course content presented in a logical manner?			
Was the course content and material complete and comprehensive?			
Will this information be useful to you in your job role?			
Relevance	Yes	No	Unsure
Do you feel you now have a better understanding of the product/procedure/training area*?			
Did the course challenge your thinking and understanding of the product/procedure/training area*?			
Do you feel the training is beneficial to your team?			
Trainer	Yes	No	Unsure
Did the trainer communicate and explain the material clearly?			
Did you feel the instructor was knowledgeable in the area covered?			
Did the trainer encourage discussions and questions?			
Comments			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: <i>Aye Swaine</i> Date: <i>28/2/2023</i>			

*Please delete as applicable

TOF 3D Introductory Competency Questions
Please use your PCs as reference

Please use your PCs as reference

1. What standard accessories are included when purchasing 2510091?

- PATIENT CABLE ✓
- STIMULATION CABLE ✓
- ACCELERATION TRANSDUCER ✓
- SEALING COVER FOR PATIENT CABLE PORT ✓
- " ————— " Drain interface between ✓
- TPU ✓
- INTERFACES ✓

2. What is a brief description of the TOF 3D Neuromuscular Monitor?

THE T03D QUANTITATIVELY MONITORS THE LEVEL OF NEUROMUSCULAR TRANSMISSION DURING SURGERY OR IN ICU AS A MEANS OF ASSESSING THE DEPTH OF MUSCLE RELAXATION ✓

3. Is the TOF 3D a subjective or quantitative monitor and what does this mean?

Quantitative - This measurement is carried out by the machine, not the operator.

4. What is the main difference between the Microstim DB3 and the TOF 3D?

The microsim was a subroutines machine & has no display

5. What part numbers does 2510101 complete patient cable assembly consist of, including part number?

- 2510102 - MAIN CABLE
- 2510103 - STIM CABLE
- 2510104 - ACCELERATION TRANSDUCER
- SHALING CABLE & PATIENTS CABLE (NOT)

FAA' memo

6. What monitoring sites can the TOF 3D be used and what accessories are required to do this?

- HAND + THUMB using HAND ADAPTATION OR THUMB ADAPTATION ✓
- TOES - using THUMB ADAPTATION ✓
- FACE - using EYE ADAPTATION ✓

7. What training can Viamed provided for the TOF 3D?

no on-site training offered, but function training has been given on a website ✓