

Distributed by

Replacement SpO₂ Probe

This compatible replacement adult oxygen saturation and pulse rate monitoring SpO₂ probe, is for use when continuous non-invasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing greater than 40kg.

For use with product compatible oximeter (See table on reverse)

Instructions for Use

An index finger is the preferred location for the sensor.

Orientate the sensor so that the fingernail is directly under the sensor. Plug the sensor into the appropriate oximeter and verify proper operation in accordance with the original manufacturers oximeter manual.

NB: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff or intravascular infusion line.

Reusable sensors may be

used on the same site for a maximum of 4 hours, provided the site is inspected routinely to ensure skin integrity and correct positioning. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

NB: If the patient wears fingernail polish or artificial fingernails choose a different probe site, or remove the fingernail polish or artificial fingernails.

Troubleshooting

The probe may not reliably track pulse and SpO₂ if the finger is incorrectly positioned within the sensor, or the sensor site is too thick, thin, or deeply pigmented to permit appropriate light transmission. If any of these situations occur, reposition the sensor.

Avoid any site that may have excessive ambient light, excessive motion, or electrical interference

Not for use within MRI

scanning suites.

Cleaning

Sensor may be surface cleaned with 70% isopropyl alcohol.

Do not immerse the sensor or connector in any liquid solution.

Do not autoclave this probe.

The SpO₂ Finger probe will not be harmed by exposure to an ethylene oxide mixture at 49° to 54°C (120° to 130°F)

Note

Refer to the instrument operator's manual for complete instructions for use of the probe with the monitor.

Accuracy

Compatible to Original manufacturers specification

warranty is to repair or replace the product.

To return a product under warranty the local distributor should be contacted for a return authorisation.

Internal components are manufactured to our proprietary design specifications. Viamed does not have a licensing agreement with any Original Oximeter manufacturer

Manufactured By MCI USA.

WARRANTY

This product is warranted for a period of 12 months from the date of purchase against defects in workmanship and materials by Viamed whose sole obligation under this