Spirolog®

de	Flowsensor
en	Flow Sensor
fr	Capteur de débit
it	Sensore di flusso
es	Sensor de flujo
pt	Sensor de fluxo
sv	Flödessensor
no	Flowsensor
da	Flowsensor
fi	Flowsensor (Virtausanturin)
nl	Flowsensor
pl	Czujnik przepływu
el	Αισθητήρας ροής
ja	フローセンサー
enUS	Flow Sensor
cs	Průtokový senzor
hu	Áramlászérző
ru	Датчик потока
zh	流量传感器
It	Tékmés jutiklis
sk	Prietokový senzor
sr	Senzor za protok
tr	Akış Sensörü

de	Kurzanleitung
en	Short Instructions
fr	Instruction d'utilisation abrégée
it	Brevi istruzioni per l'uso
es	Instrucciones resumidas
pt	Instruções curtas
sv	Snabbinstruktion
no	Kort veiledning
da	Kort vejledning
fi	Pikaopas
nl	Handleiding
pl	Krótką instrukcja
el	Σύντομος οδηγίες χρήσης
ja	簡易取扱説明書
enUS	Quick Reference Guide (for use in the USA)
cs	Stručný návod
hu	Rövid útmutató
ru	Краткая инструкция
zh	简要说明
It	Trumpa instrukcija
sk	Stručné pokyny
sr	Kratko uputstvo
tr	Kısa Talimatlar

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Spirolog®*

en Flow Sensor Short Instructions

Note the Instructions for Use of the basic equipment!

Intended use

Spirolog® Flow Sensor (not sterilisable) – hot-wire sensor for measuring the volumetric flow (Flow) and calculating the volume in Dräger anaesthesia and ventilation units.

Operation

Do not introduce any medications or other substances based on combustible solvents, such as alcohol, into the patient system. Fire risk!

Do not use explosive anaesthetic agents (e.g. ether or cyclopropane)!

Residue of nebulised medications can shorten the service life of the flow sensor. Therefore

- After using with medicament nebulisers, immediately clean and treat the flow sensor as specified
- Alternatively
- Install a filter upstream of the flow sensor

Care

Do not clean in an automatic cleaning/disinfecting machine. Do not clean with compressed air, water jet, brush etc. The thin wires in the flow sensor would be destroyed.

Regularly check the flow sensor for damage – visual inspection

- Disinfected for approx. 1 h in 70 % ethanol or isopropanol solution
 - Then allow to dry in air for at least 30 minutes
- Residual alcohol can otherwise lead to destruction of the flow sensor during calibration.

The following instrument disinfectants may be used:

- Korsalex Basic, Korsalex Extra, Korsalex Plus (Bode Chemie, Hamburg, Germany)
- Gigasept FF, Gigasept AF (Schulke & Mayr, Norderstedt, Germany)

Sterilising

The Spirolog® Flow Sensor cannot be sterilised.

Disposal

Dispose of spent sensor with infectious special waste
Low-emission incineration at more than 800 °C.

Technical data

Material Spirolog® Flow Sensor	Acrylonitrile-butadiene-styrene (ABS) with 13 µm platinum wire
Environmental conditions Operation	See Instructions for Use of basic equipment
Storage Temperature	–40 to 75 °C
Atmospheric pressure	115 to 1100 hPa
Relative humidity	5 to 95 % (no condensation)
Flow measurement Measuring principle	Hot-wire anemometry
Service life Spirolog® Flow Sensor	Can be used as long as successful calibration is possible

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© Dräger Medical AG & Co. KGaA
* Spirolog® is a registered trade mark of Dräger
Dräger Medical AG & Co. KGaA