OxiMax[®] N-65[™]

English	2
Français	3
Deutsch	4
Nederlands	5
Italiano	6
Español	7
Svenska	8
Dansk	9
Norsk	10
Suomi	11
Português	12
Русский	13
中文	14
Polski	15
Cesky	16
Magyar	17
Greek	18
TÜRKÇE	19
العربية	20
日本語	21

Oxygen Saturation Accuracy Specification Grid for OxiMAX® N-65™ Pulse Oximeter

OxiMAX® Sensor Models	SpO ₂ Range
Accuracy Specifications	
Accuracy specifications are based on controlled hypoxia studies with healthy non-smoking adult volunteers over the specified saturation SpO ₂ range(s). ² Pulse oximeter SpO ₂ readings were compared to SaO ₂ values of drawn blood samples measured by hemoximetry. All accuracies are expressed as $\pm$ "X" digits. Pulse oximeter equipment measurements are statistically distributed; about two-thirds of pulse oximeter measurements can be expected to fall in this accuracy (A _{RMS}) range.	
Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO ₂ , SaO ₂ , or both. Accordingly, observations of clinical accuracy may not achieve the same levels as those obtained under controlled laboratory conditions.	
Neonatal Accuracy	
When sensors are used on neonatal subjects as recommended, the specified accuracy range is increased by ± 1 digit , as compared to adult usage, to account for the theoretical effect on oximeter measurements of fetal hemoglobin in neonatal blood. For example, OxiCliq N accuracy on neonates is ± 3.5 digits, rather than ± 2.5 .	
Oxygen saturation accuracy can be affected by certain environmental, equipment, and patient physiologic conditions (as discussed in the operator's manual for the monitor) that influence readings of SpO ₂ , SaO ₂ , or both. Accordingly, observations of clinical accuracy may not achieve the same levels as those obtained under controlled laboratory conditions.	
Compatibility Information	
Use Nellcor sensors only with Nellcor patient monitors and monitors containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for accuracy specifications and compatibility information of particular instruments and Nellcor sensor models.	
Range of Peak Wavelength and Maximum Output Power	
Nellcor pulse oximeter sensors contain LEDs that emit red light at a wavelength of approximately 660 nm and infrared light at a wavelength of approximately 900 nm. The total optical output power of the sensor LEDs is less than 15 mw.	
OxiMAX Single Patient Use	70%–100%
MAX-A, MAX-AL	± 2
MAX-N (Adult)	± 2
MAX-P	± 2
MAX-I	± 2
MAX-FAST	± 2
SC-A (Adult)	± 2
SC-PR (Neonate)	± 3
SC-NEO (Neonate)	± 3
MAX-R ¹	± 3.5
OxiMAX OxiCliq® Single Patient Use	70%–100%
OxiCliq A	± 2.5
OxiCliq P	± 2.5
OxiCliq N (Adult)	± 2.5
OxiCliq N (Neonate)	± 3.5
OxiCliq I	± 2.5
OxiMAX Reusable Sensors	70%–100%
D-YS (Infant to Adult)	± 3
D-YS (Neonate)	± 4
D-YS with D-YSE Ear Clip	± 3.5
D-YS with D-YSPD SpotClip	± 3.5
DS-100A	± 3
OXI-A/N (Adult)	± 3
OXI-A/N (Neonate)	± 4
OXI-P/I	± 3

1. The accuracy specification has been determined between saturations of 80%–100%.

2. Subjects used to validate SpO₂ measurement accuracies were healthy and recruited from local population. Comprised of both men and women, subjects spanned a range of skin pigmentations and ranged in age from 18-50 years old.

Sensor Accuracy Grids can be viewed online at www.nellcor.com

© Copyright 2006 Nellcor Puritan Bennett Inc. All rights reserved.

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© Copyright 2006 Nellcor Puritan Bennett Inc. Tous droits réservés.
© Copyright 2006 Nellcor Puritan Bennett Inc. Alle Rechte vorbehalten.
© Copyright 2006 Nellcor Puritan Bennett Inc. Alle rechten voorbehouden.
© Copyright 2006 Nellcor Puritan Bennett Inc. Tutti i diritti riservati.
© Copyright 2006 Nellcor Puritan Bennett Inc. Todos los derechos reservados.
© Copyright 2006 Nellcor Puritan Bennett Inc. Alla rättigheter förbehålles.