

BIS and us

Pierre Kaygin
Sales and Clinical Applications

*Certified Anesthesia Nurse
MBA Int
Mail: pierre.kaygin@anandic.com*

For internal use only



History of Anandic Medical Systems AG

Our Domestic Business

- Distributor of GE over last >26 years
- 60 Employers
- Good after sales support for devices
- Selling BIS sensors 1>Mio CHF
- High quality demand in market

Export Business



- BIS sensors for Medtronic and OEM modules
- Entropy sensors for GE
- Made in China
- Warehouse in Europe, fast delivery
- And after sales support



For internal use only

Waves of EEG

Electric signals are measured and displayed



For internal use only

Spectral Entropy

Level of consciousness

Entropy parameters



RESPONSE ENTROPY (RE):

- Sensitive to activation of facial muscles
- Rapid reaction

STATE ENTROPY (SE):

- Assesses the hypnotic effect of anaesthetics on the brain
- Stable and robust

Entropy monitoring is based on the acquisition and processing of raw electroencephalogram (EEG) and frontal electromyogram (f-EMG) data.




Spectral Entropy

Level of consciousness

100	Fully awake and responsive
60	Clinically meaningful anaesthesia with low probability of consciousness
40	
0	Suppression of cortical electrical activity

Note: the data shown in the graph are a simulation of a typical case.



Spectral Entropy

Single-use patient sensors (25/box)



M1174414
EasyFit Sensor




M1038681
GE Entropy Sensor

Reusable GE entropy cable



M1050784



Entropy vs BIS

Comparison:

- **Two meaningful indexes** - State Entropy (SE) based on **1-channel raw EEG** and Response Entropy (RE) based on **EMG**. For BIS EMG is just an artifact
- **GE uses other parameters like NMT and SSI(Surgical Stress Index)**
- **More real time information** – RE reacts in few seconds, while BIS may take 30 to 120 seconds to reflect the state change⁽¹⁾ → RE detects status changes **quicker than BIS**




RE/SE Display



E-Entropy vs. E-BIS




Entropy & EasyFit



Other technologies

Masimo Sedline, USA

- **SedLine** Patient State Index (PSi) for monitoring state of the brain during anesthesia - bilateral acquisition and processing of 4-ch EEG signals. RD SedLine EEG sensor and O3 Regional Oximetry sensors. Monitor Root
- 4-ch EEG more diagnostic and improved performance during interference (electrocautery, EMG interference, low power EEG)
- Raw EEG analyzes




Medtronic BIS, USA

- BIS for depth of anesthesia in OR
- BISx4: OR - > ICU + 4-ch EEG
 - OR; prognostic value esp. cardiac surgery
 - ICU; avoid awareness/over-light sedation, detect changes and help treat status epilepticus
- Standalone monitor available
- Module integration by all major competitors (GE, Mindray, Philips etc.)
- Uses EEG waveforms in database and compares this
- Sensor can be 3 times connected, after it cannot be used




Our commitment

- Best compatible BIS and Entropy sensors
- High quality
- Affordable
- Experienced company
- Our BIS can be connected 10 times to a Medtronic device/Module
- Soft sensors, less irritation



PEOPLE WHO HEALTHCARE.

by anandic




Disposable EEG Sensors for Measurement of the Bispectral Index

Disposable EEG sensors for monitoring the depth of anesthesia with the Bispectral Index, as well as acquisition and transmission of biological EEG signals of anesthesia depth of patients.

Features:

- Compatible with all available BIS monitors and modules
- Clinically tested high measurement accuracy
- Soft and biocompatible material
- Latex free
- Shelf life from production: 1,5 year

Adult	MK-01
Children from 1 year	MK-02
Bilateral Adults	MK-04
Packaging unit	25 pieces



ANANDIC MEDICAL SYSTEMS AG
Postfach, Stadtweg 24
8245 Feuerthalen

info@anandic.com
www.anandic.com
0848 800 950

