SpiroTrue® A

Flow sensor

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en Instructions for use

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Intended use and medical purpose of SpiroTrue®A

Intended use:

SpiroTrue®A is a hot wire anemometric flow sensor, for measuring gas flow in Dräger anesthesia and ventilation equipment. SpiroTrue®A flow sensor is compatible to:

- Dräger Spirolog, REF 8403735-06

For further information on system compatibilities, take into consideration the lists of accessories or other documents provided by the manufacturer of the basic device.

Intended medical indications:

Flow sensor for measuring respiratory gas flow in Dräger anesthesia and ventilation equipment (**without RFID** functionality).

Intended user profile:

Healthcare professional, trained for the application (ventilator).

Intended patient groups:

Adult and pediatric, limited by the application (ventilator).

Clinical environment:

Professional healthcare facility.

Contraindications and other restrictions:

Not intended to be used with flammable anesthetic mixtures.

Not intended to be autoclaved or sterilized.

Manufacturer information

bluepoint MEDICAL GmbH & Co.KG
An der Trave 15
23923 Selmsdorf, Germany

Phone: +49 (38823) 5488 – 0
Fax: +49 (38823) 5488 – 29
E.Mail: info@bluepoint-medical.com
Web: www.bluepoint-medical.com

CE
0482

General warnings and safety notices

Safety symbols used in this document:

	Indicates potential harmful conditions that may lead to serious injury.
	Important information on a topic or procedure, or conditions that may lead to damage or malfunctions of the device.


General warnings:

	Do not use this medical product if it appears or is suspected to be damaged. Damage to the product may result in degraded performance.
	Do not use this product if its packaging has previously been opened. Risk of infection, cross contamination!
	No modification of this medical product is allowed.
	Strictly follow the instructions for use. Handling of the device requires full understanding and strict observation of the operating instructions of basic device (ventilator) and these instructions for use.
	The flow sensor must be calibrated in the basic device (ventilator) and system check must be performed prior to patient use. Observe the selected gas/gas-mixture type.
	Due to fire risk - do not introduce any medicaments or other substances based on combustible solvents, such as alcohol, into the patient system. Do not use explosive anesthetic agents e.g. ether or cyclopropane.

General safety notices:

	The medical device is intended only for the application described herein.
	Any serious incident that has occurred in relation to this medical device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. To report patient safety related events, send an email to the following manufacturer's address: prrc@bluepoint-medical.com with at least following information: - order number and product name - batch number of the product - date of the serious incident - description of the serious incident affecting the patient or any injuries - your contact data (institution, address, contact person/representation, title and telephone number)

SpiroTrue®A Sensor overview

1. Product name
2. Sensor tube
3. Sensor connector
4. Sensing wires
5. Sensor pins

Installing, calibrating and operating of the flow sensor

	This medical product must be connected/installed to the basic device (ventilator) according to the instructions for use of basic device. Risk for faulty operation!
	Ensure secure connection of this product to the basic device (electrically, mechanically). Risk of leakages and faulty operation!
	The flow sensor has to be securely connected to the flow sensor cable /flow sensor socket of the basic device (ventilator). The sensor cable must be loosely aligned so as to not stress the flow sensor plug during operation. Risk of faulty operation!
	The flow sensor must be calibrated in the basic device (ventilator) prior to patient use.
	The flow sensor calibration must be performed in accordance with the instructions for use of the basic device (ventilator). If a malfunction is detected, do not operate the device! Patient hazard!
	A preoperational device system check must be performed according to the instructions for use of the basic device (ventilator) prior to patient use. If a malfunction is detected, do not operate the device! Patient hazard!
	Due to fire risk – do not use this product with medicaments or other substances based on combustible solvents, such as alcohol, into the patient systems.
	Do not use explosive anesthetic agents (e.g. ether or cyclopropane).
	Life cycle of the flow sensor is limited. It can be used depending on hospital's infection control regulations and as long as calibration in the basic device (ventilator) is possible.

Cleaning and disinfecting of the flow sensor

Classification of SpiroTrue®A flow sensor for reprocessing according to the Spaulding classification system:

Semi-critical

	Reusing the flow sensor may affect patient safety due to the possible contamination of the sensor.
	Reusable products must be cleaned and disinfected after every use.
	Use only validated procedures for reprocessing of this medical product and refer to the hospital's infection control regulations.
	Use only certified cleaning agents and disinfectants with good compatibility to the materials of medical device.
	Strictly follow the manufacturer's instructions on cleaning agents and disinfectants.
	Ensure particle-free cleaning and disinfecting procedure.
	Do not clean and disinfect this medical product in automated cleaning or disinfecting equipment!
	Do not clean with compressed air, water jets, brushes etc.
	Do not use ultrasound cleaning to clean the flow sensor.
	Do not use a thermal disinfection to disinfect the flow sensor.

Validated reprocessing procedure for SpiroTrue®A:


Manual cleaning and disinfection procedure of SpiroTrue®A by immersion at room temperature:

Validated cleaning agents and disinfectants for SpiroTrue®A flow sensor			
Cleaning agent / Disinfectant	Manufacturer	Concentration	Contact time
Neodisher MediClean Forte	Dr. Weigert GmbH & Co. KG	0,5%	10 min
Korsolex Extra	Bode Chemie GmbH	5%	15 min

Manual cleaning by immersion:

Step No. /Description – Cleaning procedure

1. Prepare at least 750ml of a cleaning solution: **0,5% Neodisher MediClean Forte**
2. Immerse the flow sensor completely and bubble-free in the **cleaning solution** for at least **10 minutes**. Ensure that the cleaning solution reaches all internal and external surfaces.
3. Agitate the flow sensor in the **cleaning solution** for at least **1 minute**.
Remove and place the flow sensor on a clean surface and leave it to drain completely.

4. Rinse the flow sensor in minimum of **1L of water** (preferably demineralized water);
Swirl/move the flow sensor back and forth for at least **1 minute**.
Remove and place the flow sensor on a clean surface and leave it to drain completely.

5. Inspect the flow sensor for visible soiling and repeat the cleaning process if necessary, (until no more soiling can be seen). Inspect the flow sensor, the sensing wires and the sensor pins for visible damages and check the plastic parts for cracks and deformations.

Manual disinfecting by immersion:

Step No. / Description – Disinfecting procedure	
1. Prepare at least 500ml of a disinfectant: 5% Korsalex Extra	
2. Immerse the flow sensor completely and bubble-free in the disinfectant solution for at least 15 minutes. Ensure that the disinfectant solution reaches all internal and external surfaces.	
3. Agitate the flow sensor in the disinfectant solution for at least 1 minute. Remove and place the flow sensor on a clean surface and leave it to drain completely.	
	
4. Rinse the flow sensor in minimum of 1L of water (preferably demineralized water); Agitate the flow sensor for at least 1 minute. Remove and place the flow sensor on a clean surface and leave it to drain completely.	
	
 Repeat the rinsing process TWO (2) more times with fresh sterile water. THREE (3) separate washes are required! Discard the water after each rinse!	
5. After final rinse, remove and place the flow sensor on a clean surface and then leave it to dry in air for at least 30 minutes.	
 After cleaning and disinfection, allow the flow sensor to dry in air for at least 30 minutes. Ensure that there are no residual vapors or deposits of easily flammable disinfectants. Risk of fire!	
 Inspect the flow sensor for any visible damage or soiling. Inspect the flow sensor, the sensing wires and the sensor pins for visible damages and check the plastic parts for cracks and deformations.	

Sterilization/Steam autoclaving of the flow sensor

SpiroTrue®A flow sensor cannot be autoclaved or sterilized.

Disposal

The flow sensor must be disposed of as infectious waste.

Technical information

Measuring principle:	Hot wire anemometry
Measuring range:	0 to 200 L/min, depending on ventilator (evaluation electronics), gas type and operating conditions.
Accuracy:	DIN EN ISO 80601-2-12, DIN EN ISO 80601-2-13, Typically $\pm 8\%$ of measured value under calibration conditions depending on ventilator (evaluation electronics)

Recommended operating conditions:

Temperature:	10°C to 40°C
Relative humidity:	5% to 95%, non-condensing
Atmospheric pressure:	700hPa to 1060hPa

Recommended storage and transport conditions:

Temperature:	-20°C to 50°C
Relative humidity:	5% to 95%, non-condensing
Atmospheric pressure:	700hPa to 1060hPa

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SpiroTrue®A Explanation of symbols

	de: Katalognummer en: Catalogue number fr: Référence du catalogue es: Número de catálogo	it: Numero di catalogo el: Αριθμός καταλόγου ja: カタログ番号	ro: Număr de catalog tr: Katalog numarası ru: Номер по каталогу
	de: Chargennummer en: Batch code fr: Code de lot es: Código de lote	it: Codice lotto el: Κωδικός παρτίδας ja: バッチコード	ro: Codul lotului tr: Parti kodu ru: Номер партии
	de: Unique Device Identification (einmalige Produktkennung) en: Unique Device Identification fr: Identification unique du dispositif es: Identificación del dispositivo única	it: Identificativo unico del dispositivo el: Μοναδικό αναγνωριστικό συσκευής ja: 機器固有識別子	ro: Identificarea unică a dispozitivului tr: Benzersiz Cihaz Tanımı ru: Уникальный идентификатор прибора
	de: Medizinprodukt en: Medical device fr: Dispositif médical es: Dispositivo médico	it: Dispositivo medico el: Ιατρική συσκευή ja: 医療器具	ro: Dispozitiv medical tr: Tıbbi cihaz ru: Медицинское изделие
	de: Stückzahl en: Number of units fr: Nombre d'unités es: Número de unidades	it: Numero di unità el: Αριθμός μονάδων ja: 個数	ro: Număr de unități tr: Birim sayısı ru: Количество
	de: Herstellungsdatum en: Date of manufacture fr: Date de fabrication es: Fecha de fabricación	it: Data di produzione el: Ημερομηνία κατασκευής ja: 製造年月日	ro: Data fabricației tr: Üretim tarihi ru: Дата изготовления
	de: Ablaufdatum en: Expiry date fr: Date d'expiration es: Fecha de caducidad	it: Data di scadenza el: Ημερομηνία λήξης ja: 有効期限	ro: Data expirării tr: Son kullanma tarihi ru: Дата истечения срока годности
	de: Hersteller en: Manufacturer fr: Fabricant es: Fabricante	it: Produttore el: Κατασκευαστής ja: 製造者	ro: Producător tr: Üretici ru: Изготовитель

 0482	de: CE Europäische Konformitätskennzeichnung. Kennnummer der benannten Stelle. en: CE-European Conformity authorization mark. Identification number of the notified body. fr : Marquage d'autorisation de conformité européenne CE. Numéro d'identification de l'organisme notifié. es: Marca de autorización de Conformidad Europea CE. Número de identificación del organismo notificado.	it: Marcatura di autorizzazione CE - Conformità Europea. Numero di identificazione dell'organismo notificato. el: Σήμα ευρωπαϊκής συμμόρφωσης CE. Αριθμός αναγνώρισης κοινοποιημένου οργανισμού. ja: CE (European Conformity) 認定マーク 通知機関の ID 番号	ro: Marcaj de autorizare CE – Conformitate Europeană. Numărul de identificare al organismului notificat. tr: CE Avrupa Uygunluk yetki işareti. Onaylı kuruluşun kimlik numarası. ru: Знак соответствия нормам Евросоюза. Идентификационный номер органа технической экспертизы.
	de: Gebrauchsanweisung beachten en: Consult instructions for use fr: Consulter les instructions d'utilisation es: Consulte las instrucciones de uso	it: Consultare le istruzioni per l'uso el: Συμβουλευτείτε τις οδηγίες χρήσης ja: 取扱説明書を参照	ro: Consultați instrucțiunile de utilizare tr: Kullanım talimatlarına bakın ru: См. инструкцию по применению
	de: Die Warnhinweise in der Gebrauchsanweisung beachten! en: Observe the warnings in the instructions for use fr: Respecter les avertissements figurant dans le mode d'emploi es: Siga las advertencias de las instrucciones de uso	it: Attenersi alle avvertenze contenute nelle istruzioni per l'uso el: Λάβετε υπόψη τις προειδοποιήσεις στις οδηγίες χρήσης ja: 取扱説明書に記載されている警告を順守すること	ro: Respectați avertismentele din instrucțiunile de utilizare tr: Kullanım talimatlarındaki uyarılara dikkat edin ru: Соблюдать предупреждения в инструкции по применению
	de: Nicht steril en: Non-sterile fr: Non stérile es: No estéril	it: Non sterile el: Μη αποστειρωμένο ja: 滅菌不可	ro: Nesteril tr: Steril değildir ru: Нестерильно

	de: Bei geöffneter oder beschädigter Verpackung nicht mehr verwenden. en: Do not use if package is damaged and consult instructions for use. fr: Ne pas utiliser si l'emballage est endommagé et consulter les instructions d'utilisation. es: No utilizar si el embalaje está dañado y consultar las instrucciones de uso.	it: Non utilizzare se l'imballaggio è danneggiato e consultare le istruzioni per l'uso. el: Μην χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιές και συμβουλευτείτε τις οδηγίες χρήσης. ja: パッケージが損傷している場合は使用せず、取扱説明書を参照すること	ro: Nu utilizați dacă ambalajul este deteriorat și consultați instrucțiunile de utilizare. tr: Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına bakın. ru: Не использовать в случае повреждения упаковки и см. инструкцию по применению.
	de: Bei der Herstellung wurde kein DEHP verwendet en: Does not contain Phthalates fr: Ne contient pas de phtalates es: No contiene ftalatos	it: Non contiene ftalati el: Δεν περιέχει φθαλικές ενώσεις ja: フタル酸エステル類不使用	ro: Nu conține ftalați tr: Ftalat içermez ru: Не содержит фталаты
	de: Bei der Herstellung wurde kein Naturkautschuk-latex verwendet. en: Does not contain Latex fr: Ne contient pas de latex es: No contiene látex	it: Non contiene lattice el: Δεν περιέχει λατέξ ja: ラテックス不使用	ro: Nu conține latex tr: Lateks içermez ru: Не содержит латекс
	de: Von Hitze und Sonnenlicht fernhalten en: Keep away from heat and sunlight fr: Conserver à l'abri de la lumière du soleil et de la chaleur es: Mantener alejado de la luz solar y del calor	it: Non esporre alla luce diretta del sole e al calore el: Φυλάσσετε μακριά από τη θερμότητα και το ηλιακό φως ja: 高温および日光を避けること	ro: A se proteja de lumina solară și căldură tr: Isı ve güneş ışığından uzak tutun ru: Не подвергать воздействию тепла и солнечных лучей
	de: Lagertemperaturbegrenzung en: Storage temperature limitation fr: Limite de température de stockage es: Limitación de la temperatura de almacenamiento	it: Limite di temperatura di conservazione el: Περιορισμός θερμοκρασίας αποθήκευσης ja: 保管温度限界	ro: Limita temperaturii de depozitare tr: Depolama için sıcaklık sınırlaması ru: Предельная температура хранения

	de: Feuchtigkeitsgrenzwerte en: Humidity limitation fr: Limite d'humidité es: Limitación de la humedad	it: Limite di umidità el: Περιορισμός υγρασίας ja: 湿度限界	ro: Limită de umiditate tr: Nem sınırlaması ru: Предельная влажность
	de: Grenzwerte des atmosphärischen Drucks en: Atmospheric pressure limitation fr: Limite de pression atmosphérique es: Limitación de la presión atmosférica	it: Limite di pressione atmosferica el: Περιορισμός ατμοσφαιρικής πίεσης ja: 大気圧限界	ro: Limită de presiune atmosferică tr: Atmosferik basınç sınırlaması ru: Предельное атмосферное давление
	de: Nicht mit dem Hausmüll entsorgen en: Do not dispose of with domestic waste fr: Ne pas éliminer avec les ordures ménagères es: No desechar con la basura doméstica	it: Non smaltire insieme ai comuni rifiuti domestici el: Μην απορρίπτετε σε οικιακά απορρίμματα ja: 家庭ごみと一緒に は廃棄しないこと	ro: Nu aruncați împreună cu gunoiul menajer tr: Evsel atıklar ile birlikte bertaraf etmeyin ru: Запрещается утилизация с бытовыми отходами

SpiroTrue®A – Instructions for use

Language version: multilingual (Part 1: DE/EN/FR/ES/IT/ET/JA/RO/TR/RU)

Information within this document may be subject to changes. Amendments and alterations shall remain reserved without prior notices.

Copyright © bluepoint MEDICAL GmbH & Co. KG.

The content and works published in this document are protected by copyright.

All rights and liabilities are reserved. Without written consent of the bluepoint MEDICAL GmbH & Co. KG the information, work and intellectual property contained in this document must not under any circumstances be circulated, reproduced, duplicated, copied or translated in other languages.

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