Instructions for use: 4000 Series Pulse Oximetry Finger Sensors

DIRECTIONS FOR USE

Intended Use: Viamed 4000 Series oximetry sensors are indicated for use in continuous noninvasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters and patient monitors. Select the sensor model number that matches the brand of oximeter or patient monitor that you are using.

Precautions: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. Do not autoclave or immerse the sensor in liquid. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings: Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor every 8 hours and reposition to an alternate location every 24 hours. Avoid application of the sensors to edematous or fragile tissue. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional hemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

Application of the Sensor

1. Connect the sensor cable to the Viamed adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.
2. Insert the patient's digit into the sensor. The preferred digit to use for adults is the index finger. The digit is correctly inserted when:
 - The tip of the digit touches the rear guide posts.
 - The sensor cable extends along the top of the patient's hand.
3. Visually monitor the sensor site to ensure the integrity of the skin.
4. Secure the lead with medical tape.

Cleaning the Sensor

To clean the sensor between uses or between patients, wipe the sensor surfaces with alcohol.

Environmental Requirements:

- Operating temperature range:
10° - 40° C (50°-104° F)
- Storage and transportation temperature range
0° - 50° C (32°-122° F)
- Relative humidity:
0 - 75% non-condensing

SpO₂ Accuracy:

70 - 100%: A_{rms} <3.0%
0 - 69%: Unspecified

Accuracy test methods available from the manufacturer upon request

LED Specifications:

LED	Wavelength	Radiant Power (typical)	Brand
Red	660 nm	3.5 mW	All below
Infrared	885 nm	3.5 mW	Nellcor, CSI
	890 nm	3.5 mW	HP/Philips
	905 nm	3.5 mW	Datex, BCI, Viamed
	910 nm	3.5 mW	Nonin
	940 nm	3.5 mW	Ohmeda

Disposal:

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed.

All registered trademarks recognized.

Viamed Model	Oximeter Brand	Cable Length
SF4001VM	Nellcor	0.90 m
SF4002VM	Nellcor	2.00 m
SF4006VM	HP/Philips	2.95 m
SF4007VM	HP/Philips	2.00 m
SF4008VM	HP/Philips	2.95 m
SF4011VM	Datex	0.90 m
SF4012VM	Datex	2.95 m
SF4013VM	GE/Datex	2.95 m
SF4016VM	Ohmeda	2.95 m
SF4017VM	Ohmeda	1.00 m
SF4021VM	BCI	1.20 m
SF4026VM	Nonin	0.90 m
SF4036VM	CSI	0.90 m
SF4037VM	CSI	2.95 m
SF6500VM	Viamed	1.20 m

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Viamed Ltd