Instructions for use: 4000 Series Autoclavable Pulse Oximetry Wrap Sensors

DIRECTIONS FOR USE

Intended Use: Viamed 4000 Series oximetry sensors are indicated for use in continuous non-invasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters and patient monitors. Select the sensor model number that matches the brand of oximeter or patient monitor that you are using.

Precautions: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings: Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor every 8 hours and reposition to an alternate location every 24 hours. Avoid application of the sensors to oedematous or fragile tissue. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional haemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

Application of the Sensor:

1. Connect the sensor cable to the Bluepoint Medical adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.
2. Select the sensor site on the patient. The preferred sensor sites are the index finger for adults and paediatrics, the hand for infants and on the foot below the toes for neonates.
3. Apply the sensor on the patient:
 - Position the emitter and detector directly opposite each other.
 - Apply the sensor firmly but not too tightly.
4. Secure the sensor using an adhesive or foam wrap.
5. Visually monitor the sensor site to ensure the integrity of the skin.
6. Secure the lead with medical tape.

Adult / Paediatric

Infant

Neonatal

Cleaning the Sensor

To clean the sensor between uses or between patients, wipe the sensor surfaces with alcohol.

Autoclaving the Sensor

The sensor can be autoclaved at 134°C / 2bar and thereby effectively sterilized.

Environmental Requirements

- Operating temperature range:
10° - 40° C (50° -104° F)
- Storage and transportation temperature range:
0° - 50° C (32° -122° F)
- Relative humidity:
0 - 75% non-condensing

S_PO₂ Accuracy

70 - 100%: A_{rms} <3.0%

0 - 69%: Unspecified

Accuracy test methods available from Viamed upon request.

LED Specifications

LED	Wavelength	Radiant Power (typical)	Brand
Red	660 nm	3.5 mW	All below
Infrared	890 nm	3.5 mW	HP/Philips
	905 nm	3.5 mW	Viamed

Disposal

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed.

All registered trademarks recognized.

Viamed Model	Oximeter Brand	Cable Length
WA4007VM	HP/Philips	2.00 m
WA4008VM	HP/Philips	2.95 m
WA6500VM	Viamed	1.20 m

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Bluepoint Medical GmbH & Co. KG