Instructions for use: 4000 Series Pulse Oximetry Soft Sensors

DIRECTIONS FOR USE



Intended Use: Viamed 4000 Series oximetry sensors are indicated for use in continuous noninvasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters and patient monitors. Select the sensor model number that matches the brand of oximeter or patient monitor that you are using.

Precautions: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. Do not autoclave. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings: Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor every 8 hours and reposition to an alternate location every 24 hours. Avoid application of the sensors to edematous or fragile tissue. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional hemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

Application of the Sensor

1. Connect the sensor cable to the Viamed adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.
2. Insert the patient's digit into the sensor. The preferred digit to use is the index finger. The digit is correctly inserted when:
 - The tip of the digit touches the rear raised lip.
 - The sensor cable extends along the top of the patient's hand.
3. Visually monitor the sensor site to ensure the integrity of the skin.
4. Secure the lead with medical tape.

Adult

Paediatric

Cleaning the Sensor

To clean the sensor between uses or between patients, wipe the sensor surfaces with alcohol.

Environmental Requirements:

- Operating temperature range:
5° - 40° C (41°-104° F)
- Storage and transportation temperature range
-20° - 70° C (-4°-158° F)
- Relative humidity:
15 - 90% non-condensing

S_PO₂ Accuracy:

70 - 100%: A_{rms} <3.0%

0 - 69%: Unspecified

Accuracy test methods available from Viamed upon request.

LED Specifications

LED	Wavelength	Radiant Power (typical)	Brand
Red	660 nm	3.5 mW	All below
Infrared	885 nm	3.5 mW	Nellcor, CSI
	890 nm	3.5 mW	HP/Philips
	905 nm	3.5 mW	Datex, BCI, Viamed
	910 nm	3.5 mW	Nonin
	935 nm	3.5 mW	Nihon-Kohden

Disposal

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed.

All registered trademarks recognized.

Viamed Model - Adult	Viamed Model - Paediatric	Oximeter Brand	Cable Length
SC4001VM	SCP4001VM	Nellcor	0.90 m
SC4002VM	SCP4002VM	Nellcor	2.00 m
SC4006VM	SCP4006VM	HP/Philips	2.95 m
SC4007VM	SCP4007VM	HP/Philips	2.00 m
SC4008VM	SCP4008VM	HP/Philips	2.95 m
SC4011VM	SCP4011VM	Datex	0.90 m
SC4012VM	SCP4012VM	Datex	2.95 m
SC4013VM	SCP4013VM	GE/Datex	2.95 m
SC4021VM	SCP4021VM	BCI	1.20 m
SC4026VM	SCP4026VM	Nonin	0.90 m
SC4046VM	SCP4046VM	Nihon-Kohden	2.95 m
SC4047VM	SCP4047VM	Nihon-Kohden	1.60 m
SC6500VM	SCP6500VM	Viamed	1.20 m
SC6501VM	SCP6501VM	Viamed	0.30 m

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Viamed Ltd