SPIROQUANT H FLOW SENSOR

Instructions for Use / Mode d'emploi / Gebrauchsanweisung

FLOW SENSOR

The SpiroQuant H is a single patient use sensor used to measure the airflow in the breathing system of the patient for use with Hamilton ventilators.

Ventilator connection

To connect the Flow Sensor to the ventilator attach the sensor tubes underneath the front panel of the ventilator. The blue tube must be attached to the „PROXIMAL“ inlet of the ventilator.

Patient connection

Insert the Flow Sensor between the Y piece of the patient circuit and the patient connection. The side of the Flow Sensor equipped with the blue tube should face the patient. It is recommended that short flexible tubing is used between the sensor and the patient connection. Position the sensor with the tubing upright to prevent kinking and condensation.

Tubing clip

The tubing clip is designed to secure the Flow Sensor tubes to the ventilation tubes. Place the wide opening of the clip around the inspiratory ventilation tube. Place the Flow Sensor tube in the narrow opening of the clip, leaving a loop of tubing between the clip and the Flow Sensor to accommodate the full range of patient movement.

Calibration

Caution: The Flow Sensor must be calibrated prior to patient use.

To calibrate, follow the instructions in the appropriate ventilator manual. After successful completion of the procedure the Flow Sensor is ready for patient use.

Caution: Federal law (USA) restricts this device to sale by or on order of a licensed physician.

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Consult the ventilator's manual for detailed information!

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RETURNING THE PRODUCT


Contact the customer service department at Honeywell Healthcare Solutions GmbH to obtain an RMA number (Return Material Authorisation). Every product sent to us for checking must have been cleaned and disinfected beforehand.

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	Follow the instructions given in the operating manual
	Follow the Instructions
	Manufacturing Date
	Manufacturer
	Dispose of responsibly
	Batch Code
	Single use / Single patient use
	Medical Device
	CE-Symbol with Notified Body Number

Honeywell Healthcare Solutions GmbH

Alter Holzhafen 18, 23966 Wismar,
Germany
Phone: +49 (0)3841-360-1
Phone: +49 (0)3841-360-200
Fax: +49 (0)3841-360-222

**Honeywell Sensing &
Safety Technologies**

830 East Arapaho Road
Richardson, TX 75081
www.honeywell.com