
## Gas Sampling Line H, Disposable Airway Adapter

**en** Instructions for Use

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**Type: Adult/Pediatric, 2m/7ft**

**REF** 8090121315

**Type: Infant, 2m/7ft**

**REF** 8090121316

## Indications for Use

The Gas Sampling Line H is intended for continuous monitoring of expired and inspired respiratory gases. The Gas Sampling Line H is a single patient use device with a cumulative duration of use of up to 72 hours.

The Gas Sampling Line H is for use with sidestream respiratory gas monitors. The Gas Sampling Line H is provided in two sizes: large for use with adult/pediatric patients, and small for use with infant patients, for connection to the patient's breathing circuit.

Gas Sampling Line H is intended for continuous monitoring of Carbon Dioxide (CO<sub>2</sub>), Oxygen (O<sub>2</sub>), Sevoflurane, Isoflurane and respiratory rate in professional/hospital care environments, such as operating rooms and intensive care units (ICU).

The Gas Sampling Line H is intended for use by qualified medical personnel only.

Prescription use only.

## Warnings

**Warning:** Only for use with anaesthetic agents Sevoflurane and Isoflurane. Do not use with any other anaesthetic agents.

**Warning:** Carefully route the sampling line to reduce the risk of patient entanglement or strangulation.

**Warning:** Do not cut or remove any part of the sampling line. Cutting the sampling line could lead to erroneous readings.

**Warning:** Due to risk of cross contamination, do not re-use.

**Warning:** Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

**Warning:** Check tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate measurement.

**Warning:** Do not use adult or pediatric type sampling line configurations with infants, as the adapter adds 6.6ml dead space to the patient circuit.

**Warning:** Do not use infant sampling line configurations with adults or pediatric, as this may result in excessive flow resistance (infant adapter dead space: <1ml).

**Warning:** Do not place the airway adapter (tee) between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.

## Cautions

**Caution:** Do not attempt to clean, disinfect, sterilize or rinse any part of the sampling line with liquids.

**Caution:** Dispose of sampling lines in accordance with standard operating procedures or local regulations for the disposal of contaminated medical waste.

## Notes

**Note:** During nebulization or suction, in order to avoid moisture build up and sampling line occlusion, remove the sampling line connector from the monitor.

**Note:** For accuracy specifications and further instructions, warnings and cautions, please consult the gas monitor instructions for use.

## Contraindication

There are no known contraindications for patient monitoring with the Airway Gas Sampling Products provided that the data obtained by the gas monitoring is evaluated with consideration given to the patients clinical condition.

## Instructions for Use


Figure 1


Figure 2

1. Use the Adult/Pediatric size with endotracheal tubes of inner diameter  $> 4$  mm and the Infant size with inner diameter  $< 4$  mm.
2. Connect the Airway Adapter to the proximal end of the breathing circuit between the endotracheal tube and the ventilator Y-section as shown in Figure 1.
3. Connect the gas sample line (male Luer connector) to the respiratory gas monitor as shown in Figure 2. Rotate the Luer connector clockwise for a secure fit to ensure that there is no leakage of gases at the connection point.
4. After connection of the sampling line, check that gas values appear on the monitor display.
5. For reasons of hygiene, replace sampling line within a maximum period of 72 hours.

## Technical information

Operating temperature:  $+10^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$  ( $50^{\circ}\text{F}$  to  $104^{\circ}\text{F}$ )

Storage temperature:  $-25^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$  ( $-13^{\circ}\text{F}$  to  $140^{\circ}\text{F}$ )

Atmospheric pressure: 660 to 1060 hPa (495 to 795 mmHg)

Ambient Humidity: 10 to 95% RH (non-condensing)

For technical specification, see host device user guide or manual.

## Reporting of serious incidents

If there is a *serious incident* in connection with the use of this product, this should be reported. The incident should be reported to the manufacturer and the health authority or the competent authority for the installation site of the product. A *serious incident* is when there is a **death or a temporary or permanent serious deterioration in the health of a patient, user or other person.**

Send an email to the following manufacturer address: [prrc@bluepoint-medical.com](mailto:prrc@bluepoint-medical.com)

When doing so, please provide the following information:

- Order number and model designation of the product as indicated on the product
- Serial number/batch number of the product
- Date of the *serious incident*
- Description of the *serious incident*, including its impact on the patient or any injury
- Your contact details (institution, address, contact name (substitute), title and phone number)

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bluepoint medical GmbH & Co. KG