Sauerstoffsensoren / Oxygen Sensors

Product: Sauerstoffsensoren / Oxygen Sensors

Description:

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*The device/system encompasses the following **components**:*

- This clinical evaluation refers to the following medical O2 sensors: OOMXXX (Except. OOM109)

Intended Use:

The EnviteC Medical Oxygen Sensors are intended as oxygen-sensing component of an oxygen analyzer that measures oxygen concentration in breathing gas mixtures in the following applications:

- Sensing device for oxygen in control device of oxygen concentrators
- Sensing device for oxygen in medical ventilators
- Sensing device for oxygen in anaesthesia equipment
- Sensing device for oxygen in incubators

The use is limited to system monitoring. The sensors are not suited for breath by breath analysis of breath gases.

Functional Description:

EnviteC Medical Oxygen Sensors are based on the principle of electro-galvanic sensors. They are constructed in a plastic housing containing two electrodes - a precious metal cathode and an anode immersed in a liquid electrolyte medium. Electrically the device resembles a very low voltage battery cell. A gas permeable diffusion membrane provides the interface to the gas sample. The oxygen gas is reduced on the sensing electrode (cathode) and the anode is oxidized. The resulting current produces an external electrical voltage signal on the load resistor that is proportional to the conversion of the oxygen. The sensor signal is temperature dependent and typically compensated with an internal temperature-compensating resistor network.

1 Performance of Clinical Evaluation

The evaluation is based on the following documents:

- European Medical Device Directive 93/42/EEC (last change 2007/47/EEC) Annex X
- MEDDEV2.7.1 Evaluation of Clinical Data: A Guide for Manufacturers and Notified Bodies Coordination of Notified Bodies Medical Device (NB-MED) on Council Directives 90/385/EEC and 93/42/EEC

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The following path is chosen for the clinical evaluation:

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Table 1. Bewertungsteam / Appraisal team

Team of clinical evaluation	Qualification (Factual knowledge about topic under discussion plus length of service, if applicable)
Marko Sprössel	Mechanical Engineer, More than 17 years experience in development of medical products
Maik Lingies	Electrical Engineer, More than 17 years experience in development of medical products
Christa Dumschat	PhD in Chemistry, More than 15 years experience in development of medical products
Eleonore Helms	Quality Manager, More than 10 years experience in automotive, electronics and medical industries Certified Manager Regulatory Affairs Medical Device International Certified 6 Sigma Black Belt

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Remark: Literature search and appraisal of pre- and post market clinical data to be performed only by suitably qualified staff experienced in the design and application of the product.

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This approach of clinical evaluation is based on the compilation and appraisal of clinical data from relevant scientific literature, standards and specifications with subsequent critical evaluation.

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The objective of this evaluation is to identify comprehensive specifications for the medical device, its safe use, performance, design features and intended use as well as environmental conditions. The literature used defines as completely as possible all specifications for the safety and performance of the medical device. No known publications with regard to the use, risks, side effects of the product type nor of its intended use are withheld.

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² Clinical trials must be applied for separately with the competent authorities (DIMDI, ethical review committee etc.) and are subject to special rules for implementation and evaluation.

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Literature source No.	Aufgeführte Anwendung/Zweckbestimmung laut Literaturstelle (Indikation, Kontraindikation) <i>Listed application/intended use from literature source (indication, contra-indication)</i>	Offerierte Anwendung zur Zweckbestimmung des Produktes <i>Offered application for product's intended use</i>
I	RESPIRATORY GAS MONITOR (RGM)	Oxygen-sensing component of an oxygen analyzer that measures oxygen concentration in breathing gas mixtures for system monitoring
IV	Oxygen concentrator equipment	Sensing device for oxygen in oxygen concentrator equipment for system monitoring
X	Oxygen sensor Incubators	Sensing device for oxygen in incubators
III/VI	Oxygen sensor in Anesthesia Equipment	Sensing device for oxygen in anaesthesia equipment for system monitoring
II/IX/XI	Oxygen sensor in critical care ventilators	Sensing device for oxygen in medical ventilators for system monitoring
VII	Oxygen gas analyzer	Oxygen-sensing component of an oxygen analyzer that measures oxygen concentration in breathing gas mixtures for system monitoring

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The intended use of the product and its limitations of use are fully described in the above literature.

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Geforderte Leistungsdaten: <i>Required performance data:</i>	Anforderungen / Werte: <i>Requirements / Values:</i>	Literaturstelle Nr.: <i>Literature source No.:</i>	Leistungsdaten des Produktes: <i>Performance data of the product:</i>
Measurement accuracy O ₂	Volume fraction of 2.5% + 2.5 % of gas level	I	Meets ISO 80601-2-55 requirements
Drift of measurement accuracy	Not more than volume fraction of 2.5% + 2.5 % of gas level in 6h	I	Meets ISO 80601-2-55 requirements
Measurement accuracy of gas readings for gas mixtures	According to ISO 80601-2-55:2011	I	Meets ISO 80601-2-55 requirements
Stability of marking and labels during cleaning and disinfection	according to IEC 60601-1:2012; section 7.1.3	II	According to requirements (see VIII)
Shock and vibration for professional transportation	According to IEC 60600-2-27:2008		According to requirements (see VIII)
Measurement range	18-100% O ₂ V/V	I u. VI	0-100% O ₂ V/V
Humidity	Incubator 30 to 99% ventilator dry air	X; XI	0 – 99% RH non condensing – 0.03 % rel. per % RH at 25°C
Pressure (open circuit)	125 hPa above ambient pressure	II u. III	600 hPa – 2000 (0 – 1500 mBar O ₂)

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Evaluation of the product's known residual risks from citations

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Comparison of residual risks analyzed in the risk management file and their mitigation from the cited literature:

Con- sec. no.	Potential residual risk	Citation no.:	
		VII	VI
1.	Cross infection after use with another patient	No problem in class CCL reported	p.176
2.	Wrong measurement due to temperature or pressure outside of the specification	No problem in class CCL reported	p. 575

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The literature search proves the clinical benefit of the product. The residual risks of the product are minimized to an acceptable level (recognized state of clinical practice).

The selected literature is sufficient to evaluate the product's conformity with the Essential Requirements of the European Directive 93/42/EEC.

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One sensor type from Envitec was chosen and compared with compatible competitive sensors.

Table 6:

Product, performance data	Equivalent product I	Equivalent product II	Comments
Hersteller / Manufacturer:	City Technology	Teledyne Analytical Instruments	EnviteC
Model:	MOX-9	R24	OOM 202 (as example)
Revision der GA <i>Revision of instructions for use</i>	None found	C71746 Rev.4	Rev.1 2015
Revision der Kurzanleitung <i>Revision of brief instruction</i>	Product Data Sheet Issue 6, Feb.14	Spec Control DWG Oxygen sensor Class R24 Rev.10 (23.04.2002) retrieved at Sept. 2015	Product Data Sheet 010/2006
Zulassungs-Nr. <i>Accreditation number</i>	FDA K041773	FDA K001095	FDA K082655
Product code (UMDNS ...)			13-538
Intended use	These sensors are designed to be used to monitor the partial pressure of oxygen in anaesthesia (not including xenon), critical care, neonatal incubators and general oxygen monitors.	Teledyne medical oxygen sensors are intended to accurately measure the concentration of oxygen in a mixture of gases used in medical applications. The sensors are the oxygen-sensing component. Intended for sensor replacement of finished medical devices such as an oxygen monitor that is used to monitor the oxygen concentration in a patient's breathing environment.	The EnviteC Medical Oxygen Sensors are intended as oxygen-sensing component of an oxygen analyzer that measures oxygen concentration in breathing gas mixtures in the following applications: - Sensing device for oxygen in control device of oxygen concentrators - Sensing device for oxygen in medical ventilators - Sensing device for oxygen in anesthesia equipment - Sensing device for oxygen in incubators The use is limited to system monitoring. The sensors are not suited for breath by breath analysis of breath gases.
Patient population	The target patient population consists of those patients who require the oxygen concentration in their breathing environment to be monitored.	Not found	The target patient population consists of those patients who require the oxygen concentration in their breathing environment to be monitored.
Application	See intended use	See intended use	See intended use
...			

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Remark: *Clinical evaluation allows only approved equivalent products with comparable intended use and application range.*

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Product, performance data:	Equivalent product I	Equivalent product II	Product of clinical evaluation	Comments
Safety advice in accompanying documents:	Connection should be made via recommended mating parts only. Soldering to the sensor will damage it and invalidate the warranty. However, it is important that exposure to high concentrations of solvent vapours is avoided, both during storage, fitting into instruments and operation. When using sensors with printed circuit boards (PCBs), degreasing agents should be used before the sensor is fitted. Do not glue directly on or near the CiTiceL as the solvent may cause crazing of the plastic. The sensor is not suitable for sterilisation by steam or exposure to chemicals such as ethylene oxide or hydrogen peroxide. When installing the sensor, it must only be screwed in hand-tight and a gas tight seal ensured. Spanners and similar mechanical aids may not be used, as excessive force may damage the sensor thread. MediceLs contain metal and may be susceptible to RFI or EMI. For further information please contact City Technology.	Inspect the sensor for any physical damage or electrolyte leakage. If the sensor is defective DO NOT USE. The sensor electrolyte is caustic. Do not let it come in contact with the skin or eyes.	It is the responsibility of the user to determine the suitability for use of the sensor. Follow the instructions for use of the oxygen analyzer and for replacement of oxygen sensor. To avoid cross infection please follow strictly the instructions of the oxygen analyzer manufacturer. Refer to the oxygen analyzer operation manual to determine any needed preoperative checks. Sensor contains encapsulated by a housing, lead (Pb), lead oxide (PbO) and concentrated potassium hydroxide solution (between 2 and 5 mol / L). Lead and lead oxide are toxic and dangerous for the environment. Concentrated potassium hydroxide is corrosive (see safety data sheet). Do not open the housing or penetrate the permeable membrane. Do not touch a damaged sensor without protective gloves. In the case of leakage avoid contact with eyes. The sensor is not suited for use in a magnetic resonance imaging (MRI) environment. Cleaning/ Disinfection: The sensor membrane and the printed circuit board should not come in contact with disinfectant or cleaning agent. The other parts of the sensor can be disinfected by disinfectant wipes or with a surface disinfection agent. Follow the instructions of the producer of the disinfection material. Before insertion into the device check the sensor for mechanical damages and for humidity or crystallization of salts on the housing. Do not use a damaged sensor or a sensor with crystallization of salts outside.	

Sauerstoffsensoren / Oxygen Sensors

Safety advice on product:	Not available	On R-17 Med sensor from 2009 Pouch: The product contains a chemical known to the state of California to cause cancer or reproductive harm. On Label: Contains caustic. Causes severe burns. Maybe fatal if swallowed.	 The sensor contains caustic liquid. In case of leakage avoid contact with eyes and skin.	
Labelling, type label:	Not available	See below	See below	
...				

Label Envitec
ENVITEC

 EnviteC-Wismar GmbH
 www.envitec.com CE 0123

Medical Oxygen Sensor

 TYPE:
OOM202

SN: A246353

2014-09

Label Teledyne

TELEDYNE ANALYTICAL INSTRUMENTS CE 0086		OXYGEN SENSOR CLASS R-17MED U.S. PAT. NO. 5,714,046	
16830 CHESTNUT ST. INDUSTRY, CA 91746-1020 (626) 934-1500 USA	WARNING! CONTAINS CAUSTIC. CAUSES SEVERE BURNS! MAY BE FATAL IF SWALLOWED	S/N	350555 A9 Made in USA

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Product, performance data:	Equivalent product I	Equivalent product II	Product of clinical evaluation	Comments
Measuring range:	0-100% 0-1500 mBar O ₂	0-100%	0-100% (ppO ₂ 0 - 1250 mBar O ₂)	
Measuring accuracy:	Not specified	Within $\pm 2\%$ @ 60% O ₂ ; $\pm 3\%$ @ 100% O ₂ at constant temperature and pressure when calibrated in air	Meets ISO 80601-2-55 requirements	
Output	10-13 mV in 210 mBar O ₂	10-17 mV in air, sea level at 25° C	13-16 mV in ambient air	
Response time (T ₉₀)	15 s (air to 100% O ₂)	10 s	12 s	
Offset	<200µV	Less than 0.5% of oxygen equivalent in zero gas after 60 s	150µA	
Temperature compensation	< 2% O ₂ equivalent (0-40°C)	5% of full scale over the operating temperature range, worst case tracking error (within the first hour after a maximum temperature step) is 7.5% full scale	Between +25° C and +40° C: 3% relative error between 0° C and +50° C: 8% relative error	
Material, design:				
Standards and classifications:	ISO 60601-2-55	Not specified	ISO 60601-2-55	
EMC:	Not applicable	Not applicable	Not applicable	
Environmental conditions: Operation	-20°C to 50°C 0.5 to 2.0 Bar 0-99% RH non-condensing	0°C to 40°C 0-99% RH non-condensing	0°C to 50°C 0.6 to 2.0 Bar 0-99% RH non-condensing	
Environmental conditions: Storage, transport	-10°C to 40°C (short excursions to 50°C allowed)	0°C to 50°C	-20°C to 50°C (recommended 5°C to 15°C)	
Warranty	13 months from date of dispatch		15 months	
Expected Cell Live	0.9 x 106% O ₂ hours at 20° C 0.6 x 106% O ₂ hours at 40° C	24 months in air at 25°C and 50% R.H.	>1.000 000 % vol. oxygen hours	

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Medical O₂ sensors are produced and sold for 20 years by ENVITEC. All complaints of that product group from 2005 until May 2015 were reviewed (since that time our complaints are recorded in a database). No cases were identified that resulted in a hazardous situation or could have resulted in a hazardous situation. Due to the long-term market experience without hazardous situations, the residual risk is evaluated as acceptable.

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Based on the above clinical evaluation the Product complies with the Essential Requirements of the European Directive 93/42/EEC. The Product is to be considered safe in the scope of the Directive for the above application and the intended use as described. The residual risks attributed to the Product are acceptably low compared with the Product's clinical benefit.

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