

# **ELECTRONIC TEMPERATURE MONITORS**

**(INCLUDING CALIBRATOR)**

---

## **SERVICE INFORMATION**

---

**TYPES:** ET 100 SINGLE CHANNEL  
ET 200 TWIN CHANNEL  
TC 100 CALIBRATOR

**CONTENTS:** INTRODUCTION  
CIRCUIT DIAGRAMS  
COMPONENT LAYOUT  
COMPONENT LIST  
CALIBRATION PROCEDURE



**LIBRA MEDICAL**


**138 Overdown Road, Reading, Berkshire RG3 6NJ  
Telephone: 0734 25235**

## INTRODUCTION


These instruments, in single and twin channel versions, have been designed to provide accurate and reliable long term temperature monitoring and display in any clinical situation. They have been designed to accept thermistor temperature probes with a specification similar to the Yellow Springs 400 series. The information contained in this brochure should enable any technical service department to carry out repairs or recalibration if required. Calibration can be checked with the in built calibration facility, at 40°C, and an external calibrator covering 22 – 50°C in six steps, is available for routine checks, Part No. TC 100.

## CIRCUIT DIAGRAMS AND COMPONENT LAYOUT


1 & 2 CHANNEL: DISPLAY BOARD




1 CHANNEL: ANALOGUE BOARD




1 CHANNEL: ANALOGUE BOARD LAYOUT



1 & 2 CHANNEL: DISPLAY BOARD LAYOUT



2 CHANNEL: ANALOGUE BOARD



## COMPONENT LIST

### DISPLAY BOARD

R26	390K	1%	.25W	DISP 1	LCD 3½ digit
R27	47K	1%	.25W		.5in high.
C6	47pF	ceramic		19993	Knob
C7, 8	.47µF	met.poly.		19994	Case – white
C9	.22µF	met.poly.		19995	ET100 front panel
IC2	ICL 7126 CPL			19996	ET200 front panel

### SINGLE CHANNEL ANALOGUE BOARD

R1	6K34	.5%	.25W	PR1, 2, 3, 4	1K	preset
R2, 4, 7	100K	1%	.25W	C1		.47µF met.poly.
R3	1K	1%	.25W	C2		2.2µF tant.
R5	3K3	1%	.25W	C3		.1µF met.poly.
R6	470R	1%	.25W	C4		100µF 16v al. elect.
R8	5K9	1%	.25W	D1		ICL 8069 DCZR
R9	330R	1%	.25W	Q1		2N6027
R10	18K	1%	.25W	SW1		switch: 6way 2pole
R11	1K2	1%	.25W			
R12	1M	1%	.25W			
R13	11K3	1%	.25W			
R15	270K	1%	.25W			
R16	4K7	1%	.25W			
R14	470K	1%	.25W			

### TWIN CHANNEL ANALOGUE BOARD

R1	6K34	.5%	.25W	PR1, 2, 3,	1K	preset
R2, 4, 7	100K	1%	.25W	4, 6, 7		
R3	1K	1%	.25W	PR5	22K	preset
R5	3K3	1%	.25W	C1, 2		.47µF met.poly.
R6	470R	1%	.25W	C3		100µF 16V al.elect.
R8	5K9	.5%	.25W	C4		.1µF met.poly.
R9	330R	1%	.25W	C5		2.2µF tant.
R10	18K	1%	.25W	IC1		ICL 7611 DCPA
R11, 12	1K2	1%	.25W	D1		ICL 8069 DCZR
R13	15K	1%	.25W	Q1		2N6027
R14	2K7	1%	.25W	SW1		switch: 6way 2pole
R15, 16,						
17, 18	1M	1%	.25W			
R19, 20,						
21	4K7	1%	.25W			
R22	270K	1%	.25W			
R24,25	11K3	.5%	.25W			
R23	470K	1%	.25W			

## **CALIBRATION PROCEDURE**

<b>CALIBRATION VALUES (°C)</b>						
Calibration	22	26	32	38	44	50
Thermometer display	22.0	25.9	31.9	38.0	43.9	49.8
Critical values				●		

### **SINGLE CHANNEL**

1. Plug calibration box into probe input. Set all pots (PR1 – 4) halfway. Switch thermometer to USE and set box to 22.0
2. Adjust Zero (PR3) so display reads 22.0
3. Set box to 50° and adjust Max (PR1) so display reads 49.8
4. Repeat 2. & 3. until satisfactory.
5. Set box to each position and check display against table above. Make final adjustments if necessary.
6. Disconnect calibration box. Switch thermometer to CAL.
7. Adjust Cal (PR2) so display reads 40.0
8. Switch thermometer to USE. Adjust 0° (PR4) so display reads 0.0

### **TWIN CHANNEL**

1. Plug calibration boxes into T1 & T2 inputs. Set all pots halfway. Switch thermometer to T1 and set T1 calibration box to 22.0
2. Adjust T1 zero (PR3) so display reads 22.0
3. Set box to 50° and adjust Max (PR1) so display reads 49.8
4. Repeat 2 & 3 until satisfactory.
5. Set box to each position and check display against table above. Make final adjustments if necessary.
6. Switch thermometer to T2 and set T2 calibration box to 22.0
7. Adjust T2 zero (PR4) so display reads 22.0
8. Repeat 5 but only adjust PR4.
9. Disconnect calibration boxes. Switch thermometer to CAL.
10. Adjust CAL (PR2) so display reads 40.0
11. Switch thermometer to T1. Adjust T1 0° (PR6) so display reads 0.0
12. Switch thermometer to T2. Adjust T2 0° (PR7) so display reads 0.0
13. Plug in calibration boxes, set to 38°. Switch thermometer to  $\Delta T$ . Adjust  $\Delta T$  zero (PR5) until display reads 0.0
14. Leave T2 box at 38°. Set T1 box to each position and check that display error is less than 0.2°. Set T1 box to 38° and repeat for T2.

### **BATTERY LOW VOLTAGE: BOTH MODELS**

1. Remove PP3 battery from thermometer and replace with bench type DC power supply. Also connect digital voltmeter across supply for more accurate reading. Set supply to 9V and switch on. Switch thermometer to CAL.
2. Slowly reduce supply voltage until display flashes. This should occur in the range 7.2 – 6.6V.