**VIAMED SMARTSAT®
DISPOSABLE PULSE OXIMETRY SENSORS**

Intended Use: SMARTsat® Disposable Pulse Oximetry Sensors, models 10-AP-VM, 10-PP-VM, 10-IP-VM, 10-NP-VM, are intended for use with Viamed monitoring systems that use SMARTsat® technology.

The SMARTsat® Disposable Pulse Oximetry Sensors are indicated for use in non-invasive continuous monitoring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate on adults to neonates.

The use environment may include: hospitals, hospital-type facilities and intra-hospital transport.

The use environment is dependent upon the use environments specified for the monitoring system utilized in conjunction with the sensor.

Note: This Instruction for Use is intended for the medical professional.

Caution: The SMARTsat® Disposable Pulse Oximetry Sensors are for prescription use only.

Caution: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter.

Warnings

1. A functional tester (like Index II or equivalent) cannot be used to assess the accuracy of pulse oximeter probes.
2. Connect the sensor only to a SMARTsat® technology monitor or module. A mismatch may cause patient injury, inaccurate measurements or equipment damage.
3. Check the patient skin condition at least every 4 hours and reposition to an alternative location if skin integrity changes.
4. Do not reuse disposable sensors as this can cause cross infection. In addition the probe function can be out of specification upon reuse.
5. A misuse of a pulse oximeter probe with excessive pressure and with prolonged periods can cause pressure injury.
6. Avoid application of the sensors to oedematous or fragile tissue.
7. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure.
8. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional haemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement.
9. Route all cables carefully to reduce the possibility of patient entanglement or strangulation.
10. Note: The materials used in the manufacture of the sensor contain no natural rubber latex protein. The materials have undergone extensive biocompatibility testing. Further information is available upon request.
11. Do not use this sensor with MRI scanning. Burns may be caused by conducted current.

Application of the Sensor

1. Select the sensor site on the patient. The preferred sensor sites are the index finger for adults and paediatrics, the big toe for infants and on the foot below the toes for neonates.
2. Remove the protective film from the adhesive surface of the sensor **(A)**.
3. Apply the sensor on the patient **(B)**:
 - Position the emitter and detector directly opposite each other
 - Apply the sensor firmly, but not too tightly
4. Visually monitor the sensor site to ensure the integrity of the skin.
5. Secure the lead with medical tape.
6. Connect the sensor cable to the adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.

Environmental Requirements

- Operating temperature range:
+5°C to +40°C (+41°F to +104°F)
- Storage and transportation temperature range:
-20°C to +60°C (-4°F to +140°F)
- Storage relative humidity range:
15 to 95% non-condensing

SpO₂ Accuracy:

60 – 100%:	A _{rms} < 3%
90 – 100%:	A _{rms} < 2%
80 – 90%:	A _{rms} < 2%
70 – 80%:	A _{rms} < 3%
60 – 70%:	A _{rms} < 3%
<60%:	Unspecified

Pulse Rate Accuracy:

20 – 300 BPM: A_{rms} < 2 BPM
BPM: Beats per minute

Pulse oximeter equipment measurements are statistically distributed. Only about two-thirds of measurements can be expected to fall within ± Arms of the measured values by a co-oximeter. To verify the function of pulse oximeter probes, a functional tester (like Index II or equivalent) can be used.

SpO₂ accuracy was validated by clinical tests, where the measured values of the sensors were compared with those of the arterial co-oximetry in healthy adult subjects over the specified functional oxygen saturation range. These values were validated on the following monitor: VM-2160 / Oxytrue®A (SMARTsat®OEM III). A functional tester may not be used to validate SpO₂ accuracy.

LED Specifications*

*This information is especially useful for clinicians

LED	Wavelength	Radiant Power (typical)
Red	660 nm (Peak)	2.5 - 3.5 mW
Infrared	905 nm (Centroid)	2.5 - 3.5 mW

Disposal

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed Ltd.

REF	Model	Viamed Type	Intended patient
6020131204	10-AP-VM	0015010	Adult
6020131207	10-PP-VM	0015011	Paediatric
6020131209	10-IP-VM	0015012	Infant
6020131211	10-NP-VM	0015013	Neonatal

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© 2018 bluepoint medical GmbH & Co. KG
SMARTsat® and OxyTrue® are registered trademarks of bluepoint medical GmbH & Co. KG.
bluepoint medical GmbH & Co. KG, An der Trave 15, 23923 Selmsdorf, Germany