
SMARTsat®

Reusable Pulse Oximetry Sensors


en Instructions For Use

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**SC7500VM, SCM7500VM,
SCP7500VM**

SF7500VM


W7500VM


**SMARTsat®
Technology**


(A)


SC7500VM, SCM7500VM
(>20kg)


SCP7500VM (10 - 20kg)


SF7500VM (>20kg)

W7500VM (>20kg)


W7500VM (10-20kg)


W7500VM (<3kg)

(B)


SC7500VM, SCM7500VM,
SF7500VM


SCP7500VM, W7500VM

(C)


SC7500VM,
SCM7500VM,
SCP7500VM

(D)


W7500VM

(E)


SC7500VM, SCM7500VM, SCP7500VM

(F)


SC7500VM, SCM7500VM, SCP7500VM

(G)

(H)

XT6500VM/ XT6501VM → ← SC7500VM/ SCM7500VM/ SCP7500VM/ SF7500VM/ W7500VM

IFU-01-02-0004 Rev.0 2018-06

SMARTsat® and OxyTrue® are registered trademarks of bluepoint medical GmbH & Co. KG.

™* Third party brands are trademarks of their respective owners.

© 2018 bluepoint medical GmbH & Co. KG.

bluepoint medical GmbH & Co. KG

An der Trave 15

D-23923 Selmsdorf, Germany

www.bluepoint-medical.com

Viamed SMARTsat®

Reusable Pulse Oximetry Sensors

Directions for Use

Note: This Instruction for Use is intended for the medical professional. For use in the home care environment the user may use the device after instruction by, or under supervision of, a trained healthcare professional.

Indications for Use

The SMARTsat® Reusable Pulse Oximetry Sensors Model SC7500VM, SCM7500VM, SCP7500VM, SF7500VM and W7500VM are intended for use with Viamed monitoring systems that use SMARTsat® technology.

The SMARTsat® Reusable Pulse Oximetry Sensors are intended for use in non-invasive continuous or spot-check monitoring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate.

The use environment may include: hospitals, hospital-type facilities, intra-hospital transport, mobile emergency medical applications and home use.

Transport environments include intra-hospital transport and ground emergency transport (road ambulances). The use environment is dependent upon the use environments specified for the monitoring system utilized in conjunction with the sensor.

The SMARTsat® Reusable Pulse Oximetry Sensors are for prescription use only.

Disposal

Do not dispose of with domestic waste. Follow local regulations when disposing or recycling this device and its component parts. A list of component materials is available by contacting Viamed Ltd. (see back page).

Instructions for Use

Sensor application:

1. *Apply SC7500VM, SCM7500VM, SCP7500VM, SF7500VM:*

Insert the patient's finger into the sensor. The preferred finger to use is the index finger (A). Alternative sites may be used depending on finger sizes (B).

Apply W7500VM:

Fix a Wrap-Tape to the sensor (D) and apply the sensor to the selected sensor site. The preferred sensor sites are the index finger for adults and paediatrics, the hand for infants and on the foot below the toes for neonates.

2. The finger is correctly inserted when

SC7500VM, SCM7500VM, SCP7500VM, SF7500VM:

- a. The tip of the finger touches the sensor's inside end (C).
- b. The sensor cable extends along the patient's hand (A).

NOTE: Ensure the sensor is not placed backwards (E), or with the finger extending (F).

W7500VM:

- a. The emitter and detector are positioned directly opposite each other.
- b. The sensor is applied firmly, but not too tight.

3. Connect the sensor to the monitor. Use an extension cable if needed (G), (H).
4. Switch on the oximeter and verify proper operation. (Ensure that the oximeter is providing a reading.)
5. Visually monitor the sensor site at least every 6 hours to ensure integrity of the skin and that no pressure injury is developing. Reposition to an alternative location if skin integrity changes.

Sensor removal:

1. Open the sensor by pressing on the sides of the sensor (SC7500VM, SCM7500VM, SCP7500VM) or at the front part of the housing (SF7500VM). Remove the W7500VM by opening the velcro fastener. Remove the sensor from the finger.
2. Disconnect the sensor from the monitor and follow the cleaning instructions, if necessary.
3. Store the sensor away for its next use.

Sensor expected life is 2 year.

WARNINGS

WARNING: Pulse oximeter probes and extension cables are designed for use with specific monitoring systems. Incompatible components can result in degraded performance and electromagnetic incompatibility, or patient injury. The operator is responsible for checking compatibility prior to use.

WARNING: A functional tester cannot be used to assess the accuracy of pulse oximeter probes.

WARNING: Check the patient skin condition at least every 6 hours and reposition to an alternative location if skin integrity changes.

WARNING: Misapplication of a pulse oximeter probe with excessive pressure (e.g. sensor too small or applied too tight) and for prolonged periods can cause pressure injury.

WARNING: Avoid application of the sensors to areas with poor skin integrity.

WARNING: Do not sterilize by irradiation, steam, or ethylene oxide – refer to cleaning and disinfection instructions. Use of agents other than specified may damage the sensor.

WARNING: Do not use the sensor or extension cable if it is damaged. Use of a damaged sensor or extension cable could cause patient injury or equipment failure.

WARNING: Excessive patient motion, excessive ambient light, electromagnetic interference (e.g. stacking and location to other medical equipment), dysfunctional haemoglobin, intravascular dyes, fingernail polish and long or artificial fingernails may affect the sensor performance and the accuracy of the measurement.

WARNING: Route all cables carefully to reduce the possibility of patient entanglement or strangulation.

WARNING: Do not alter or modify the sensor or extension cable. Alterations or modifications may affect performance or accuracy.

WARNING: Do not use the SMARTsat® sensor, or extension cable or other oximetry sensors during MRI scanning. Conducted current may cause burns.

CAUTIONS

CAUTION: Avoid possible interference from sources of electromagnetic interference such as, but not limited to: Cellular phones, radio transmitters, motors, telephones, lamps, electrosurgical units, defibrillators, and other devices. Interference may cause inaccurate measurements. If you are unsure if your device is working properly, contact your clinician.

CAUTION: Sensors or extension cables may work with different monitors. Refer to the specific monitor instructions for use for essential performance, accessories, compliance information (EMC) and electrosurgical applications.

Cleaning

It is necessary to clean and disinfect the sensor before attaching it to a new patient. Routine cleaning or low-level disinfection is recommended when using the sensor on the same patient.

Cleaning SC7500VM, SCM7500VM, SCP7500VM, W7500VM:

To clean the sensor with 70% Isopropyl alcohol:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Rinse the sensor and cable with distilled water for at least 2 minutes.
CAUTION: Do not immerse the connector of the sensor cable in liquid.
3. Clean the surface of the cable by wiping it with a cotton pad soaked in 70% Isopropyl alcohol.
4. Immerse the sensor housing in 70% Isopropyl alcohol for 5 minutes and let it air dry for at least 10 minutes.
5. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the cleaning process.

To disinfect (low-level) with 1:10 bleach dilution:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Rinse the sensor and cable with distilled water for at least 2 minutes.
CAUTION: Do not immerse the connector at the sensor cable in liquid.
3. Prepare 1:10 bleach dilution in distilled water. Wipe the cable with a cotton pad soaked with the solution of diluted bleach.
4. Immerse the sensor housing in the solution of diluted bleach for 5 minutes.
5. Rinse the sensor and cable by immersing them in distilled water for 5 minutes.
6. Allow 10 minutes for the sensor, cable, and connector to dry after rinsing.
7. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the disinfection process.

To disinfect (high-level) - CIDEX OPA™*:

1. Rinse sensor to remove surface dirt – the sensor housing may be fully immersed in liquid.
CAUTION: Do not immerse the connector at the sensor cable in liquid.
2. Clean the sensor and patient contact surfaces. Prolystica™* is recommended as a cleaning agent. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution and wipe the inner and outer surfaces with a soft brush or cloth to remove any visible soil. Rinse sensor with water following cleaning and wipe dry prior to disinfection. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the cleaning process.
3. CIDEX OPA™* is recommended for high-level disinfection. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution for 12 minutes and remove.
4. Thoroughly rinse the sensor with water following disinfection (5 minutes minimum distilled water rinse/soak).
5. Manually dry the sensor and verify the sensor, cable, and connector are dry prior to use.

Cleaning SF7500VM:

To clean the sensor with 70% Isopropyl alcohol:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Clean the surface of the cable and sensor head at the inside and outside by wiping it with a cotton pad soaked in 70% Isopropyl alcohol.
3. Allow 10 minutes for the sensor, cable, and connector to dry after cleaning.
4. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the cleaning process.

Specifications	LED Wavelength	Output Power ³
	Red: 660nm Peak	< 5 mW
	Infrared: 905nm Centroid	< 5 mW
Accuracy	SpO ₂ Accuracy: ⁴	Pulse Rate Accuracy: ⁵
	60 – 100%: $A_{rms} \leq 2\%$	20 – 300 BPM.: $A_{rms} < 2$ BPM. (refer to monitor's operating manual for more information) <i>BPM: Beats per minute</i>
	60 – 70%: $A_{rms} \leq 3\%$	
	70 – 80%: $A_{rms} \leq 2\%$	
	80 – 90%: $A_{rms} \leq 2\%$	
	90 – 100%: $A_{rms} \leq 2\%$	
	<60%: unspecified	
Water spray protection	SC7500VM, SCM7500VM, SCP7500VM, W7500VM	IP33: Protected from water spray less than 60 degrees from vertical
	SF7500VM	IP22: Protected from water drops less than 15 degrees from vertical

Note: Because pulse oximeter measurements are statistically distributed, only about two-thirds of measurements can be expected to fall within $\pm A_{rms}$ of the measured values by a co-oximeter. To verify the function of pulse oximeter probes a functional tester (Index 2 or equivalent) can be used.

³⁾ This information can be especially useful to clinicians

⁴⁾ Validated by clinical tests, where the measured values of the sensors were compared with those of the arterial co-oximetry in healthy adult subjects over the specified functional oxygen saturation range. These values were validated on the following monitor: VM-2160 / Oxytrue® A (SMARTsat® OEM III). A functional tester may not be used to validate SpO₂ accuracy.

⁵⁾ Pulse rate accuracy was verified by simulated bench tests to ensure entire range was verified. These values were verified on the following monitor: VM-2160 / Oxytrue® A (SMARTsat® OEM III).

Biocompatibility

The materials that come into contact with the patient have undergone extensive biocompatibility testing. Further information is available on request.

Environmental Specification

Operational temperature range:	<i>SC7500VM, SCM7500VM, SCP7500VM, W7500VM:</i> +0°C to +40°C (32°F to 104°F) <i>SF7500VM:</i> +5°C to +40°C (41°F to 104°F)
Transient operating conditions:	Sensor will operate for 20 minutes at -20°C and +50°C temperature conditions.
Storage, and transport temperature range:	<i>SC7500VM, SCM7500VM, SCP7500VM, W7500VM:</i> -40°C to +70°C (-40°F to 158°F) <i>SF7500VM:</i> -20°C to +70°C (-4°F to 158°F)
Relative humidity range:	15 - 95% non-condensing
Atmospheric pressure range:	620 hPa to 1060 hPa
Stabilization Time (from storage to operational conditions):	Up to 20 minutes

Note: These environmental specifications are specific to the sensor. For information on the environmental specification of the compatible monitoring systems, please see the Operator's Manual for that system.

Warranty

To obtain information about a warranty, if any, for this product, contact Viamed Ltd. (See back page).

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