VIAMED SMARTSAT® REUSABLE EAR PROBE PULSE OXIMETRY SENSORS

Intended Use: SMARTsat® Reusable Ear probe Pulse Oximetry Sensors, model EP7500VM, are intended for use with Viamed monitoring systems that use SMARTsat® technology.

The SMARTsat® Reusable Ear probe Pulse Oximetry Sensors are indicated for use in non-invasive continuous monitoring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate on adult and paediatric (excluding infant and neonatal) patients weighing greater than 20 kg.

The use environment may include: hospitals, hospital-type facilities, intra-hospital transport, mobile emergency medical applications and home use.

Transport environments include intra-hospital transport and ground emergency transport (road ambulances). The use environment is dependent upon the use environments specified for the monitoring system utilized in conjunction with the sensor.

Note: This Instruction for Use is intended for the medical professional. For use in the home care environment the user may use the device after instruction by, or under supervision of, a trained healthcare professional.

Caution: The SMARTsat® Reusable Ear probe Pulse Oximetry Sensors are for prescription use only.

Caution: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter.

Warnings

1. A functional tester (like Index II or equivalent) cannot be used to assess the accuracy of pulse oximeter probes.
2. Connect the sensor only to a SMARTsat® technology monitor or module. A mismatch may cause patient injury, inaccurate measurements or equipment damage.
3. Check the patient skin condition at least every 4 hours and reposition to an alternative location if skin integrity changes.
4. A misuse of a pulse oximeter probe with excessive pressure and with prolonged periods can cause pressure injury.
5. Refer to cleaning and disinfection instructions.
6. Avoid application of the sensors to oedematous or fragile tissue.
7. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure.
8. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional haemoglobin, low perfusion, intravascular dyes, may affect the sensor performance and the accuracy of the measurement.
9. Route all cables carefully to reduce the possibility of patient entanglement or strangulation.
10. Note: The materials used in the manufacture of the sensor contain no natural rubber latex protein. The materials have undergone extensive biocompatibility testing. Further information is available upon request.
11. Do not use this sensor with MRI scanning. Burns may be caused by conducted current.

Application of the Sensor

The use of the ear sensor is contraindicated for patients whose ears are pierced at the indicated site.

1. The preferred site is the patient's ear lobe, securing the sensor using the cable support. Alternatively, position the ear sensor on the ear pinna and secure using medical tape **(A)**.
2. Wipe the patient's ear clean with alcohol solution.
3. Insert the ear sensor cable into the cable support, leaving at least 10 cm residual cable to connect the sensor to the fleshy portion of the ear lobe **(B)**.
4. Hook the cable support over the patient's ear **(B)**.
5. Clip the sensor to the fleshy portion of the ear lobe, with the light source on the outside of the ear **(B)**.
If applying the sensor to the pinna, clip the sensor to the pinna with the light source on the outside of the ear, and ensure that the cable runs behind the ear and down the side of the patient's body.

6. Secure the sensor cable with medical tape.
7. Connect the sensor cable to the adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.

Cleaning or Disinfecting the Sensor

Clean or disinfect the sensor before attaching it to a new patient.

Disconnect the sensor from the monitor before cleaning or disinfecting.

Clean the sensor and patient contact surfaces with a soft cloth moistened in water.

To disinfect the sensor, wipe the sensor in 70% isopropyl alcohol.

Environmental Requirements

- Operating temperature range:
+5°C to +40°C (+41°F to +104°F)
- Storage and transportation temperature range
-30°C to +70°C (-22°F to +158°F)
- Storage relative humidity range:
10 to 95% non-condensing

SpO₂ Accuracy:

60 – 100%:	$A_{rms} < 2.5\%$
90 – 100%:	$A_{rms} < 2\%$
80 – 90%:	$A_{rms} < 2\%$
70 – 80%:	$A_{rms} < 3\%$
60 – 70%:	$A_{rms} < 4\%$
<60%:	Unspecified

Pulse Rate Accuracy:

20-300 BPM:	$A_{rms} < 2 \text{ BPM}$
BPM:	Beats per minute

Pulse oximeter equipment measurements are statistically distributed. Only about two-thirds of measurements can be expected to fall within ± 2 Arms of the measured values by a co-oximeter. To verify the function of pulse oximeter probes, a functional tester (like Index II or equivalent) can be used.

SpO₂ accuracy was validated by clinical tests, where the measured values of the sensors were compared with those of the arterial co-oximetry in healthy adult subjects over the specified functional oxygen saturation range. These values were validated on the following monitor: VM-2160 / Oxytrue® A (SMARTsat® OEM III). A functional tester may not be used to validate SpO₂ accuracy.

LED Specifications*

*This information is especially useful for clinicians

LED	Wavelength	Radiant Power (typical)
Red	660 nm (Peak)	2.5 - 3.5 mW
Infrared	905 nm (Centroid)	2.5 - 3.5 mW

Disposal

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed Ltd.

REF	Model	Viamed Type	Cable Length
6020132264	EP7500VM	0014850	1.20 m

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© 2018 bluepoint medical GmbH & Co. KG
SMARTsat® and OxyTrue® are registered trademarks of bluepoint medical GmbH & Co. KG.
bluepoint medical GmbH & Co. KG, An der Trave 15, 23923 Selmsdorf, Germany