
RD900 Neopuff[™] Infant Resuscitator



RD900 Neopuff[™]Infant Resuscitator

The RD900 Neopuff[™] Infant Resuscitator is an easy to use manually operated, gas powered resuscitator which provides controlled and accurate resuscitation of newborn babies in delivery suites, nurseries and neonatal intensive care units.

PRECISE CONTROL OF AIRWAY PRESSURE

- Once setup, the Neopuff[™] Infant Resuscitator will safely inflate the babys' lungs and provide optimum oxygenation without the risks associated with under or over inflation at uncontrolled pressures.
- The desired peak inspiratory and positive end expiratory pressures are set using simple fingertip controls.

FAST & ACCURATE MANOMETER

 Provides constant monitoring and reassurance of the delivered peak inspiratory pressure and positive end expiratory pressure.

LIGHTWEIGHT PATIENT T-PIECE

 Provides simple, effortless breath by breath resuscitation by occluding the patient T-piece aperture with your thumb or finger.









 The lightweight patient T-piece and patient supply line are especially useful over extended resuscitation periods.

INCREASED SAFETY

 The maximum pressure relief valve provides increased safety by allowing you to preset the maximum possible peak inspiratory pressure.

RD900 Neopuff[™] Infant Resuscitator



Accessories shown are optional. Please refer to product catalogue for part numbers.

Fisher Paykel

RD900 Neopuff[™]Infant Resuscitator

PRODUCT SPECIFICATIONS

Manometer range	-20 to 80 cm H ₂ O (-2.0 to 7.8 kPa)	
Maximum pressure relief	@ 10 LPM	3 to 80 cm H ₂ O (0.3 to 7.8 kPa)
Peak inspiratory pressure (PIP)	@ 10 LPM	2 to 80 cm H ₂ O (0.2 to 7.8 kPa)
Positive end	@ 5 LPM	1 to 5 cm H ₂ O
expiratory pressure	@ 10 LPM	2 to 15 cm H ₂ O
(PEEP)	@ 15 LPM	3 to 25 cm H ₂ O
Gas inlet	5 LPM (min) to 15 LPM (max)	
flow range	If the gas inlet flow rate increases from 5 to 15	
	LPM, the peak inspiratory pressure typically	
		proximately 8 cm H ₂ O (0.8 kPa).
Delivered oxygen concentration	Up to 100% depending on gas supply	
Operating time	@ 5 LPM	80 minutes
(400 L cylinder)	@ 10 LPM	40 minutes
	@ 15 LPM	26 minutes
		5

RD900 NEOPUFFTM INFANT RESUSCITATOR

Height	25 cm (9.8") 20 cm (7.9")	
Width		
Depth	10 cm (3.9°)	
Weight	2.1 kg (4.6 lb) excluding accessories Up to 10 kg (22 lb)	
Recommended body weight		
Operating and storage limits	90% relative humidity	

For full product specifications please contact your Fisher & Paykel supplier.

Certification



ISO 8382

Type Certified No. U0111

CITECH

MODEL SHOWN

RD900 Neopuff™ Infant Resuscitator.

ALSO AVAILABLE

Installed in the CosyCot™ Infant Warmer.

ACCESSORIES SHOWN

The accessories shown are intended only to illustrate a typical configuration for resuscitation. The 900 series Neopuff™ Infant Resuscitator and accessories can be customised to suit your requirements.

OTHER MODELS & ACCESSORIES

For more information please refer to our product catalogue.

ORDERING INFORMATION

The 900 series Neopuff™ Infant Resuscitator and accessories are available separately or as a complete system. Please refer to our product catalogue or contact your local Fisher & Paykel supplier for further information and assistance.

Fisher & Paykel have a policy of continuous product improvement and reserve the right to alter specifications without notice.







INTERNATIONAL AUSTRALIA UK / IRELAND* US FRANCE / BENELUX 25 Carbine Road, PO Box 14-348, Panmure, Auckland, New Zealand. Tel: 64 9 574 0100. Fax: 64 9 574 0158. Email: warmers@tphcare.tp.co.nz 36-40 New Street, Ringwood, Victoria 3134, Australia. Tel: 61 3 9879 5022. Fax: 61 3 9879 5232

Unit 16 Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK. Tel: 44 1628 626 136. Fax. 44 1628 626 146. *EU authorised representative. 22982 Alcalde Drive, Suite 101, Laguna Hills CA 92653, USA. Tel: 1 949 470 3900 or 1 800 446 3908. Fax: 1 949 470 3933 10 Avenue du Quebec, SILIC S12-Villebon, 91946 COURTABOEUF CEDEX, France. Tel: 33 164 465 201. Fax: 33 164 465 221