Final Q.A. of TT460.

Attach O2 hose to flowmeter on TT460.

Set flowrate to 10 Lpm.

Read Flow through Flowmeter CE079

Centre of ball on CE0xx should be >9 and < 11

Gentlt shake Flowmeter the ball should bounce then return to 10 Lpm

Attach tubing of Test Guage (CE078)

Meter should read 45 Cm WG +/- 5 on test gauge.

remove and reattach tubing several times readings should be consistant.

Use Snoop (Liquid leak detector) on O2 hose and right angle to cheack for leaks.

Check Flowmeter tube for scraches
Check All Lables are attached
CE Lable
S/N Lable
O2 Flow Lable

Fill in and sign Q.A. Form:

Final Q.A. TT495.



Attach O2 hose to the O2 Flowmeter on TT495 and set flow to 10 Lpm. attach tubing of Test Flowmeter. The centre of the ball on the Test Flowmeter (**CE079**) should be between 9 and 11.

Gently shake the TT495 so the Flowmeter Ball bounces it should return to the 10Lpm mark

Attach Air hose to the Air Flowmeter on TT495 and set flow to 10 Lpm. attach tubing of Test Flowmeter. The centre of the ball on the Test Flowmeter should be between 9 and 11.

Gently shake the TT495 so the Flowmeter Ball bounces it should return to the 10Lpm mark

Attach output of TT495 to Test Gauge (CE078 Manometer) tubing using 15mm Adaptor Flowrate set at 10Lpm Set adjustable valve to lowest setting.

Test Gauge should read 6 or less While attached to the Test Guage increase the adjustable valve to the Maximum setting

Test Gauge should move up the scale to >42 and < 50(+/- 5 CmWg)

Block all holes in adjustable valve. Using thumb and finger.

Safety valve should exhaust at >45 and <50 (+/- 5 CmWg)

Set adjustable valve to read 10cm H2O on the pressure gauge and note the readings obtained from the Manometer.

Check Flowmeter tubes for scratches
Check All Labels are attached
CE Label
S/N Label
Thumb Vent Label
Reduce pressure Label

Fill in and sign Q.A. form: