DIRECTIONS FOR USE


Intended Use: Viamed 4000 Series oximetry sensors are indicated for use in continuous noninvasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters and patient monitors. Select the sensor model number that matches the brand of oximeter or patient monitor that you are using.

Precautions: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. Do not autoclave or immerse the sensor in liquid. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings: Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor every 8 hours and reposition to an alternate location every 24 hours. Avoid application of the sensors to edematous or fragile tissue. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. Excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional hemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

Application of the Sensor

1. Connect the sensor cable to the Viamed adaptor cable that fits your oximeter or directly to the patient monitor. Then switch on the oximeter and verify proper operation.
2. Insert the patient's digit into the sensor. The preferred digit to use for adults is the index finger. The digit is correctly inserted when:
 - The tip of the digit touches the rear guide posts.
 - The sensor cable extends along the top of the patient's hand.
3. Visually monitor the sensor site to ensure the integrity of the skin.
4. Secure the lead with medical tape.

Cleaning the Sensor

To clean the sensor between uses or between patients, wipe the sensor surfaces with alcohol.

Environmental Requirements:

- Operating temperature range:
10° - 40° C (50°-104° F)
- Storage and transportation temperature range
0° - 50° C (32°-122° F)
- Relative humidity:
0 - 75% non-condensing

SpO₂ Accuracy:

70 - 100%: A_{rms} <3.0%
0 - 69%: Unspecified

Accuracy test methods available from the manufacturer upon request.

LED Specifications

	Wave Length	Radiant Power (typical)	
Red	660nm	3.5 mW	
Infrared	885nm	3.5 mW	Nellcor, CSI
	905nm	3.5 mW	Datex, BCI, Viamed
	910nm	3.5 mW	Nonin
	940nm	3.5 mW	Ohmeda, Novamatrix

Disposal:

Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed.

All registered trademarks recognized.

Viamed Model	Oximeter Brand
SC4001	Nellcor
SC4002	Nellcor
SC4006	HP/Phillips
SC4007	HP/Phillips
SC4011	Datex
SC4012	Datex
SC4016	Ohmeda
SC4017	Ohmeda
SC4021	BCI
SC4026	Nonin
SC4036	CSI
SC4037	CSI
SC6500 VM	Viamed

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Viamed Ltd