Manufactured for
VIAMED LTD – United Kingdom

Flow sensor

Type: APC, 4310002

REF 7030131001

EN Instructions for use

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bluepoint MEDICAL GmbH & CO. KG
An der Trave 15
23923 Selmsdorf, Germany
www.bluepoint-medical.com

CE 0482

EN Instructions for use

Flow sensor APC, 4310002

Intended Use

Hot wire anemometric flow sensor, for measuring gas flow in Dräger anaesthesia and ventilation equipment. The sensor is compatible to:

- **Spirolog - REF 8403735-06**
- **Spirolife - REF MK01900**

Operation

- Due to fire risk - do not introduce any medicaments or other substances based on combustible solvents, such as alcohol, into the patient systems.
- Do not use explosive anaesthetic agents (e.g. ether or cyclopropane).
- After using with medicament nebulisers, immediately clean and treat the flow sensor as specified. Alternatively, install a filter upstream of the flow sensor.

Care

- Alignment, operation and sensor change must be carried out in accordance with the relevant device-specific instruction manual.
- Do not clean in automated cleaning or disinfection equipment. Do not clean with compressed air, water jets, brushes etc.
- Reusing the flow sensor may affect patient safety due to the possible contamination of the sensor.
- **Life cycle of the sensor is limited. It can be used as long as calibration in the anaesthesia or ventilation equipment is possible.**

Cleaning and Disinfection

The following cleaning agents and disinfectants can be used:

- Korsorex Basic/Extra/Plus (Bode Chemie, Hamburg, Germany)
- Gigasept FF/AF (Schülke&Mayr, Norderstedt, Germany)

Strictly follow the manufacturer's instructions on the disinfectant.

- Disinfect the flow sensor by immersion.
- Inspect the flow sensor for soiling and damage.
- Rinse the flow sensor under water (preferably sterile water) and then allow to dry in air for at least 30 minutes.

Sterilization

The flow sensor can be autoclaved (134°C).

- Use a steam sterilizer with fractional vacuum processes (e.g. Webeco A30-B). Follow the instructions for use of the sterilizer.
- Pay attention to the water quality for steam generation.
- The flow sensor must be cleaned and disinfected before sterilization.
- The flow sensor must not be covered with other items to be sterilized.
- The flow sensor must not touch the chamber walls of the sterilizer.
- Steam sterilization parameters:
134 °C with pre-vacuum and an exposure time of 3.5 minutes

Disposal

The flow sensor must be disposed of as infectious waste.

Measuring range: 0 to 200 L/min, depending on the gas type and operating conditions

Accuracy: DIN EN ISO 80601-2-12,
DIN EN ISO 80601-2-13
Typically ±8% of measured value under calibration conditions depending on Ventilator(evaluation electronics)

Recommended operating conditions

- Temperature: 10°C to 40°C
- Relative humidity: 5% to 95%, non-condensing
- Atmospheric pressure: 700hPa to 1060hPa

Recommended storage and transport conditions

- Temperature: -20°C to 50°C
- Relative humidity: 5% to 95%, non-condensing
- Atmospheric pressure: 700hPa to 1060hPa

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