Medical Oxygen Sensor / Analyzer Warning Information

Teledyne medical oxygen sensing products are of the highest quality and are designed, manufactured and tested to provide reliable and accurate service when used properly and in accordance with the manufacturer's instructions.

However, certain risks will always be present as is the case with many products used in the medical arena. The role of this document is to help, you, the user minimize such risks.

The risk of improper output:

- The analyzer and/or sensor uses electrochemical technology and may be subject to drift over time. Accurate and frequent calibration of the sensor/analyzer with a known reference gas or gas mixture is essential. 100% oxygen is the preferred and recommended calibration gas followed by check in ambient air (20.9% oxygen). If ambient air is used as a calibration gas, only oxygen concentration up to 50% can be accurately measured. Also, make sure that the air has sufficient circulation such that the concentration in the room is certain to be 20.9%. Typically sensors and analyzers are capable of achieving accuracy in excess of 2% oxygen.
- If the output of the sensor is erratic and/or unstable when exposed to a gas mixture do not attempt to use the sensor or analyzer. The sensor/ analyzer should be checked out by a qualified technician.
- Do not attempt to use the sensor/analyzer outside the maximum ambient conditions (temperature and pressure) specified by the manufacturer. Always allow the sensor/analyzer to come to equilibrium with the local environment before use.
- Do not allow moisture to accumulate on the sensor's sensing surface. When used under high humidity conditions, ensure that the sensor face points downward and that the sensor is not at the lowest point in the circuit.
- Sensors are pressure sensitive and respond to pressure changes. Many sensors are supplied with a shorting device which is used to lessen the initial stabilizing period when the sensor is first installed. If for some reason the shorting device is not effective, the stabilizing time may be longer than normal.

The risk of contamination:

- The oxygen sensor as well as the batteries which power the analyzer contain a caustic solution. If any leakage from the sensor or batteries appears to be present, handle with extreme care and do not attempt to use the product. Special care should be taken so as not to puncture the sensing membrane.
- Please follow the correct disposal procedures in the disposing of sensors and batteries; consult the relevant Material safety Data Sheet (MSDS) for disposal and electrolyte exposure information.
- Sensors/analyzers are not shipped in a sterile condition.
- Sterilize appropriate items when exposed to infections or contaminated persons, environments or surfaces. Follow the sterilization procedure listed in the user's manual.

The risk due to mishandling:

- The sensor/analyzer should be handled with care at all times to prevent damage, possible inaccurate measurements, or electrolyte leakage.

If further information is needed contact Teledyne Analytical Instruments (TAI) or your nearest TAI distributor.

TELEDYNE
ELECTRONIC TECHNOLOGIES
Analytical Instruments
16630 CHESTNUT STREET, CITY OF INDUSTRY, CA 91748-1020
P.O.BOX 1580, CITY OF INDUSTRY, CA 91749-1580
TELEFONO: (626) 934-1500, (626) 961-9021
FAX: (626) 961-2138

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European Master Distributor
Viamed, Ltd.
Kiehl, Reino Unido

DECLARATION OF CONFORMITY

APPLICATION OF COUNCIL DIRECTIVE : 89/336/EEC

STANDARDS TO WHICH
CONFORMITY IS DECLARED : EN60601-1-2

MANUFACTURER'S NAME : Teledyne Brown

MANUFACTURER'S ADDRESS : 16830 Chestnut Street
City of Industry, Ca. 91748
USA
(818) 961-9221

TYPE OF EQUIPMENT : Oxygen Analyzer-Medical
Equipment

EQUIPMENT CLASS : Medical Electrical
Equipment

MODEL NUMBER : TED60

I THE UNDERSIGNED, HEREBY DECLARE THAT THE EQUIPMENT
SPECIFIED ABOVE, CONFORMS TO THE ABOVE DIRECTIVE(S)
AND STANDARD(S).

PLACE: City of Industry, CA

SIGNATURE: [Signature]

FULL NAME: STEPHEN BROY

POSITION: ENGINEERING MANAGER