

14 Russell Square, London, WC1B 5EP Telex 883669 DHSSHO G Fax 01-637 8990 3103/3317 Telephone 01-636 6811 Ext.

Your reference:

Our reference: STD/PG4e SCB/NHR

Date:

10 June 1988

#### LIMITED HAZARD NOTICE

#### VIAMED - PHOTOTHERAPY LIGHT SHIELD

### Summary

The Department has received a report that these tinted perspex shields which are used with headboxes during phototherapy may not be effective in protecting the eyes of the babies from the ultra violet radiation and bright light emitted by the phototherapy equipment.

#### 1. Background

The traditional practice during infant phototherapy is to cover the eyes of babies to protect them from the light and ultra violet radiation.

- As an alternative to this Viamed introduced in 1984 a tinted perspex sheet, shaped in the form of an inverted U which fits over a standard incubator head box and protects the babies eyes from radiation.
- The Department has some limited and so far uncorroborated evidence that these shields may not be effective in stopping harmful radiation.
- This HN is published as an early warning, and as more information becomes available, supplemented by measurements commissioned by the Department, a further publication will be issued.
- As the location of all these Phototherapy Light Shields is known with certainty, it is possible to send a Limited Hazard Notice directly to known users.

#### Action

- Users of Perspex Phototherapy Light Shields, supplied by Viamed, 15 Station Road, Cross Hills, Keighley, West Yorkshire and shaped in the form of an inverted U, some 9" x 8" x 6" are requested to withdraw these from use immediately, and revert to the traditional method of eye protection.
- 7. Further information will be sent to you in due course.
- Enquiries to the Department may be addressed to Mr S C Bishop or Mr N H Richardson at the address given at the top of this Notice.

35900/MV/1

# SAFETY ACTION BULLETIN

No 45

DEPARTMENT OF HEALTH
SCOTTISH HOME AND HEALTH DEPARTMENT
WELSH OFFICE
DEPARTMENT OF HEALTH AND SOCIAL SERVICES
(NORTHERN IRELAND)

SAB(89)19 March1989

#### PHOTOTHERAPY LIGHT SHIELDS

#### SUMMARY

An effective material called Amber 300 has been identified which overcomes the problems discussed in the Limited Hazard Notice issued in June 1988 to the 55 hospitals having defective light shields. At that time the Department recommended the withdrawal from use of these dark grey tinted headbox shields because they were suspected of transmitting more ultra violet and blue light during phototherapy than was safe. Headboxes constructed of Amber 300 are now available.

#### ACTION

This information should be brought to the attention of all who need to know or be aware of it. This will include safety officers, safety representatives, supplies officers, nursing and midwifery staff in neonatal intensive care units, SCBU's and maternity departments.

It is recommended that tinted headbox shields made only from Amber 300 should be used in conjunction with headboxes as eye protection during phototherapy.

Plain perspex headboxes and coloured perspex phototherapy shields made of Amber 300 can be obtained from either of the following:

Bio-Med Engineering Ltd Hendlands Road Liversedge Yorks WF15 7NT (Tel: 0924 407811) Viamed 15 Station Road Cross Hills Keightley West Yorks BD20 7DT (Tel: 0535 34542).

#### BACKGROUND

- 1. The result of light transmission tests of the origional dark grey tinted headbox shields by the Department of Clinical Physics, Glasgow and the National Physical Laboratory confirmed the excessive transmission of UV and blue light through the material.
- 2. Following discussions with ICI a perspex material called Amber 300 was identified and tested. It was evident from the results that a shield made from this coloured material would prevent any transmission of UV or blue light over the test range of 300 to greater than 500 nm. This includes the range of radiation considered to be injurous to the infants' retinas.
- 3. Being of a bright colour these shields enable a rather better view of the baby's head than the originnal dark grey tinted shields which of course should not be used.



Our reference: Your reference:

Door Ar Lamb.

13 Sept 85

Herewith to got top which we discussed this morning

Jon kindly offered to carry out your proposed muelification on this one as a Rample prior to madifying the remainder in service at various Dospitals.

Hank you for your cooperation in the matter and I lope we can quickly bring it to a catefactory conclusion.

your Sincorely

C.J. GADD STBID.

## **ENQUIRIES**

4. 7	rechnical enq	uiries regarding	his bulletin	should be	addressed	as follows:
------	---------------	------------------	--------------	-----------	-----------	-------------

England

Mr S C Bishop or Mr N H Richardson Department of Health NIIS Procurement Directorate 14 Russell Square London WC1B 5EP (Tel: 01 636 6811 Ext 3103/3317

for health authority use as instructed locally.						
TITLE: HUMAN DURA AND CREUTZFELDT-JAKOB DISEASE	Hazard(89)11 17 March 1989					
Sent to department responsible for action:	Date					
Action completed: SignatureLocation	Date					
Comments on action taken:						
Return to after completed action:	Date					