

## VM3COP64.03 Aii-2000PALM O2 Calibration Test

Before using this calibration test procedure, the user should refer to the appropriate Aii-2000PALM O2 user manual (Revision 10).


Required equipment: Time Electronics LTD Microcal 1030 – Calibrated Voltage source (CE076) - Calibration traceable to UKAS calibration standards. The specifications can be found on Document ID 30969.

Test lead – Molex to 4 mm banana plugs.


- 1) Remove the 4 case screws on the back of the Aii-2000PALM.



- 2) Open the case and disconnect the Oxygen sensor from the PCB.



- 3) Connect the Microcal to the Aii-2000PALM O2R analyser using the test lead as shown.



- 4) Set the Microcal to 100 mV scale and the polarity switch to NORM.



- 5) Set the output adjustment dial to 35.0mV.
- 6) Power up the Aii-2000PALM O2 analyser.
- 7) Press the CAL button on the front panel to calibrate the analyser.



- 8) At the end of the calibration procedure the analyser should display 100.0. If the Aii-2000PALM O2 analyser fails to calibrate, check the connections and repeat the test.

The Microcal is simulating the output of an Oxygen sensor reading a 100% Oxygen concentration. By decreasing the output of the Microcal to set values the accuracy of the Aii-2000PALM O2 analyser can be determined across a range of simulated Oxygen concentrations.





- 9) In turn, dial in the input voltages from the table below. Record the values on the certificate (0180030).
- 10) Once the calibration test has been performed, disconnect the test lead and reconnect the Oxygen sensor.
- 11) Carefully align the upper and lower case and replace the case screws.
- 12) Save the document and upload to the Intrastats calibration index.
- 13) Print a copy on Viamed letter headed paper.

Microcal output	Aii-2000PALM O2 expected reading ( $\leq 2\%$ relative)
35.0 mV	100% (98.0% - 102.0%)
31.5 mV	90% (88.2% - 91.8%)
28.0 mV	80% (78.4% - 81.6%)
24.5 mV	70% (68.6% - 71.4%)
21.0 mV	60% (58.8% - 61.2%)
17.5 mV	50% (49.0% - 51.0%)
14.0 mV	40% (39.2% - 40.8%)
10.5 mV	30% (29.4% - 30.6%)
7.3 mV	20.9% (20.5% - 21.3%)
6.3 mV	18% (17.6% - 18.4%)