
**ARISTO**<sup>TM</sup>  
ARISTOMEDICAL PRODUCTS, INC.


## Aristo's Comp,Sat Neonatal Disposable Sensors

LOT: 0038-1


### Instructions for Use

1. Select a suitable site for the sensor. The ball of the foot is the preferred site. The heel of the hand may be used as an alternate site.

2. Remove the "Pull Here" adhesive liner and position the detector on the sole of the foot. Align the light source on the top of the foot *directly opposite* the detector. Wrap tape around foot. **Do not stretch or pull tape during application.**

3. Route the cable along the top side of the foot and secure cable with tape if necessary.
- 3a. If using heel of hand, route the cable along side the arm and secure cable with tape if necessary.

4. Connect the sensor cable into the patient cable and verify proper operation as described in the instrument's operator manual.

### Aristo Comp,Sat Neonatal Disposable

|   | Manufacturer            | Cat#  |
|---|-------------------------|-------|
|   | BCI <sup>®</sup>        | 221-4 |
|   | CSI <sup>®</sup>        | 281-4 |
|   | Datex <sup>®</sup>      | 251-4 |
|   | Nellcor <sup>®</sup>    | 311-4 |
|   | NK <sup>®</sup>         | 201-4 |
| ✓ | Ohmeda <sup>®</sup>     | 241-4 |
|   | Novamatrix <sup>®</sup> | 261-4 |
|   | Aristo <sup>®</sup>     | 606   |

International distributor:

**ARISTO**  
ARISTOMEDICAL PRODUCTS, INC.

Manufacturer:

Opto Sensors (M) Sdn. Bhd.  
No 8, Jalan Firma2/2  
Tebrau Industrial Estate 1  
81100 Johor Bahru, Malaysia  
Tel : (607) 3537008  
Fax : (607) 3537010

EC Representative:

Rapiscan Security Products Ltd  
Unit B1, The Fleming Centre  
Fleming Way, Crawley  
West Sussex  
England RH102NN  
Tel : 44-1293-540601  
Fax : 44-1293-516945

CE  
0123

## Indications

The Comp./Sat™ Neonatal Disposable Sensor is indicated for use when continuous or spot monitoring of non-invasive arterial oxygen saturation is required for patients weighing less than 3 kg. The sensor is indicated for single patient use only.

## Contraindications

The Comp./Sat™ Disposable Sensor is contraindicated for use on patients who exhibit allergic reaction to the adhesive tape.

NOTE: The materials used in the manufacture of the sensor contain **NO NATURAL LATEX PROTEIN**. The sensor's materials have undergone extensive biocompatibility testing. Further information is available upon request.

## Warnings and Precautions

1. Federal law (US) restricts the use or sale of this device by or on the order of a physician.
2. The Comp./Sat™ Disposable Sensor is intended for use only with the designated oximeter monitors or modules. Use of this sensor with instruments other than the designated compatible monitors or modules may result in improper performance.
3. This product should not be used in the presence of flammable anesthetics.
4. Operations may be affected in the presence of strong electromagnetic sources such as electrosurgery equipment.
5. The product is not recommended for use in the presence of imaging equipment such as magnetic resonance imaging (MRI) and computed tomography (CT) devices.
6. Do not autoclave or immerse in liquid of any kind including cleaning solutions. Do not attempt to sterilize.
7. If the sensor is damaged in any way, discontinue use immediately.
8. Significant levels of dysfunctional hemoglobin such as carboxyhemoglobin or methemoglobin will affect the accuracy of the oxygen measurement.
9. Operation may be affected in the presence of high ambient light. Shield the sensor area with a surgical towel if necessary.
10. Check the sensor every 4 to 8 hours and reposition to an alternate location within a maximum of 24 hours. Application of the sensor to edematous or fragile tissue must be avoided.
11. Failure to apply the sensor or the replacement tape properly may cause incorrect measurement.
12. Intravascular dyes may interfere with performance of the sensor and the monitor which may cause inaccurate measurements.
13. The performance of the sensor and the monitor may be compromised by excessive motion.
14. A long or artificial fingernail may interfere with proper optical alignment. Also avoid finger nails with protective or colored polish. Relocate the sensor to another recommended site.
15. Poor plethysmograph wave forms, or inaccurate saturation measurements may be caused when the sensor is applied too tightly. This should be avoided.
16. For optimal performance keep the sensor site at the level of the heart.
17. Placing the sensor on any extremity with an arterial catheter, blood pressure cuff, or venous infusion line should be avoided.
18. High oxygen levels may present risk of retinopathy in the premature infant. In such cases, set the upper limit of SpO<sub>2</sub> in accordance with the acceptable clinical standards.

## SpO<sub>2</sub> Accuracy

70%-100% SaO<sub>2</sub> +/- 3 digits  
50%-70% SaO<sub>2</sub> unspecified

## Environmental Requirements

Operating temperature range: 50-104 F (10-40 C)  
Storage temperature range: 32-122 F (0-50 C)  
Relative humidity: 0-75% non-condensing

US Patients 5,217,012

Follow local greening ordinances and recycling plans regarding disposal or recycling of device components. A list of component material is available by contacting the manufacturer.

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