SpiroTrue® APC

Flow sensor

REF 4310002

EN

Instructions for use

DE

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IT

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bluepoint MEDICAL GmbH & Co KG
An der Trave 15
23923 Selmsdorf, Germany
www.bluepoint-medical.com

CE 0482

EN Instructions for use

SpiroTrue® APC – Flow sensor

Intended Use

Hot wire anemometric flow sensor, for measuring gas flow in Dräger anaesthesia and ventilation equipment. The sensor is compatible to:

- Spirolog - REF 8403735-06
- Spirolife - REF MK01900

Operation

- Due to fire risk - do not introduce any medicaments or other substances based on combustible solvents, such as alcohol, into the patient systems.
- Do not use explosive anaesthetic agents (e.g. ether or cyclopropane).
- After using with medicament nebulisers, immediately clean and treat the flow sensor as specified. Alternatively, install a filter upstream of the flow sensor.

Care

- Alignment, operation and sensor change must be carried out in accordance with the relevant device-specific instruction manual.
- Do not clean in automated cleaning or disinfection brushes etc.
- Reusing the flow sensor may affect patient safety due to the possible contamination of the sensor.
- Life cycle of the sensor is limited. It can be used as long as calibration in the anaesthesia or ventilation equipment is possible.

Disinfection

The following instrument disinfectants can be used:

- Korsorex Basic/Extra/Plus (Bode Chemie, Hamburg, Germany)
- Gigasept FF/AF (Schülke&Mayr, Norderstedt, Germany)
- Incidin Liquid (Ecolab Deutschland GmbH, Germany)

Strictly follow the manufacturer's instructions on the disinfectant.

- Disinfect the flow sensor by immersion.
- Inspect the flow sensor for soiling and damage.
- Rinse the flow sensor under water (preferably sterile water) and then allow to dry in air for at least 30 minutes.

Sterilization

- The flow sensor can be autoclaved (134°C).

Disposal

The flow sensor must be disposed of as infectious waste.

Recommended storage and transport conditions

- Temperature: -20°C to 50°C
- Relative humidity: 5% to 95%, non-condensing
- Atmospheric pressure: 700hPa to 1060hPa

Recommended operating conditions

- Temperature: 10°C to 40°C
- Relative humidity: 5% to 95%, non-condensing
- Atmospheric pressure: 700hPa to 1060hPa

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VIAMED LTD - UK