

COMPANY OPERATING PROCEDURES

Teledyne R21A (PN 0110121),
O₂ Sensor Test Procedure

VM3/COP/37.36

Date: 3-Feb-03

Revision date: 12-Jul-18

Issue: 3

Test & Equipment Required

Functional test required using R21 tester and R21 test adapter, (digital multimeter and adapter).

Output specification as in document: (VM3COP37.00, column (i)).

Method.

1. Open cell packaging and remove O₂ sensor.
2. Check the O₂ sensor for damage and signs of leaking electrolyte.
3. Connect the O₂ sensor to the R21 tester, select 200mV range (marked with arrows) and wait a few seconds for meter reading to stabilise.
4. Check for an output Ensure that meter reading lies within the output specification range,

7mV - 14mV 9-14 mV
* Scan into Intrastats as passed with output on Q/A Results
5. If no, or very low output, or sensor is returned as faulty leave the O₂ sensor exposed to air for 1 Hour Minimum.

6. If the O₂ sensor fails the output specification test, it should be left unsealed and boxed, placed in a ~~docket~~ ^{and} docket and put on the goods in desk with a short note explaining why it failed, ready to be booked in with a SRN ~~the~~

ON Q/A Bench ~~test~~ labeled Return to Supplier
Each sensor should have its barcode attached with the output/ problem on it.
7. Disconnect the O₂ sensor from the R21 tester and label with Viamed stickers if appropriate. ^{Box and add Viamed outer box labels.}
8. Reseal in serial numbered packet, integrity seal the edge of the packet with initialed and dated Viamed sticker and re-box.

Also needs to be noted as a fail with output on Intrastat system on Q/A Results.

cg 23.8.18