
SMARTsat®


Disposable SpO₂ Sensors

de	Gebruiksaanweisung	4
en	Instructions For Use	9
fr	Mode d'emploi	14
bg	Инструкции за употреба	19
cs	Návod k použití	25
da	Brugsanvisning	30
el	Οδηγίες χρήσης	35
es	Instrucciones de uso	40
et	Kasutusjuhend	45
fi	Käyttöohjeet	50
hr	Upute za upotrebu	55
hu	Használati útmutató	60
it	Istruzioni per l'uso	65
ja	取扱説明書	70
ko	사용 지침	75
lt	Naudojimo instrukcija	79
lv	Lietošanas instrukcija	84
nl	Gebruiksaanwijzing	89
no	Bruksanvisning	94
pl	Instrukcja użycia	99
pt	Instruções de utilização	104
pt-br	Instruções de uso	109
ro	Instrucţiuni de utilizare	114
ru	Инструкция по применению	119
sk	Návod na použitie	125
sl	Navodila za uporabo	130
sr	Uputstvo za upotrebu	135
sv	Bruksanvisning	140
tr	Kullanım Talimatları	145
zh	使用说明	150


AP


PP


IP


NP


Y7500


(A) AP / PP / IP / NP


(B) AP / PP


(C) IP




(D) NP

(E) Y7500 + HGS7500

(F) Y7500 + HGP7500**(G)****(H)**

XT6500/ XT6501 → ← AP / PP / IP / NP / Y7500

#	REF
10-AP  24	6020131194
10-PP  24	6020131197
10-IP  24	6020131199
10-NP  24	6020131201
XT6500	1020132275
XT6501	1020122058
Y7500  24	12020132001
HGP7500  48	12020121002
HGS7500  24	12020121001

#	KX
10-AP / 10-PP	0.45m / 1.5ft
10-IP / 10-NP	0.9m / 3ft
Y7500	1.2m / 4ft
XT6500	1.2m / 4ft
XT6501	2.4m / 8ft

SMARTsat®**Pulsoximetrie-Einweg-Sensoren**

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Disposable Pulse Oximetry Sensors

1 Directions for use

Note: This Instruction for Use is intended for the medical professional. Before use, carefully read the pulse oximetry sensor Instructions for Use, including all warnings, cautions, and instructions.

1.1 Intended purpose

Pulse oximeters including accessories provide measurements of arterial blood oxygen saturation (SpO₂) and pulse rate.

1.2 Indications for use

The SMARTsat® Disposable Pulse Oximetry Sensors are intended for use in non-invasive continuous and/ or spot-check monitoring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate on adult, paediatric, infant and neonatal patients, depending on the sensor style used.

The use environment may include: hospitals, hospital-type facilities, intra-hospital transport and mobile emergency medical applications.

Transport environments include intra-hospital transport and ground emergency transport (road ambulances). The use environment is dependent upon the use environments specified for the compatible monitoring system in combination with the sensor used.

The SMARTsat® Disposable Pulse Oximetry Sensors are intended for use with monitoring systems with integrated SMARTsat® technology. For monitor/ sensor compatibility, refer to the monitors instructions for use. The SMARTsat® Disposable Pulse Oximetry Sensors are for prescription use only.

1.3 Contraindications

There are no known contraindications to the use of the sensors.

1.4 Warnings and cautions

- Pulse oximeter probes and extension cables are designed for use with specific monitoring systems. Incompatible components can result in degraded performance and electromagnetic incompatibility, or patient injury. The operator is responsible for checking compatibility prior to use.
- A functional tester cannot be used to assess the accuracy of pulse oximeter probes.
- To prevent skin erosion and pressure necrosis, check the patient skin condition at least every 4 hours and reposition to an alternative location if skin integrity changes. Check the measurement site more frequently with patients which may have low tissue perfusion or edema at the measurement site.
- Do not apply the sensors to areas with poor skin integrity to prevent damage to the skin.
- Do not reuse disposable sensors as this can cause cross infection, risk of electric shock and inaccurate readings.
- Do not strap the sensor tape, wrap or the medical tape too tight to prevent blood flow restriction, inaccurate readings and pressure injuries to the skin beneath the sensor.
- During measurement, the sensor generates red and infrared light with specific fixed wavelengths as listed in the specifications below. Consider that these wavelengths might influence diagnostic parameters of other optical applications.
- Keep the SpO₂ sensor out of the radiation field during computed tomography or full-body radiation to avoid possible inaccurate readings for the duration of the active irradiation period.
- Do not sterilize, clean, disinfect or immerse in liquid of any kind to prevent product damage which may result in patient injury.

- Do not use the sensor or extension cable if it is damaged. Use of a damaged sensor or extension cable could cause patient injury or equipment failure.
- Excessive patient motion, excessive ambient light, electromagnetic interference (e.g. stacking and location to other medical equipment), blood pressure cuff on sensor site, venous congestion, abnormal venous pulsations, abnormal pulse rhythms due to physiological conditions or induced through external factors, dysfunctional haemoglobin, intravascular dyes, fingernail polish and artificial fingernails may affect the sensor performance and the accuracy of the measurement.
- Route all cables carefully to reduce the possibility of patient entanglement or strangulation.
- Do not alter or modify the sensor or extension cable. Alterations or modifications may affect performance or accuracy.
- Properly apply sensors according to sensor's directions for use. If the sensor is not applied properly or partially gets loose, this may cause no or inaccurate readings.
- Select a well perfused site for monitoring, very low perfusion at the monitored site may result in no or inaccurate readings.
- Remove the SMARTsat[®] sensor, extension cable or other oximetry sensors from the patient during magnetic resonance imaging MRI. Conducted current may cause burns.
- If SpO₂ values indicate hypoxemia, a laboratory blood sample should be taken to confirm the patient's condition.
- Avoid possible interference from sources of electromagnetic interference such as, but not limited to: Cellular phones, radio transmitters, motors, telephones, lamps, electrosurgical units, defibrillators, and other devices. Interference may cause inaccurate measurements. If you are unsure if your device is working properly, contact your clinician.

2 Instructions for use

2.1 Sensor selection

Select the sensor and application site on the patient. The preferred sensor sites are the index finger for adults and paediatrics **(B)**, the big toe for infants **(C)** and on the foot below the toes for neonates **(D)**.

Sensor	AP	PP	IP	NP	Y7500 *
Intended patient	Adult	Paediatric	Infant	Neonatal	All
	> 30 kg	10 to 50 kg	3 to 20 kg	< 3kg (foot/ Hand) > 30kg (Finger/Toe)	> 1kg
Application site	Finger or toe	Finger	Thumb or great toe	Foot or hand Thumb or great toe	Neonatal: hand or foot Adult: Finger or toe

* Not available in all countries

2.2 Sensor application

2.2.1 Adhesive sensors

1. Remove the protective film from the adhesive surface of the sensor **(A)**.
2. Apply the sensor on the patient: AP and PP **(B)**; IP **(C)**; NP **(D)**
 - Position the emitter and detector directly opposite each other
 - Apply the sensor firmly, but not too tight
3. Visually monitor the sensor site to ensure the integrity of the skin.
4. If needed, secure the cable with medical tape. Ensure the tape securing the cable does not restrict blood flow.

2.2.2 Y7500

1. Position the Hydrogel stripes HGS7500 (E) or Hydrogel pads HGP7500 (F) at the application site.
2. Place the Y7500 with emitter and detector directly opposing each other on top of the Hydrogel.
3. Apply a Wrap-Tape to fixate the sensor; ensure it is not applied too tight.

2.3 Monitoring

1. Connect the sensor to the monitor. Use an extension cable if needed (G, H)

NOTE: Not suitable for 90-degree sensors plugs

2. Switch on the oximeter and verify proper operation. (Ensure that the oximeter is providing a reading.)
3. Visually monitor the sensor site at least every 4 hours to ensure integrity of the skin and that no pressure injury is developing. Reposition to an alternative location if skin integrity changes.

3 Specifications

3.1 Accuracy Specifications

SpO ₂ Accuracy ¹⁾	Pulse Rate Accuracy
70 – 100 %: $A_{rms} \leq 2.5\%$ (no motion, incl. low perfusion ⁴⁾) ²	$A_{rms} \leq 2\text{bpm}$ (no motion, incl. low perfusion ⁴⁾) ³
60 – 80 %: $A_{rms} \leq 3\%$ (no motion, incl. low perfusion ⁴⁾) ²	
70 – 100 %: $A_{rms} \leq 3\%$ (motion condition) ⁵	$A_{rms} \leq 3\text{bpm}$ (motion condition) ⁵
<60 %: unspecified	

Note: A functional tester may not be used to validate SpO₂ accuracy.

¹⁾ Arterial functional oxygen saturation

²⁾ SpO₂ accuracy is validated by clinical accuracy studies on healthy adult male and female test subjects of age 21 to 32 with skin pigmentation ranging from light to dark over the specified functional oxygen saturation range. Validation was performed on the following monitor: OxyTrue® A with SMARTsat® OEM III

³⁾ Pulse rate accuracy was verified by simulated bench tests to ensure that the entire range was verified. These values were verified on the following monitor: OxyTrue® A with SMARTsat® OEM III

⁴⁾ Tested with Fluke ProSim 8 Oximeter tester at infrared percentage modulation PI: 0.7% to 0.1%

⁵⁾ Tested with all Fluke Index II Oximeter tester motion patterns with pattern specific motion frequency of 0.5Hz to 6Hz at perfusion PI: 0.65% to 5% including non-repetitive motion and motion repeating every 0.5Hz.

3.2 Technical Specifications

Optical	LED	Wavelength	Output Power ⁶
	Red	650 to 670nm	< 10 mW
	Infrared	880 to 930nm	< 10 mW
Materials	The materials that come into contact with the patient have undergone extensive biocompatibility testing to meet ISO 10993. The materials used in the manufacture of the sensor do not contain natural rubber latex. Further information is available on request.		
Water spray protection	IPx2 (Protected from water drops less than 15 degrees from vertical)		

⁶⁾ This information can be especially useful to clinicians.

3.3 Environmental Specification

Operational temperature range:	+0°C to +40°C (32°F to 104°F)
Storage, and transport temperature range:	-20°C to +60°C (-4°F to 140°F)
Relative humidity range:	15 - 95% non-condensing
Atmospheric pressure range:	620 hPa to 1060 hPa (620 hPa corresponding to an altitude of 4000m)
Stabilization Time (from storage to operational conditions):	Up to 20 minutes

Note: These environmental specifications are specific to the sensor. For information on the environmental specification of the compatible monitoring systems, please see the Operator's Manual for that system.

4 Disposal

Do not dispose of with domestic waste. Follow local regulations when disposing or recycling this device and its component parts. A list of component materials is available by contacting bluepoint medical GmbH & Co. KG.

5 Warranty

To obtain information about a warranty, if any, for this product, contact bluepoint medical's Technical Services Department, or your local bluepoint medical representative.

6 Reporting of serious incidents

Report any serious incident that occurs in connection with the use of this product to the competent authority in your country and to the manufacturer. Serious incidents associated with the use of this product are low oxygen saturation, high oxygen saturation, electric shock and death.

Contact manufacturer bluepoint medical directly, or send an email to the following address: prcc@bluepoint-medical.com.

For reporting always provide the following information:

- Order number and model as indicated on the product
- Serial number/batch number
- Date and description of the serious incident, including its impact on the patient or any injury
- Your contact details (name and title, institution, address and phone number)

7 Clinical Studies

In the clinical accuracy studies the measured values of the sensors (SpO_2) were compared with those of the arterial co-oximetry (SaO_2).

The accuracy plots at the end of this booklet show the collected data points plotted as SaO_2 vs the measurement error ($SpO_2 - SaO_2$) together with the information listed below:

- | | |
|-------------|---|
| (1) Y-axis | Measurement error ($SpO_2 - SaO_2$) |
| (2) X-axis | Reference CO-Oximeter reading (SaO_2) |
| (3) — . . — | Upper 95% Limit of Agreement |
| (4) | Mean |
| (5) — . . — | Lower 95% Limit of Agreement |
| (6) | Measured A_{rms} values |

The following monitor/ sensor combinations were tested:

Monitor	Sensor
OxyTrue® A with SMARTsat® OEM III	10-AP ⁷⁾
	10-PP
	10-NP

⁷⁾ Results apply to equivalent sensor 10-IP

Below the plots the measured Arms values per SaO₂ range are listed. Note that the overall claimed accuracy values are listed in section "Accuracy".

Pulse oximeter measurements are statistically distributed, only about two-thirds of measurements can be expected to fall within $\pm$ Arms of the measured values by a co-oximeter.

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© 2017 bluepoint medical GmbH & Co. KG.
SMARTsat® and OxyTrue® are registered trademarks of bluepoint medical GmbH & Co. KG.
®* TM* Third party brands are trademarks of their respective owners.