
SMARTsat®


Reusable SpO₂ Sensors

de	Gebrauchsanweisung	4
en	Instructions For Use	11
fr	Mode d'emploi	17
bg	Инструкции за употреба	24
cs	Návod k použití	31
da	Brugsanvisning	37
el	Οδηγίες χρήσης	44
es	Instrucciones de uso	51
et	Kasutusjuhend	58
fi	Käyttöohjeet	64
hr	Upute za upotrebu	70
hu	Használati útmutató	76
it	Istruzioni per l'uso	83
ja	取扱説明書	90
ko	사용 지침	96
lt	Naudojimo instrukcija	102
lv	Lietošanas instrukcija	108
nl	Gebruiksaanwijzing	114
no	Bruksanvisning	121
pl	Instrukcja użycia	127
pt	Instruções de utilização	134
pt-br	Instruções de uso	141
ro	Instrucțiuni de utilizare	148
ru	Инструкция по применению	155
sk	Návod na použitie	162
sl	Navodila za uporabo	169
sr	Uputstvo za upotrebu	176
sv	Bruksanvisning	182
tr	Kullanım Talimatları	188
zh	使用说明	195


SoftCap®


SoftWrap®

SoftFlap®


EP7500

BF7500


(A)


SoftCap®

SC (>30kg), SCM (10 - 50kg)

SoftCap®


SCP (3 - 10kg)


SoftFlap®

(>30kg)

(B)


SoftCap® SC / SCM
SoftFlap®, SoftWrap®


SoftCap® SCP, SoftWrap®

(C)


SoftCap® SC / SCM / SCP

(D)


SoftCap® SC / SCM / SCP

(E)


SoftCap® SC / SCM / SCP


(F) SoftWrap® + WT


WT






SoftWrap® (>20kg)

SoftWrap® (<3kg)

SoftWrap® (<3kg)

(G) EP7500 (>30kg) + CS**(H) EP7500 (>30kg)****(I) BF7500 (≤5kg)****(J) XT6500 / XT6501****(K)**

XT6500/ XT6501 → ← SoftCap® / SoftFlap® / SoftWrap®

#	REF	KX	#	REF	KX
SC7500	6020132004	1.2m / 4ft	W7500	6020132006	1.2m / 4ft
SC7502	15020132001	0.2m / 0.66ft	W7501	13020132001	0.11m / 0.36ft
SCM7500	6020132010	1.2m / 4ft	W7502	15020132004	0.2m / 0.66ft
SCM7502	15020132002	0.2m / 0.66ft	EP7500	6020132254	1.2m / 4ft
SCP7500	6020132300	1.2m / 4ft	SF7500	6020132002	1.2m / 4ft
SCP7502	15020132003	0.14m / 0.46ft	BF7500	14020132001	1.2m / 4ft
XT6500	1020132275	1.2m / 4ft	XT6501	1020122058	2.4m / 8ft
CS	3020663004	-	WT	15020621001	0.12m / 0.4ft

SMARTsat®

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Reusable Pulse Oximetry Sensors

1 Directions for use

Note: This Instruction for Use is intended for the medical professional. Before use, carefully read the pulse oximetry sensor Instructions for Use, including all warnings, cautions, and instructions. For use in the home care environment the user may use the device after instruction by, or under supervision of a trained healthcare professional.

1.1 Intended purpose

Pulse oximeters including accessories provide measurements of arterial blood oxygen saturation (SpO₂) and pulse rate.

1.2 Indications for use

The SMARTsat® Reusable Pulse Oximetry Sensors are intended for use in non-invasive continuous and/or spot-check monitoring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate on adult, paediatric, infant and neonatal patients, depending on the sensor style used.

The use environment may include: hospitals, hospital-type facilities, intra-hospital transport, mobile emergency medical applications and home use.

Transport environments include intra-hospital transport and ground emergency transport (road ambulances). The use environment is dependent upon the use environments specified for the compatible monitoring system in combination with the sensor used.

The SMARTsat® Reusable Pulse Oximetry Sensors are intended for use with monitoring systems with integrated SMARTsat® technology. For monitor/sensor compatibility, refer to the monitors instructions for use. The SMARTsat® Reusable Pulse Oximetry Sensors are for prescription use only.

1.3 Contraindications

There are no known contraindications to the use of the sensors.

1.4 Warnings and cautions

- Pulse oximeter probes and extension cables are designed for use with specific monitoring systems. Incompatible components can result in degraded performance and electromagnetic incompatibility, or patient injury. The operator is responsible for checking compatibility prior to use.
- A functional tester cannot be used to assess the accuracy of pulse oximeter probes.
- To prevent skin erosion and pressure necrosis, check the patient skin condition at least every 6 hours and reposition to an alternative location if skin integrity changes. Check the measurement site more frequently with patients which may have low tissue perfusion or edema at the measurement site.
- Misapplication of a pulse oximeter probe with excessive pressure (e.g. sensor too small or applied too tight) and for prolonged periods can cause pressure injury and inaccurate readings.
- Do not apply the sensors to areas with poor skin integrity to prevent damage to the skin.
- Routinely clean the sensor according to the instructions and before attaching it to a new patient to prevent cross infection.
- During measurement, the sensor generates red and infrared light with specific fixed wavelengths as listed in the specifications below. Consider that these wavelengths might influence diagnostic parameters of other optical applications.
- Keep the SpO₂ sensor out of the radiation field during computed tomography or full-body radiation to avoid possible inaccurate readings for the duration of the active irradiation period.

- Do not clean the device or sensors whilst in use or connected to a patient.
- Do not sterilize by irradiation, steam, or ethylene oxide – refer to cleaning and disinfection instructions. Use of agents other than specified may damage the sensor.
- Do not use the sensor or extension cable if it is damaged. Use of a damaged sensor or extension cable could cause patient injury or equipment failure.
- Excessive patient motion, excessive ambient light, electromagnetic interference (e.g. stacking and location to other medical equipment), blood pressure cuff on sensor site, venous congestion, abnormal venous pulsations, abnormal pulse rhythms due to physiological conditions or induced through external factors, dysfunctional haemoglobin, intravascular dyes, fingernail polish and artificial fingernails may affect the sensor performance and the accuracy of the measurement.
- Route all cables carefully to reduce the possibility of patient entanglement or strangulation.
- Do not alter or modify the sensor or extension cable. Alterations or modifications may affect performance or accuracy.
- Properly apply sensors according to sensor's directions for use. If the sensor is not applied properly or partially gets loose, this may cause no or inaccurate readings.
- Select a well perfused site for monitoring, very low perfusion at the monitored site may result in no or inaccurate readings.
- Remove the SMARTsat® sensor, extension cable or other oximetry sensors from the patient during magnetic resonance imaging (MRI). Conducted current may cause burns.
Note: Excluded from this warning are W7501 sensors, which are validated for use during MRI and are labelled as such.
- If SpO₂ values indicate hypoxemia, a laboratory blood sample should be taken to confirm the patient's condition.
- Avoid possible interference from sources of electromagnetic interference such as, but not limited to: Cellular phones, radio transmitters, motors, telephones, lamps, electrosurgical units, defibrillators, and other devices. Interference may cause inaccurate measurements. If you are unsure if your device is working properly, contact your clinician.

2 Instructions for use

2.1 Sensor application

2.1.1 SoftCap® and SoftFlap®

Insert the patient's finger into the sensor **(A)**. The preferred finger to use is the index finger **(B)**. Alternative sites may be used for the SoftCap® models, depending on finger sizes **(B)**. See table below for recommended finger digit height **(d)** for each SoftCap® style.

Sensor model		
SC	SCM	SCP
12.5 – 25.5 mm	9 – 19 mm	7 – 9 mm

The sensor is applied correctly when:

- a. The tip of the finger touches the sensor's inside end **(C)**.
- b. The sensor cable extends along the patient's hand **(A)**.

NOTE: Ensure the sensor is not placed backwards **(D)**, or with the finger extending **(E)**.

Open the sensor by pressing at the sides of the SoftCap® housing or at the front part of the SoftFlap® housing, to remove it.

2.1.2 SoftWrap®

Fix a Wrap-Tape to the SoftWrap® sensor **(F)** and apply the sensor to the selected sensor site. The preferred sensor sites are the index finger for adults and paediatrics, the hand for infants and on the foot below the toes for neonates.

The sensor is applied correctly when:

- The emitter and detector are positioned directly opposite each other.
- The sensor is applied firmly, but not too tight.

Open the Velcro® fastener at the Wrap Tape (WT) to remove the SoftWrap® sensor from the finger.

2.1.3 EP7500

The preferred site is the patient's ear lobe, securing the sensor using the cable support **(G)**.

Insert the ear sensor cable into the cable support, leaving at least 10 cm residual cable **(G.1)** to connect the sensor to the fleshy portion of the ear lobe.

Hook the cable support over the patient's ear and clip the sensor to the fleshy portion of the ear lobe **(G.2)**.

Alternatively, position the ear sensor on the ear pinna for spot-check measurement **(H)**.

The sensor is applied correctly when:

- The emitter and detector are completely covered by the ear/earlobe.
- No tensile force is exerted on the sensor cable.

To remove, press the back housing of the EP7500 and detach the sensor together with the cable support (CS) from the ear.

2.1.4 BF7500

Position the sensor between the thumb and index finger and position the sensor on the foot below the toes **(I.1)**.

Apply slight pressure to both sides of the sensor to ensure that the detector and emitter are in good contact with the skin during the spot-check measurement **(I.2)**.

The sensor is applied correctly when the sensor is applied firmly, but not too tight during the spot-check measurement.

2.2 Monitoring

Connect the sensor to the monitor. Use an extension cable if needed **(J)**, **(K)**.

NOTE: Not suitable for 90-degree sensors plugs

Switch on the oximeter and verify proper operation. (Ensure that the oximeter is providing a reading.)

Visually monitor the sensor site at least every 6 hours to ensure integrity of the skin and that no pressure injury is developing. Reposition to an alternative location if skin integrity changes.

2.3 Sensor removal

After monitoring, remove the sensor according to the instructions above.

Disconnect the sensor from the monitor and follow the cleaning instructions, if necessary.

Store the sensor away for its next use.

3 Cleaning and disinfection

It is necessary to clean and disinfect the sensor before attaching it to a new patient. Routine cleaning or low-level disinfection is recommended when using the sensor on the same patient.

3.1 SoftCap®, SoftWrap®, BF7500

3.1.1 To clean the sensor with 70% Isopropyl alcohol:

- Disconnect the sensor from the monitor before cleaning or disinfecting.

- b. Rinse the sensor and cable with distilled water for at least 2 minutes.
Caution: Do not immerse the connector of the sensor cable in liquid.
- c. Clean the surface of the cable by wiping it with a cotton pad soaked in 70% Isopropyl alcohol.
- d. Immerse the sensor housing in 70% Isopropyl alcohol for 5 minutes and let it air dry for at least 10 minutes.
- e. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the cleaning process.

3.1.2 To disinfect (low-level) with 1:10 bleach dilution

- a. Disconnect the sensor from the monitor before cleaning or disinfecting.
- b. Rinse the sensor and cable with distilled water for at least 2 minutes.
Caution: Do not immerse the connector at the sensor cable in liquid.
- c. Prepare 1:10 bleach dilution in distilled water. Wipe the cable with a cotton pad soaked with the solution of diluted bleach.
- d. Immerse the sensor housing in the solution of diluted bleach for 5 minutes.
- e. Rinse the sensor and cable by immersing them in distilled water for 5 minutes.
- f. Allow 10 minutes for the sensor, cable, and connector to dry after rinsing.
- g. Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition is reached.

3.1.3 To disinfect (high-level) - CIDEX OPA™*

- a. Rinse sensor to remove surface dirt – the sensor housing may be fully immersed in liquid.
Caution: Do not immerse the connector at the sensor cable in liquid.
- b. Clean the sensor and patient contact surfaces. Prolystica™* is recommended as a cleaning agent. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution and wipe the inner and outer surfaces with a soft brush or cloth to remove any visible soil. Rinse sensor with water following cleaning and wipe dry prior to disinfection. Visually inspect the sensor for cleanliness. If visible soil is present, repeat the cleaning process.
- c. CIDEX OPA™* is recommended for high-level disinfection. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution for 12 minutes and remove.
- d. Thoroughly rinse the sensor with water following disinfection (5 minutes minimum distilled water rinse/soak).
- e. Manually dry the sensor and verify the sensor, cable, and connector are dry prior to use.

3.2 SoftFlap®, EP7500

3.2.1 To clean the sensor with 70% Isopropyl alcohol

- a. Disconnect the sensor from the monitor before cleaning or disinfecting.
- b. Clean the surface of the cable and sensor head at the inside and outside by wiping it with a cotton pad soaked in 70% Isopropyl alcohol.
- c. Allow 10 minutes for the sensor, cable, and connector to dry after cleaning.
- d. Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition is reached.

4 Specifications

4.1 Accuracy Specifications

SpO ₂ Accuracy ¹⁾	Pulse Rate Accuracy
70 – 100 %: $A_{rms} \leq 2\%$ (no motion, incl. low perfusion ⁴⁾ ²⁾	$A_{rms} \leq 2\text{bpm}$ (no motion, incl. low perfusion ⁴⁾ ³⁾
60 – 80 %: $A_{rms} \leq 2.5\%$ (no motion, incl. low perfusion ⁴⁾ ²⁾	$A_{rms} \leq 3\text{bpm}$ (motion condition) ⁵⁾
70 – 100 %: $A_{rms} \leq 3\%$ (motion condition) ⁵⁾	
<60 %: unspecified	

Note: A functional tester may not be used to validate SpO₂ accuracy.

¹⁾ Arterial functional oxygen saturation

²⁾ SpO₂ accuracy is validated by clinical accuracy studies on healthy adult male and female test subjects of age 21 to 32 with skin pigmentation ranging from light to dark over the specified functional oxygen saturation range. Validation was performed on the following monitor: OxyTrue® A with SMARTsat® OEM III

³⁾ Pulse rate accuracy was verified by simulated bench tests to ensure that the entire range was verified. These values were verified on the following monitor: OxyTrue® A with SMARTsat® OEM III

⁴⁾ Tested with Fluke ProSim 8 Oximeter tester at infrared percentage modulation PI: 0.7% to 0.1%

⁵⁾ Tested with all Fluke Index II Oximeter tester motion patterns with pattern specific motion frequency of 0.5Hz to 6Hz at repeating PI: 0.65% to 5% including non-repetitive motion and motion repeating every 0.5Hz.

4.2 Technical Specifications

Optical	LED	Wavelength	Output Power ⁶⁾
	Red	650 to 670nm	< 10 mW
	Infrared	880 to 930nm	< 10 mW
Materials	The materials that come into contact with the patient have undergone extensive biocompatibility testing to meet ISO 10993. The materials used in the manufacture of the sensor do not contain natural rubber latex. Further information is available on request.		
Water spray protection	IP34 (Protected from water spray from any direction)		
Expected service life	2 years		

⁶⁾ This information can be especially useful to clinicians.

4.3 Environmental Specification

Operational temperature range:	+0°C to +40°C (32°F to 104°F)
Transient operating conditions:	Sensor will operate for 20 minutes at -20°C and +50°C temperature conditions.
Storage, and transport temperature range:	-40°C to +70°C (-40°F to 158°F)
Relative humidity range:	15 - 95% non-condensing
Atmospheric pressure range:	620 hPa to 1060 hPa (620 hPa corresponding to an altitude of 4000m)
Stabilization Time (from storage to operational conditions):	Up to 20 minutes

Note: These environmental specifications are specific to the sensor. For information on the environmental specification of the compatible monitoring systems, please see the Operator's Manual for that system.

5 Disposal

Do not dispose of with domestic waste. Follow local regulations when disposing or recycling this device and its component parts. A list of component materials is available by contacting bluepoint medical GmbH & Co. KG.

6 Warranty

To obtain information about a warranty, if any, for this product, contact bluepoint medical's Technical Services Department, or your local bluepoint medical representative.

7 Reporting of serious incidents

Report any serious incident that occurs in connection with the use of this product to the competent authority in your country and to the manufacturer. Serious incidents associated with the use of this product are low oxygen saturation, high oxygen saturation, electric shock and death.

Contact manufacturer bluepoint medical directly, or send an email to the following address: prcc@bluepoint-medical.com.

For reporting always provide the following information:

- Order number and model as indicated on the product
- Serial number/batch number
- Date and description of the serious incident, including its impact on the patient or any injury
- Your contact details (name and title, institution, address and phone number)

8 Clinical Studies

In the clinical accuracy studies the measured values of the sensors (SpO_2) were compared with those of the arterial co-oximetry (SaO_2).

The accuracy plots at the end of this booklet show the collected data points plotted as SaO_2 vs the measurement error ($SpO_2 - SaO_2$) together with the information listed below:

- | | |
|-------------|---|
| (1) Y-axis | Measurement error ($SpO_2 - SaO_2$) |
| (2) X-axis | Reference CO-Oximeter reading (SaO_2) |
| (3) — . . — | Upper 95% Limit of Agreement |
| (4) | Mean |
| (5) — . . — | Lower 95% Limit of Agreement |
| (6) | Measured A_{rms} values |

The following monitor/ sensor combinations were tested:

Monitor	Sensor
OxyTrue® A with SMARTsat® OEM III	SoftCap® SC7500 ⁷⁾
	SoftFlap® SF7500
	SoftWrap® W7500 ⁸⁾
	Ear Probe EP7500

⁷⁾ Results apply to equivalent sensors SoftCap® with medium and small housing including BF7500.

⁸⁾ Results apply to equivalent SoftWrap® sensor models (W7500, W7501, W7502)

Below the plots the measured A_{rms} values per SaO_2 range are listed. Note that the overall claimed accuracy values are listed in section "Accuracy".

Pulse oximeter measurements are statistically distributed, only about two-thirds of measurements can be expected to fall within $\pm A_{rms}$ of the measured values by a co-oximeter.

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© 2017 bluepoint medical GmbH & Co. KG.
SMARTsat® and OxyTrue® are registered trademarks of bluepoint medical GmbH & Co. KG.
®* TM* Third party brands are trademarks of their respective owners.