

# Temperature Probes Revision Questions (Edition 1.1)

Please review the Introduction (DOCID 109087) and Technical (DOCID 111263) training notes in the document index and answer the following questions.

**Please be aware this is edition 1.1, revised on 10th September 2024.**

Email \*

viamed.cathy.green@gmail.com

**What is the purpose of monitoring patient temperature during surgery? \***

1 point

- ☒ To detect potential health issues like hypothermia or fever
- ☐ To track blood pressure fluctuations
- ☐ To monitor oxygen saturation
- ☐ To check heart rate during anaesthesia

**Which one of these do not measure core body temperature? \***

1 point

- ☒ Nasopharynx
- ☐ Oesophagus
- ☐ Skin
- ☐ Rectum

**What are the two main types of Bluepoint 400-AC temperature probes? \***

1 point

- ☒ Oesophageal/Rectal and Skin/Surface
- ☐ Nasopharyngeal and Surface
- ☐ Disposable and Reusable
- ☐ Handheld and Digital

**Which of the following patient sizes are available for the Bluepoint 400-AC temperature probes?**

\* 1 point

- ☒ Adult, Paediatric, and Neonatal
- ☐ Adult, Infant, and Senior
- ☐ Infant, Child, and Adult
- ☐ Newborn, Child, and Adolescent

**What material are Bluepoint 400-AC temperature probes made from to prevent latex allergies?**

\* 1 point

- ☐ Acrylic
- ☒ Santoprene
- ☐ Silicone
- ☐ Polycarbonate

**How should 400-AC temperature probes not be cleaned between uses? \***

1 point

- ☐ Hand wash only
- ☐ Steam sterilization using an autoclave
- ☐ Rinsed with alcohol
- ☒ Soaked in bleach

**Why is monitoring core temperature critical during anaesthesia? \***

1 point

- ☐ To check if patient is fully anaethetised
- ☒ Anaesthesia makes it hard for patients to maintain core temperature.
- ☐ To check SpO2 levels are accurate
- ☐ To check heart rate

**Which leading medical device manufacturers are Bluepoint 400-AC probes compatible with?**

\* 1 point

- ☒ YSI, Siemens, GE Healthcare, Dräger
- ☐ YSI, IBM, Samsung, Siemens
- ☐ Dräger, Siemens, Microsoft, Apple
- ☐ GE Healthcare, Apple, Philips, Hewlett Packard

**What is the typical warranty period for Bluepoint 400-AC temperature probes? \***

1 point

- ☐ 3 months
- ☐ 6 months
- ☒ 12 months
- ☐ Lifetime

**Which of the following is untrue? \***

1 point

- ☐ They are steam autoclavable
- ☐ They can be cleaned specified cleaning agents
- ☒ They are disposable
- ☐ They can be mechanically cleaned in an automated washer/disinfector

This content is neither created nor endorsed by Google.

Google Forms