

Neonatal Resuscitation Limbs

PIP (cmH ₂ O)	² NeoFlow	0.0	3.6	5.3	11.8	14.7	19.8	5.0
	³ NeoPuff	0.0	0.0	0.0	3.6	15.1	20.2	18.4
	⁴ Intersurgical	0.0	2.9	5.7	13.3	15.0	20.0	19.7
		0	25%	50%	75%	100%	Cap Fully CLOSED and Bleed Hole Fully Occluded	Cap Fully OPEN and Bleed Hole Fully Occluded
¹ PIP cap closure (cmH ₂ O) (0=fully open 100%=fully closed)								

¹positive inspiratory pressure

Test Details Test date: 16th March 2015 Test conducted: Armstrong Medical Limited Delivery device: Neopuff™ RD900 (Fisher Paykel Healthcare, Auckland, NZ) Oxygen line pressure: 58psi Flow rate: 10L.min Maximum PIP setting: 20cmH ₂ O Manometer: GiO1 Digital Pressure Gauge, (GiO Healthcare, Taiwan, ROC)	Breathing Circuit ² NeoFlow AMVC1871/167 (Armstrong Medical, Coleraine, NI) ³ Neopuff™ 900RD010 (Fisher Paykel Healthcare, Auckland, NZ) ⁴ Intersurgical® 6433000 (Wokingham, UK)
---	--

Neonatal Resuscitation Limbs

Effect of Cap Rotation on Positive Inspiratory Pressure

