

Operating Instructions



INFANT T-PIECE RESUSCITATOR 900 SERIES

Intended Use

The Fisher & Paykel Healthcare Neopuff™ Infant Resuscitator is an easy-to-use manually operated, gas-powered resuscitator which provides controlled and accurate resuscitation of newborn babies in delivery suites, nurseries and neonatal intensive care units.

Symbol	Definition
--------	------------




Attention: Consult the Operating Instructions.

Warnings:


- Please read and understand the instructions fully before using the Neopuff™ Infant Resuscitator and related accessories. The Neopuff Infant Resuscitator is to be used only by persons trained in infant resuscitation. It is recommended that users refer to ILCOR/AHA/ERC guidelines to determine the suitability of different types of resuscitator for use.
- It is the responsibility of the purchaser to ensure that all users of this device have been adequately trained in resuscitation techniques.
- The Neopuff Infant Resuscitator must only be used after checking that correct pressures will be delivered to the baby.
- Ensure no smoking, naked flames or sources of ignition are present while the unit is in use.
- For connection to flow-regulated oxygen or oxygen/air mixture only.
- Input gas flow rate 5 to 15 L/min. Recommended operating gas flow rate 8 L/min. Do not attempt to use a higher flow than 15 L/min. Input flow ranges are circuit specific, refer to circuit User Instructions.
- The Max Pressure Relief can be adjusted up to a nominal 80 cmH₂O [mbar], and should only be done in exceptional circumstances by persons trained in infant resuscitation. Do not attempt to set the Max Pressure Relief above 80 cmH₂O [mbar].
- Use only recommended Fisher & Paykel Healthcare Infant Resuscitator accessories.
- Use only a Fisher & Paykel Healthcare Gas Supply Line or approved equivalent.
- Ensure all oxygen and air supplies are turned off and disconnected from the Neopuff before performing cleaning procedures. Explosion and fire hazards can exist when performing cleaning procedures in an oxygen-enriched environment.
- The black Test Lung contains natural rubber latex which may cause allergic reactions.
- US Federal law restricts this device to sale by or on the order of a licensed physician.
- An alternative means of resuscitation must be available.

SET-UP




The following procedure should be carried out prior to every use of the Neopuff™ to ensure that the device is functioning correctly.

1. Check manometer reads zero with no gas flow. If not, the manometer requires calibration (refer to section 3.3.4 of the Technical Manual).




2. Connect Gas Supply

Connect an oxygen or blended oxygen/air supply to the gas inlet port using the Gas Supply Line.



3. Connect T-Piece Circuit

- Connect the T-Piece Circuit to the gas outlet port.
- Leave blue cap in place on the T-Piece Circuit or connect Test Lung to T-Piece Circuit (before use, inspect Test Lung for signs of damage such as discoloration).




4. Check Settings

Adjust the gas supply to desired flow rate between 5 and 15 L/min.

Note:


Ensure the oxygen concentration of an oxygen/air supply is either monitored using an oxygen analyzer or preset using oxygen/air flow rate graphs.



5A. To Check Max Pressure

Occlude PEEP¹ cap and turn PIP² control fully clockwise, until the knob doesn't turn anymore.


1. Positive End Expiratory Pressure 2. Peak Inspiratory Pressure



5B. Adjust max pressure control knob clockwise or counterclockwise to set desired max pressure.

Note:


- The factory setting of the Max Pressure Relief is 40 cmH₂O [mbar].
- The Max Pressure Relief valve acts as an overall limit on the achievable circuit pressure. Resuscitation above 40 cmH₂O [mbar] cannot be achieved unless the Max Pressure Relief valve is adjusted.



6. To Set PIP*

While still occluding the PEEP cap, turn PIP control knob counterclockwise until the desired peak inspiratory pressure is set.

**Neopuff™ units manufactured from June 2010 have a faster PIP adjustment.*




7. To Set PEEP


Adjust PEEP cap to the desired PEEP level.

Turn off gas supply and remove Test Lung from T-Piece. Ensure that the rigid plastic connector of the Test Lung is also removed from the T-Piece before attempting to connect a mask or endotracheal tube. Failing to do so may cause unacceptable delays during patient resuscitation.


TO RESUSCITATE



1. Adjust gas supply to the desired flow rate.



2. Fit T-Piece to neonatal resuscitation mask and place over the baby's mouth and nose.



OR Fit T-Piece to the endotracheal tube.

3. Resuscitate by placing and removing thumb over the PEEP¹ cap to allow inspiration and expiration.

1. Positive End Expiratory Pressure

CLEANING AND SERVICING

- Clean external surfaces of the Neopuff™ Infant Resuscitator and Gas Supply Line using a damp cloth and mild soapy water or Isopropyl Alcohol.
- Dry all surfaces after cleaning with a clean soft cloth or paper towel.
- The Neopuff should require minimal servicing or maintenance when used under normal conditions.
- For more information on cleaning and maintenance of the Neopuff Infant Resuscitator, Test Lung and other accessories, please refer to the Technical Manual (Part No. 185041597).

PERFORMANCE SPECIFICATIONS

Recommended Body Weight Range	Up to 10 kg
Manometer Range	-10 to 80 cmH ₂ O [mbar]
Peak Inspiratory Pressure (PIP)	@ 5 L/min approx. 2 to 70 cmH ₂ O [mbar]
	@ 8 L/min approx. 3 to 72 cmH ₂ O [mbar]
	@ 10 L/min approx. 4 to 73 cmH ₂ O [mbar]
	@ 15 L/min approx. 8 to 75 cmH ₂ O [mbar]
Positive End Expiratory Pressure (PEEP)	@ 5 L/min approx. 1 to 6 cmH ₂ O [mbar]
	@ 8 L/min approx. 1 to 10 cmH ₂ O [mbar]
	@ 10 L/min approx. 2 to 15 cmH ₂ O [mbar]
	@ 15 L/min approx. 4 to 17 cmH ₂ O [mbar]
Gas Inlet Flow Range	5 L/min (min) to 15 L/min (max)
Operating Time (400 L cylinder)	50 minutes (typical value based on a gas flow rate of 8 L/min)
NOTE: All performance figures listed above are representative only. PEEP values stated are based on typical clinical PIP settings. Higher PEEP values can be achieved if higher PIP values are set.	

For more information please contact
your local Fisher & Paykel Healthcare representative

Manufacturer 

Fisher & Paykel Healthcare Ltd
15 Maurice Paykel Place
East Tamaki, Auckland 2013

PO Box 14 348, Panmure
Auckland 1741
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.com
Web: www.fphcare.com

Australia

Fisher & Paykel Healthcare Pty Limited
36-40 New Street, PO Box 167
Ringwood, Melbourne
Victoria 3134, Australia

Tel: +61 3 9879 5022
Fax: +61 3 9879 5232

Austria

Tel: 0800 29 31 23
Fax: 0800 29 31 22

Benelux

Tel: +31 40 216 3555
Fax: +31 40 216 3554

China

Tel: +86 20 3205 3486
Fax: +86 20 3205 2132

France

Tel: +33 1 6446 5201
Fax: +33 1 6446 5221

Germany

Tel: +49 7181 98599 0
Fax: +49 7181 98599 66

India

Tel: +91 80 4284 4000
Fax: +91 80 4123 6044

Irish Republic

Tel: 1800 409 011

Italy

Tel: +39 06 7839 2939
Fax: +39 06 7814 7709

Spain

Tel: +34 902 013 346
Fax: +34 902 013 379

Sweden

Tel: +46 8 564 76 680
Fax: +46 8 36 63 10

Switzerland

Tel: 0800 83 47 63
Fax: 0800 83 47 54

Taiwan

Tel: +886 2 8751 1739
Fax: +886 2 8751 5625

Turkey

Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK 

Fisher & Paykel Healthcare Ltd
Unit 16, Cordwallis Park
Clivemont Road, Maidenhead
Berkshire SL6 7BU, UK

Tel: +44 1628 626 136
Fax: +44 1628 626 146

USA/Canada

Tel: +1 800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001