Instructions for use: 4000 Series Disposable Pulse Oximetry Sensors

DIRECTIONS FOR USE

Intended Use: Viamed 4000 Series oximetry sensors are indicated for use in continuous noninvasive monitoring of arterial oxygen saturation and pulse rate. They are original or fully compatible replacement sensors intended for use with major brands of pulse oximeters or patient monitors. The table at the back, lists the Viamed 4000 Series disposable sensors available for the following four categories of patients: adult (weight > 30 kg), paediatric (10–50 kg weight), infant (3–20 kg weight), and neonatal (weight <3 kg). Matching the brand of oximeter you are using, select the model number corresponding to the patient and tape type.

Precautions: For further instructions, warnings and precautions refer to the operating instructions provided with your oximeter. Do not autoclave or immerse the sensor in liquid. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings: Match the sensor model with your brand of oximeter. A mismatch may cause patient injury, inaccurate measurements or equipment damage. Check the sensor and reposition to an alternative location every 4 hours. Avoid application of the sensors to: edematous or fragile tissue, excessive patient motion, excessive ambient light, electromagnetic interference, dysfunctional hemoglobin, low perfusion, intravascular dyes, finger nail polish and long or artificial finger nails which may affect the sensor performance and the accuracy of the measurement. Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure. The sensor is contraindicated for use on patients who exhibit allergic reaction to the adhesive tape. Note: The materials used in the manufacture of the sensor contain no natural latex protein. The materials have undergone extensive bio-compatibility testing. Further information is available upon request.

Application of the Sensor


1. Select the sensor site on the patient. The preferred sensor sites are the index finger for adults and paediatrics, the great toe and on the foot below the toes for neonates.
2. Remove the plastic backing from the sensor adhesive wrap.
3. Apply the sensor on the patient:
 - Position the emitter and detector directly opposite each other.
 - Apply the sensor firmly but not too tightly.
4. Connect the sensor cable to the adaptor cable that is suitable for your oximeter or patient monitor. Then switch on the oximeter and verify proper operation
5. Visually monitor the sensor site to ensure the integrity of the skin.
6. Secure the lead with medical tape.

Adult/Paediatric

Infant

Neonatal

Environmental Requirements:

- Operating temperature range:
10° - 40° C (50°-104° F)
- Storage and transportation temperature range
0° - 50° C (32°-122° F)
- Relative humidity:
0 - 75% non-condensing

SpO₂ Accuracy:

70 - 100%: A_{rms} <3.0%

0 - 69%: Unspecified

Accuracy test methods available upon request.

LED Specifications

	Wave Length	Radiant Power (typical)	Brand
Red	660nm	3.5 mW	All below
Infrared	885nm	3.5 mW	Nellcor, CSI
	905nm	3.5 mW	Datex, BCI, Viamed
	910nm	3.5 mW	Nonin
	940nm	3.5 mW	Ohmeda, Novamatrix

Disposal: Follow local regulations governing disposable recycling plans regarding disposal or recycling of device components. A list of component materials is available by contacting Viamed.

All registered trademarks recognized.

<i>Nellcor®</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	1-AF	1-PF	1-IF	1-NF
Transpore	1-AT	1-PT		
Plaster	1-AP	1-PP	1-IP	1-NP

<i>Ohmeda®</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	2-AF	2-PF	2-IF	2-NF
Transpore	2-AT	2-PT		
Plaster	2-AP	2-PP	2-IP	2-NP

<i>Datex</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	3-AF	3-PF	3-IF	3-NF
Transpore	3-AT	3-PT		
Plaster	3-AP	3-PP	3-IP	3-NP

<i>CSI</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	4-AF	4-PF	4-IF	4-NF
Transpore	4-AT	4-PT		
Plaster	4-AP	4-PP	4-IP	4-NP

<i>BCI®</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	5-AF	5-PF	5-IF	5-NF
Transpore	5-AT	5-PT		
Plaster	5-AP	5-PP	5-IP	5-NP

<i>Novamatrix®</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	6-AF	6-PF	6-IF	6-NF
Transpore	6-AT	6-PT		
Plaster	6-AP	6-PP	6-IP	6-NP

<i>Nonin</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	7-AF	7-PF	7-IF	7-NF
Transpore	7-AT	7-PT		
Plaster	7-AP	7-PP	7-IP	7-NP

<i>Viamed</i>	<i>Adult</i>	<i>Paediatric</i>	<i>Infant</i>	<i>Neonatal</i>
Micro Foam	8-AF	8-PF	8-IF	8-NF
Transpore	8-AT	8-PT		
Plaster	8-AP	8-PP	8-IP	8-NP

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
Viamed	Adult	Paediatric	Infant	Neonatal
Micro Foam	8-AF	8-PF	8-IF	8-NF
Transpore	8-AT	8-PT		
Plaster	8-AP	8-PP	8-IP	8-NP

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Viamed Ltd