

P 969E8 (CRITICARE) EXT.


MALE END. (6 PIN LETO)

STEP 1- CUT 8' OF CABLE, PLACE CLEAR HEAT SHRINK OVER CABLE & STRAIN RELIEF

STEP 2- STRIP OUTER & INNER CABLES BACK TO AFX 1" REMOVED

STEP 3- PUT BOTH SHIELDS TOGETHER THEN PUT CABLE GRIP OVER TOP OF SHIELD (MAKE SURE SHIELD HAS FULL CONTACT WITH GRIP)

STEP 4- STRIP AND TIAV ENDS OF WIRES THEN SOLDER INTO PIN AS FOLLOWS -



1= YELLOW
2= BLACK
3= WHITE
4= RED
5= ORANGE