






Document Title: **Work Instructions, MaxO2ME Analyzer Assembly**Document Number: **WI-4066**Revision Number: **M**



## 1.0 SCOPE

- 1.1 This document shall provide the guidelines and documentation requirements for the final assembly of the MaxO2ME analyzers.






## 2.0 MANUFACTURING INSTRUCTIONS

ALERTS:					
 QUALITY	 SAFETY	 ACCEPT / CORRECT	 REJECT / INCORRECT	 RECORD / SIGN	
<b>SAFETY EQUIPMENT REQUIRED:</b> <ol style="list-style-type: none"> <li>1. Safety glasses</li> <li>2. Smock</li> <li>3. Electrostatic Discharge (ESD) wrist strap and work mat.</li> </ol>					
<b>MATERIALS REQUIRED:</b> <ol style="list-style-type: none"> <li>1. Material Pick List Components</li> <li>2. Syringe filled with RP39P28 Gasket material</li> <li>3. Cotton Swabs</li> </ol>					
<b>EQUIPMENT AND TOOLS REQUIRED:</b> <ol style="list-style-type: none"> <li>1. Pneumatic Torque Driver @ 60 in-oz</li> <li>2. Manual Torque Driver @ 50 in-oz</li> <li>3. Soldering Station</li> <li>4. Clippers</li> <li>5. Speaker Soldering Fixture MT390</li> <li>6. Monitor Test System MT367</li> <li>7. 4x "AA" Batteries</li> <li>8. Tweezers</li> <li>9. Pliers</li> <li>10. X-Acto Knife</li> <li>11. T-8 Torx Driver</li> <li>12. #1 Phillips Head Driver</li> </ol>					









TASK	INSTRUCTION
SETUP	  <p>ENSURE all components diagramed above are accounted for prior to assembly.</p>






TASK	INSTRUCTION
<b>1</b>	<p>Print labels according to LCS.</p> <p> <b>Mark completion on in-process form</b></p>
<b>2</b>	<p>Inspect labels according to LCS</p> <p> <b>Mark completion on in-process form</b></p> <p> Note: personal inspecting labels must be different from the one who printed the labels.</p>
<b>3</b>	<p>SEPARATE Printed Circuit Boards (PCBs) from panel with pliers.</p> <p> REMOVE excess board material. Use caution not to damage any board traces.</p> <div data-bbox="940 940 1479 1717">  </div>



TASK	INSTRUCTION
<b>4</b>	<p>TURN <b>ON</b> soldering station.</p> <p>SET temperature to approximately 800°F.</p> <p> TIN soldering tip as necessary.</p> <p> CAUTION: HOT</p> 
<b>5</b>	<p>PLACE Speaker leads into PCB.</p> <p>ENSURE <b>+</b> and <b>-</b> leads match with corresponding marks on PCB. And that the speaker sits flat.</p> <p>PLACE PCB/Speaker on Soldering Fixture MT390.</p>   





TASK	INSTRUCTION
<b>6</b>	<p>TOUCH Solder wire to speaker lead and PCB joint.</p> <p>TOUCH Iron to tip of solder wire.</p> <p>TWIST Iron ¼ turn to deposit solder.</p> <p>REPEAT on remaining lead.</p> <p>CLEAN solder joints with alcohol.</p> <p>REMOVE PCB from fixture.</p> <p> <b>Mark completion on in-process form</b></p> 
<b>7</b>	<p>REMOVE center foam from Photocell gasket. Leave adhesive backing on until ready to place on the PCBA</p> <p>PLACE foam gasket over LED.</p> <p>REMOVE adhesive backing with tweezers and PRESS gasket into place.</p> 

Document Title: **Work Instructions, MaxO2ME Analyzer Assembly**Document Number: **WI-4066**Revision Number: **M**




TASK	INSTRUCTION
<b>8</b>	<p>REMOVE inner foam section from PCB cover gasket</p> <p>REMOVE adhesive backing</p> <p>PLACE foam gasket over switch</p> <p>TOGGLE switch to the <b>LEFT</b></p> <p><input checked="" type="checkbox"/> Mark completion on in-process form</p>









TASK	INSTRUCTION
<p><b>9</b></p>	<p>UNWRAP Enclosures</p> <p>PLACE Spring in Battery Door</p> <p>THREAD T-Screw into Battery Door</p> <p> <b>Mark completion on in-process form</b></p> <div data-bbox="911 331 1511 722"> <p>REAR ENCLOSURE</p> <p>SPRING</p> </div> <div data-bbox="920 737 1503 1161"> <p>BATTERY DOOR &amp; KICKSTAND</p> <p>T-SCREW</p> </div>
<p><b>10</b></p>	<p>SEPARATE Rear Enclosure and Battery Door from Front Enclosure</p> <p>APPLY Buzzer Gasket.</p> <p>APPLY Battery Foam to door</p> <p> <b>Mark completion on in-process form</b></p> <div data-bbox="937 1249 1487 1602"> </div>








TASK	INSTRUCTION
<b>11</b>	<p>REMOVE adhesive backing and Rear Screen Protector Film from membrane switch.</p> 
<b>12</b>	<p>INSERT ribbon cable through front enclosure slot.</p> <p>PRESS the membrane switch into recessed area all around the perimeter.</p> <p> <b>VERIFY</b> no bubbles are present on the membrane switch. <b>REPLACE</b> if bubbles persist.</p> 








TASK	INSTRUCTION
<b>13</b>	<p>FILL syringe with gasket material</p> <p>CUT tip of syringe</p> <p>APPLY a bead of gasket material to the identified area behind ribbon cable.</p> <p> WIPE off any excess material that could interfere with PCB screen.</p> <p> <b>Mark completion on in-process form</b></p> 
<b>14</b>	<p>REMOVE PCB screen protector</p> 



TASK	INSTRUCTION
<b>15</b>	<p>PLACE PCB onto housing.</p> <p>SECURE PCB to housing with 4 Phillip-head screws.</p> <p>TIGHTEN to 50 in-oz with manual torque driver</p> <p> CHECK LCD for discoloration. If discoloration occurs back off screws by ¼ turn each and RECHECK. If discoloration continues reject.</p> 
<b>16</b>	<p> Record the board serial number on the Unit to Board Serial Number form.</p>  



TASK	INSTRUCTION
<b>17</b>	<p data-bbox="305 331 613 363">CONNECT ribbon cable.</p> <p data-bbox="305 720 651 751">PUSH tabs to lock in place.</p> <p data-bbox="305 758 358 804"></p> <p data-bbox="362 835 776 867">ENSURE ribbon cable is straight.</p> <p data-bbox="305 947 786 1035"> <b>Mark completion on in-process form</b></p> <div data-bbox="971 331 1453 667"></div> <div data-bbox="963 699 1412 993"></div> <div data-bbox="982 1031 1380 1325"></div>



TASK	INSTRUCTION
<b>18</b>	<p>APPLY red Loctite #262 to bottom 1/3<sup>rd</sup> of sensor jack threads.</p> <p>PLACE metal washer onto threads</p> <p>THREAD nut on until it bottoms out.</p> <p> WIPE off any excess Loctite</p> <p> This may be performed as a batch process during material setup.</p> <div data-bbox="971 327 1451 697"></div> <div data-bbox="971 743 1451 1138"></div>



**19**

CONNECT Jumper to PCB


CONNECT Sensor Jack to PCB

PLACE Sensor Jack into housing


CONNECT Battery Cable to PCB



INSERT Battery cable contacts into enclosure.







 Contact wires face **DOWN**

PRESS contacts into place with flat edge of tweezers.






 MOVE wires so they do not get pinched during rear enclosure assembly.

 **Mark completion on in-process form**









TASK	INSTRUCTION
<b>20</b>	<p>CLIP off protrusion from battery contacts.</p> <p> This can be performed as a batch process during material setup.</p> <p>INSERT <b>1x Left</b> and <b>2x Right</b> Battery contacts.</p> <p> Clipped ends face <b>DOWN</b></p> <p>PRESS contacts into place with flat edge of tweezers.</p> <p>INSERT Power Jack Cover</p> <p> <b>Mark completion on in-process form</b></p> <div></div>



TASK	INSTRUCTION
<b>21</b>	<p data-bbox="305 449 781 558">APPLY continuous bead of gasket material to identified inside edges of REAR enclosure.</p> <p data-bbox="305 911 735 953"> WIPE off any excess material.</p> <p data-bbox="305 1262 808 1335">INSERT rubber plug into the slot of the FRONT enclosure.</p> <p data-bbox="305 1415 789 1499"> <b>Mark completion on in-process form</b></p> <div data-bbox="873 331 1544 911"></div> <div data-bbox="898 968 1450 1451"></div>



TASK	INSTRUCTION
<b>22</b>	<p>Prior to placing the rear enclose onto the front enclosure, place the unit label on the rear encloser in the upper left corner of the white/grey plastic area leaving a small gap from the blue colored portion.</p> <p> <b>Record the Unit Number on the Unit to Board Serial Number form matching it to the Board Serial Number.</b></p> 
<b>22</b>	<p>PLACE rear enclosure onto front enclosure.</p> <p>FASTEN with 6 screws. USE pneumatic torque driver set to 60+/-10 in-oz</p> <ul style="list-style-type: none"> <li> <b>TIGHTEN</b> screw over sensor jack housing <b>first</b> to hold in place.</li> <li> If holes become stripped replace enclosure.</li> <li> A manual torque driver can also be used if necessary.</li> <li> <b>WIPE</b> off any excess gasket material that squeezes out after assembly.</li> </ul> 