

Risk Assessment iaw EN ISO 14971:2000 Annex D : Possible hazards with medical devices.

Ref.	Hazard.	Related part / Component posing risk.	Sev of Haz.	Like of Haz.	Risk.	Solution.	Document referenced.	Sev of Haz.	Like of Haz.	Risk.
D.2 Energy hazards and contributory factors										
D.2.1	Electricity	N/A	1	1	1			1	1	1
D.2.2	Heat	If battery polarity reversed	3	1	3	Protection diode added in circuit	(YZ) Design File	1	1	1
D.2.3	Mechanical force	N/A	2	1	2			2	1	2
D.2.4	Ionising radiation	N/A	1	1	1			1	1	1
D.2.5	Non ionising radiation	N/A	1	1	1			1	1	1
D.2.6	Moving parts	N/A	1	1	1			1	1	1
D.2.7	Unintended motion	N/A	1	1	1			1	1	1
D.2.8	Suspended masses	N/A	1	1	1			1	1	1
D.2.9	Patient support failure	N/A	1	1	1			1	1	1
D.2.10	Pressure (vessel rupture)	N/A	1	1	1			1	1	1
D.2.11	Acoustic pressure	N/A	1	1	1			1	1	1
D.2.12	Vibration	N/A	1	1	1			1	1	1
D.2.13	Magnetic fields (eg. MRI)	N/A	2	2	4		F. User Instructions	2	1	2
D.3 Biological hazards and contributory factors										
D.3.1	Bio-contamination	Instrument & Leads	2	2	4	Easy to clean – cleaning instructions in user instructions	F. User Instructions	2	1	2
D.3.2	Bio-incompatibility	Electrodes	1	1	1		Manufacturers Data	1	1	1
D.3.3	Incorrect formulation (chemical composition)	N/A	1	1	1			1	1	1
D.3.4	Toxicity	N/A	1	1	1		Manufacturers Data	1	1	1

T:\Kevin Rush\CE Files\CE files current\Microstim\CE File_Microstim_2004_PP'd\E. Risk Analysis\EN 14971 Annexe D\E. Risk Management
BS14971 Annex D_Microstim.doc Page 1

19/03/04

Severity of hazard: 1 = Insignificant, 2 = Tolerable, 3 = Critical, 4 = Intolerable. Probability of event: 1 = Improbable, 2 = Occasional, 3 = Likely, 4 = Highly likely.

Risk calculated as severity of hazard x probability of event, 1 – 16. Further explanation of risk management policy – see Risk Management Policy & Definitions.

D.3.5	Allergenicity	N/A	1	1	1		Manufacturers Data	1	1	1
D.3.6	Mutagenicity	N/A	1	1	1		Manufacturers Data	1	1	1
D.3.7	Oncogenicity	N/A	1	1	1		Manufacturers Data	1	1	1
D.3.8	Carcinogenicity	N/A	1	1	1		Manufacturers Data	1	1	1
D.3.9	Re and/or cross infection	Instrument & Leads	1	1	1		Manufacturers Data	1	1	1
D.3.10	Pyrogenicity	N/A	1	1	1		Manufacturers Data	1	1	1
D.3.11	Inability to maintain hygienic standards	Instrument & Leads	2	2	4	Easy to clean – cleaning instructions in user manual	F. User Instructions	2	1	2
D.3.12	Degradation	Instrument & Leads	2	1	2	Care Instructions given in user manual	F. User Instructions	1	1	1
D.4 Environmental hazards and contributory factors										
D.4.1	Electromagnetic fields	N/A	1	1	1			1	1	1
D.4.2	Susceptibility to electromagnetic interference	N/A	1	1	1			1	1	1
D.4.3	Emissions of electromagnetic interference	N/A	1	1	1			1	1	1
D.4.4	Inadequate supply of power	N/A	1	1	1			1	1	1
D.4.5	Inadequate supply of coolant	N/A	1	1	1			1	1	1
D.4.6	Storage / operation outside prescribed environmental conditions	Instrument	2	1	2	No specific environmental storage / operating conditions imposed		2	1	2
D.4.7	Incompatibility with other devices with which the product is intended to be used	ECG Pads	1	1	1	Standard sized ECG press studs on leads		1	1	1

T:\Kevin Rush\CE Files\CE files current\Microstim\CE File_Microstim_2004_PP'd\E. Risk Analysis\EN 14971 Annexe D\E. Risk Management BS14971 Annex D_Microstim.doc Page 2

19/03/04

Severity of hazard: 1 = Insignificant, 2 = Tolerable, 3 = Critical, 4 = Intolerable. Probability of event: 1 = Improbable, 2 = Occasional, 3 = Likely, 4 = Highly likely.

Risk calculated as severity of hazard x probability of event, 1 – 16. Further explanation of risk management policy – see Risk Management Policy & Definitions.

D.4.8	Accidental mechanical damage	Instrument & Leads	2	1	2	Relatively robust material used. If damaged, user to assess level of damage – return for repair or spares parts available	F. User Instructions	1	1	1
D.4.9	Contamination due to waste products and/or device disposal	Instrument & Leads	1	1	1	No special disposal required. Alkaline batteries used		1	1	1
D.5 Hazards resulting from incorrect output of energy and substances										
D.5.1	Electricity	N/A	1	1	1			1	1	1
D.5.2	Radiation	N/A	1	1	1			1	1	1
D.5.3	Volume	N/A	1	1	1			1	1	1
D.5.4	Pressure	N/A	1	1	1			1	1	1
D.5.5	Supply of medical gases	N/A	1	1	1			1	1	1
D.5.6	Supply of anaesthetic agents	N/A	1	1	1			1	1	1
D.6 Hazards related to the use of the medical device and contributory factors										
D.6.1	Inadequate labelling	Label	2	1	2	Product easy to use. Label	F. Labels	1	1	1
D.6.2	Inadequate operating instructions	User Instructions	2	1	2	Product easy to use. User manual	F. User Instructions	1	1	1
D.6.3	Inadequate specification of accessories	N/A	1	1	1			1	1	1
D.6.4	Inadequate specification of pre-use checks	User Instructions / Labels	2	2	4	Product easy to use. User manual / labels	F. User Instructions F. Labels	2	1	2
D.6.5	Over-complicated operating instructions	User Instructions	2	1	2	Product easy to use. User manual	F. User Instructions	1	1	1

D.6.6	Inadequate specification of service and maintenance	N/A	1	1	1	Instrument set up at manufacturer. Does not require regular service. Maintenance only required if performance degrades		1	1	1
D.6.7	Use by unskilled / untrained personnel	Instrument & Leads	3	2	6	Product easy to use. Warning in user manual – not for use by untrained personnel	F. User Instructions	2	1	2
D.6.8	Reasonable foreseeable misuse	Instrument & Leads	3	2	6	Product easy to use. Warning in user manual – not for use by untrained personnel	F. User Instructions	2	1	2
D.6.9	Insufficient warning of side effects	N/A	1	1	1			1	1	1
D.6.10	Inadequate warnings of hazards likely with re-use of single use devices	N/A	1	1	1			1	1	1
D.6.11	Incorrect measurement and other metrological aspects	N/A	1	1	1			1	1	1
D.6.12	Misrepresentation of results	Instrument & Leads	1	1	1	Only for use by qualified & trained personnel – warning & in-depth explanation in user manual	F. User Instructions	1	1	1
D.6.13	Incompatibility with consumables / accessories / other devices	Battery & Leads	1	1	1	Standard battery used. Leads colour coded – spares available	F. User Instructions	1	1	1

T:\Kevin Rush\CE Files\CE files current\Microstim\CE File_Microstim_2004_PP'd\E. Risk Analysis\EN 14971 Annexe D\E. Risk Management BS14971 Annex D_Microstim.doc Page 4

19/03/04

Severity of hazard: 1 = Insignificant, 2 = Tolerable, 3 = Critical, 4 = Intolerable. Probability of event: 1 = Improbable, 2 = Occasional, 3 = Likely, 4 = Highly likely.

Risk calculated as severity of hazard x probability of event, 1 – 16. Further explanation of risk management policy – see Risk Management Policy & Definitions.

D.6.14	Sharp edges or points	Instrument & Leads	2	2	4	Only if damaged, user to assess level of damage – return for repair or spares parts available	F. User Instructions	1	1	1
D.7 Inappropriate, inadequate or overcomplicated user interface (man/machine communication)										
D.7.1	Mistakes & judgement errors	N/A	1	1	1	No complex user interface		1	1	1
D.7.2	Lapses and cognitive recall errors	N/A	1	1	1	No complex user interface		1	1	1
D.7.3	Slips & blunders (mental or physical)	N/A	1	1	1	No complex user interface		1	1	1
D.7.4	Violation or abbreviation of instructions, procedures etc	N/A	1	1	1	No complex user interface		1	1	1
D.7.5	Complex or confusing control system	N/A	1	1	1	No complex user interface		1	1	1
D.7.6	Ambiguous or unclear device state	N/A	1	1	1	No complex user interface		1	1	1
D.7.7	Ambiguous or unclear presentation of settings, measurement, or other information	N/A	1	1	1	No complex user interface		1	1	1
D.7.8	Misrepresentation of results	N/A	1	1	1	No results presented		1	1	1
D.7.9	Insufficient visibility, audibility or tactility	N/A	1	1	1	No complex user interface		1	1	1
D.7.10	Poor mapping of controls to action or of displayed information to actual state	N/A	1	1	1	No complex user interface		1	1	1

T:\Kevin Rush\CE Files\CE files current\Microstim\CE File_Microstim_2004_PP'd\E. Risk Analysis\EN 14971 Annexe D\E. Risk Management BS14971 Annex D_Microstim.doc Page 5

19/03/04

Severity of hazard: 1 = Insignificant, 2 = Tolerable, 3 = Critical, 4 = Intolerable. Probability of event: 1 = Improbable, 2 = Occasional, 3 = Likely, 4 = Highly likely.

Risk calculated as severity of hazard x probability of event, 1 – 16. Further explanation of risk management policy – see Risk Management Policy & Definitions.

D.7.11	Controversial modes or mappings as compared to existing equipment	N/A	1	1	1	No complex user interface		1	1	1
D.8 Hazards arising from function failure, maintenance and ageing and contributory factors										
D.8.1	Erroneous data transfer	N/A	1	1	1			1	1	1
D.8.2	Lack of, or inadequate specification for maintenance including post maintenance functional tests	N/A	1	1	1			1	1	1
D.8.4	Inadequate maintenance	N/A	1	1	1	Instrument set up at manufacturer. Does not require regular service. Maintenance only required if performance degrades		1	1	1
D.8.5	Lack of adequate determination of end of device life	Instrument	1	1	1	User decision based on damage or cost effectiveness of repair		1	1	1
D.8.6	Loss of electrical integrity	Instrument	2	2	4	Failure of any electrical component will result in loss of audio or visual output indication		1	1	1
D.8.7	Loss of mechanical integrity	Instrument	1	1	1	User decision based on damage or cost effectiveness of repair		1	1	1
D.8.8	Inadequate packaging (contamination and / or deterioration of the device)	Instrument & Leads	2	2	4	Boxed in protective plastic case. Bubble wrap / cardboard box through mail / courier	M. Packaging	1	1	1

T:\Kevin Rush\CE Files\CE files current\Microstim\CE File_Microstim_2004_PP'd\E. Risk Analysis\EN 14971 Annexe D\E. Risk Management BS14971 Annex D_Microstim.doc Page 6

19/03/04

Severity of hazard: 1 = Insignificant, 2 = Tolerable, 3 = Critical, 4 = Intolerable. Probability of event: 1 = Improbable, 2 = Occasional, 3 = Likely, 4 = Highly likely.

Risk calculated as severity of hazard x probability of event, 1 – 16. Further explanation of risk management policy – see Risk Management Policy & Definitions.

D.8.9	Re-use and/or improper re-use	Instrument & Leads	1	1	1	User decision on suitability for next use based on pre-use check & tolerable damage to casing		1	1	1
D.8.10	Deterioration in function (gradual occlusion of fluid / gas path or change in resistance to flow, electrical conductivity) as a result of repeated use	Instrument & Leads	1	1	1	Degradation of electrical components will eventually result in loss of audio or visual output indication		1	1	1