

Steve Hardaker <viamed.steve.hardaker@gmail.com>

Re: Ceratherm 600-3 Mains Failure Alarm

1 message

Andrew Bullock <andrew.bullock@mbht.nhs.uk>

To: "steve.hardaker@viamed.co.uk" <steve.hardaker@viamed.co.uk>

22 November 2022 at 08:25

Morning Steve,

Yes it worked well. The pins on the cap needed realigning slightly to fit the board but nothing drastic. So far so good.

Thanks

Andrew

Andrew Bullock

Senior Medical Engineer

University Hospitals of Morecambe Bay NHS Foundation Trust

Phone: 01229 614115

Email: andrew.bullock@mbht.nhs.uk

Click HERE to report any NON-URGENT faults with medical equipment.

From: Main Account <viamedinbox@gmail.com> on behalf of Steve Hardaker <office@viamed.co.uk>

Sent: 17 November 2022 17:19

To: Andrew Bullock <andrew.bullock@mbht.nhs.uk> Subject: Re: Ceratherm 600-3 Mains Failure Alarm

Hi Andrew,

A couple of months ago we discussed replacing the Super Cap on a radiant warmer, I just wanted to check with you that this worked as expected?

Any feedback would be appreciated and will help other customers should this issue arise again, thanks.

Best Regards,

Steve Hardaker

Technical Support Manager

Viamed Ltd.

Please note: Viamed is now enacting hybrid working to allow sales and admin staff to work remotely and in the office. Telephone calls to the main office will be answered, but please continue to use email where possible.

Website: www.viamed.co.uk

email: steve.hardaker@viamed.co.uk

Tel: +44 (0)1535 634542 Fax: +44 (0)1535 635582

Sign up to our mailing list, go to www.viamed.co.uk/sign-up

Twitter: twitter.com/ViamedLtd

Facebook: facebook.com/ViamedLtdUK Company registered in England, No. 1291765.

E-mail Disclaimer

The information in this email is confidential and may be legally privileged.

It is intended solely for the addressee. Access to this email by anyone else is unauthorized.

If you are not the intended recipient, any disclosure, copying, distribution or any action taken or omitted to be taken in reliance on it, is prohibited and may be unlawful.

If you feel you have received this email in error or would like to be removed from our mailing list, please reply with delete in the subject line or call +44