

Compatible SpO₂ Finger Probes Cross Reference

Original Equipment Manufacturer	Cross Reference of Model and Part Numbers	Length *Metres	Part Number
Baxter	Explorer Monarch	0.9	P856
Baxter	600333-003 ASAT	3.7	P860RA
Biochem (BCI)	3044 3300 3301 3302	0.9	P861RA
Bruker	Bruker	0.9	P861RA
Corometrics	556	0.9	P856RA
Criticare	511SD/934SD 1100 1100-3 501 503 504 507 508 581 Poet Plus	3.0	P862RA
Criticare	511SD/934SD 1100 1100-3 501 503 504 507 508 581 Poet Plus (Requires P968E8 or P969E8)	0.9	P869RA
Critikon	9700 Series	0.9	P856RA
Critikon	9700 Series	3.7	P858RA
Critikon	9084 8700 Series	0.9	P865RA
Critikon	8997 Dinamap Oxytrak	0.9	P877RA
Datascope	0600-00-0026-02 Passport (Requires P963E10)	0.9	P863RA
Datascope	0600-00-0026-01 Accusat Accutor	3.7	P864RA
Datex	Fits Most Datex Models	3.7	P872RA
Datex	Fits Most Datex Models (Requires P973E10)	0.9	P873RA
Dräger	Vitalert 1000 Narkomed Cicero	0.9	P856RA
Elmed	6006194 PSO1000	0.9	P861RA
Entri Medical	Propaq 102 104 106	0.9	P856RA
Fukuda Denshi	HG302 Module	0.9	P856RA
Gambro	BP806	0.9	P856RA
Invivio	4500 Plus 2	0.9	P856P
Invivio	9383 Scout ES 4500 Scout Scout + Scout +1 Scout +2 Scout +3 Scout +4	3.0	P878K
Ivy Medical	405A	0.9	P856RA
Kontron	7840	3.7	P886RA
Marquette	Eagle 3100	0.9	P856RA
Marquette	Tram System Modular System	3.7	P867RA
Medical Data	2800-50 OP11A OP30 Escort 100 Escort 11	0.9	P856RA
Mennen	800-060-730 Horizon XL (Nellcor)	0.9	P856RA
Mennen	800-060-540 Horizon XL (Ohmeda)	3.7	P867RA
Nellcor	C-13-200M N-200 DS100A (Compatible with Nellcor derivatives. Requires P956E4 or P956E8)	0.9	P856RA
Nellcor	N6000 N180E N180	3.7	P858RA
Nonin	8500 8604K (Requires P971E3 or P971E10)	0.9	P871RA
Novamatrix	8660 500 Series	3.0	P875RA
Novamatrix	8774 8776 505 511 515 520 7000 Oxypleth	3.0	P876RA
Ohmeda	3700 4700 Oxystat Plus 4700 380-1000-042 6051-0000-049 Oxystat Modular CD Anae	3.7	P867RA
Pacetek	520	3.7	P859RA
Propaq	Propaq	0.9	P856RA
S & W Vickers	4621643 4700708 DSNT	0.9	P856RA
Sensormedics	766721 767367 Oxyshuttle OMNI-SAT	0.9	P856RA
Siemens	1281	0.9	P856RA
Simed	S-100c	3.7	P860RA
Spacelabs	015-0150-00 90465/66/67/89 90651A-08 690-0003-00 (Nellcor). Requires P926E10	0.9	P856I
Spacelabs	90651A-08 90419-01 90485 90489 90467 90466 90465 90351-D/G	3.0	P857RA
Spacelabs	369-083-001 90419 90501 90651A-03 (Novamatrix)	3.0	P875RA
	Extension Cables		
Critikon/Sensormedics	9084 8700 Series **(P865RA)	3.0	P965E10
Criticare	518LD 1100 1100-3 Poet Plus 503 504 507 508 581 (DB9 to DB9) **P869RA	2.4	P968E8
Criticare	518DD 503 Spot (Lemo to DB9) **P869RA	2.4	P969E8
Datascope	012-00-0516-12 Passport (Switchcraft to DB9) **P863RA	3.0	P963E10
Datex	Nicolet to DB9 **P873RA	3.0	P973E10
Nellcor	EC-4 N180 N180E N6000 (DB9 to DB9) **P856RA	1.2	P956E4
Nellcor	EC-4 N180 N180E N6000 (DB9 to DB9) **P856RA	2.4	P956E8
Nonin	85001 (DB9 to DB9) **P871RA	0.9	P971E3
Nonin	8604x-10 (DB9 to DB9) **P871RA	3.0	P971E10
Spacelabs	175-0646-00 90465/66/67/89 90651A-08 90351-D/G (7 Pin Hypertronic to 12 Pin) **P876RA	0.6	P925E2
Spacelabs	175-0706-00 90465/66/67/89 90651A-08 690-0003-00 90351-DG (7 Pin Hypertronic to 12 Pin) **P876RA	3.0	P926E10

VIAMED

15 Station Road, Cross Hills, Keighley,
West Yorkshire Bd20 7DT UK
Tel: 01535 634542 Fax: 01535 635582
www.Viamed.co.uk Email: info@viamed.co.uk

CE0086

For any manufacturer not mentioned please contact us.
All Manufacturers' trade marks and registered names recog

* Lengths are approximate
** Indicates suitable for use with stated finger probe.

Distmci97 20/7/99

7 May 1998

Part number	Compatible manufacturer	Dot colour	Pad seat colour	Pad colour	Cable			Comments
					Colour & material	Length		
						Feet (ft)	Metres (m)	
P856RA	Nellcor	Blue	Blue	White	White PVC	3	0.90	
P857RA	Spacelabs	White	Blue	White	White PVC	10	3.00	
P858RA	Nellcor	Blue	Blue	White	White PVC	12	3.65	
P859RA	Pace Tech	White	Blue	White	White PVC	12	3.65	
P860RA	Baxter	White	Blue	White	White PVC	12	3.65	
P861RA	BCI	Orange	Orange	White	White PVC	3	0.90	
P862RA	Criticare	Purple	Purple	White	White PVC	10	3.00	
P863RA	Datascope	White	Green	Black	White PVC	3	0.90	
P864RA	Datascope	White	Green	Black	White PVC	12	3.65	
P865RA	Sensormedics	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P867RA	Ohmeda	Green	Green	Black	White PVC	12	3.65	
P869RA	Criticare	Purple	Purple	White	White PVC	3	0.90	
P870RA	Critikon	White	White	White	White PVC	3	0.90	Discontinued, replaced by P865RA.
P871RA	Nonin	White	Tan	White	White PVC	3	0.90	
P872RA	Datex	Tan	Tan	White	White PVC	12	3.65	
P873RA	Datex	Tan	Tan	White	White PVC	3	0.90	
P875RA	Novamatrix	Yellow	Yellow	Black	White PVC	10	3.00	
P876RA	Novamatrix	Yellow	Yellow	Black	White PVC	10	3.00	
P877RA	Critikon	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P878RA	Invivo	White	Blue	White	White PVC	10	3.00	
P886RA	Kontron	?	?	White	White PVC	12	3.65	Not released yet due to unavailability of connectors.
P887RA	Spacelabs	Yellow	Yellow	Black	White PVC	12	3.65	Unissued part number.

7 May 1998

Part number	Connector				Supplied by MCI	Derivation	Comments
	Type	Number of contacts	Gender	Comments			
P856RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P857RA	Hypertronics	12	M	Pins 5,6,7,9 & 10 missing.		P858	
P858RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P859RA	Switchcraft	5	F			P858	
P860RA	Hypertronics	9	M	Pins 2 & 5 missing?		P858	
P861RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P862RA	Lemo	5	M		✓	Original	
P863RA	D	9	M	Pin 8 missing.		P867	
P864RA	DIN	8	M			P867	
P865RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P867RA	Hypertronics	9	M	Pins 3 & 6 missing.	✓	Original	
P869RA	D	9	M			P862	
P870RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P871RA	D (Amp)	9	M			P873?	
P872RA	DPK1	10	F	Pins 2,3 & 0 missing.	✓	Original	No connector supplied
P873RA	D	9	M	Pins 2 & 8 missing.	✓	Original	
P875RA	Hypertronics	7	M	Pin 7 missing.		P876	Wire cross over
P876RA	Hypertronics	7	M	Pin 7 missing.	✓	Original	
P877RA	D	9	M	Pins 7, 8 & 9 missing		Kit	Sensormedics LED's & cable.
P878RA	Lemo	7	M			P858	Query cable length
P886RA					✓	Original	
P887RA	Hypertronics	12	M	Pins 5,6,7,8 & 12 missing.		P876	

7 May 1998

Part number	Connection		Extension cable part number
	Monitor	Extension cable	
P856RA	✓	✓	Nellcor 958E8 or P958E4. Spacelabs P926E10
P857RA	✓		
P858RA	✓	✓	
P859RA	✓		
P860RA	✓		
P861RA		✓	
P862RA	✓		
P863RA		✓	P963E10
P864RA	✓		
P865RA		✓	
P867RA	✓		
P869RA		✓	P968E8 or P969E8
P870RA		✓	
P871RA		✓	P971E3 or P971E10
P872RA	✓		
P873RA		✓	P973E10
P875RA	✓	✓	Spacelabs P925E2
P876RA	✓	✓	Spacelabs P925E2
P877RA		✓	
P878RA	✓		
P886RA	✓		
P887RA	✓		

MCI. NO.	MANUFACTUR	Identification	LENGTH_FT
P856RA	Baxter	Explorer Monarch	3.00
P856RA	Corometrics	556.00	3.00
P856RA	Critikon	9700 series	3.00
P856RA	Draegar	Vitalert 1000	3.00
		Narkomed Cicero	
P856RA	3.00	BP806	3.00
P856RA	Entri Medical	Propaq 102 104 106	3.00
P856RA	Fukuda Denshi	HG302 Module	3.00
P856RA	Gambro	BP806	3.00
P856RA	Hewlett Packard	Miracle Adult system	3.00
		Merlin M1020A	
		Neonatal	
P856RA	Invivo	4500+2	3.00
P856RA	Ivy Medical	405A	3.00
P856RA	Marquette	Eagle 3100	3.00
P856RA	Medical Data	2800-50 OP11A	3.00
		OP30 Escort 100	
		Escort II	
P856RA	Mennen	800-060-730 Horizon	3.00
		XL	
P856RA	Nellcor	C-13-200M N-200	3.00
		DS100A	
P856RA	3.00	1281.00	3.00
P856RA	Propaq	Propaq	3.00
P856RA	S & W Vickers	4621643 4700708	3.00
		DSNT	
P856RA	Siemens	1281.00	3.00
P856RA	Spacelabs	015-0150-00 90489	
		90465 90466 90467	
		90651A-08	
	3.00		
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	3.00	N6000,N180E,N180	
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	Nellcor	N6000 N180E N180	0.00
P860RA	Baxter	600333-003,ASAT	
	11.00		

P859RA	Pacetek	520.00	0.00
P860RA	Baxter	600333-003 ASAT	
P861RA		0.00 3044,3300,3301,3302	
	P860RA Simed	S-100e	
		0.00	
P861RA	Biochem (BCI)	3044 3300 3301 3302	3.00
	3044.00 Elmed	6006194,PSO1000	
		3.00	
P861RA	Bruker	Bruker	3.00
P861RA	Elmed	6006194 PSO1000	
		3.00	
P862RA	Criticare	511SD/934 1100 1100-3 501 503 504 507 508 581 Poet Plus	12.00
P863RA	Datascope	0600-00-0026-02 Passport	
		12.00	
P864RA	Datascope	0600-00-0026-01 Accustat	12.00
P865RA	Sensor medics	766721 767367 Oxyshuttle Omnistat	
P867RA	Marquette	Tram system Modular Sys	12.00
P867RA	Mennen	800-060-540,Horizon XL	
		3.00	
P867RA	Marquette	Tram system Modular Sys	12.00
P867RA	Mennen	800-060-540 Horizon XL	
	4700.00	12.00	
P867RA	Ohmeda	380-1000-042 6051-0000-049 Oxystat Modular CD Anaes 3700 4700	
,Oxystat Plus			
	4700.00	12.00	
P870RA	Critikon	9084,8700 series	
		12.00	
P871RA		12.00 8604K,Onyx	0.00

P870RA	Critikon	9084 8700 series	
,CH-2S,SATLITE,ULTS00,OSCAR II,M-ESTPR,CH-RS,ULTSVO,CARDIO CAP II		12.00	
P871RA	Nonin	8604K Onyx	0.00
P872RA	Datex	878579-1 SAS-F4 CG-2CS CH-S CH-RS ULTSIO ULTSVI AS/3 OSP-200,	
CH-2S SATLITE ULTS00 OSCAR II M-ESTPR CH-RS ULTSVO CARDIO CAP II		12.00 369-083-001,90419,90501,90651A-03	
P875RA	Novametrix	8660 500 series	0.00
P875RA	Spacelabs	369-083-001 90419 90501 90651A-03	12.00
P876RA		0.00 015-0130-00,90651A-08,90419-01,90489,90467,90466,90465,90351-D/G	
P876RA	Novametrix	8774 8776 505 511 515 520 7000 Oxypleth	12.00
P876RA	Spacelabs	015-0130-00 90651A-08 90419-01 90489 90467 90466 90465 90351-D/G	12.00
P877RA	Critikon	8997,Dinamap Oxytrak	0.00
P878RA		12.00 9383,Scout ES,Scout +4,4500Scout,Scout +,Scout+1,Scout+2,Scout +3	
P877RA	Critikon	8997 Dinamap Oxytrak	0.00

P878RA	Invivo	9383 Scout ES Scout
		+4 4500Scout Scout
		+ Scout+1 Scout+2
		Scout +3
		12.00
		12.00
P886RA	Kontron	78400 7840
		12.00

MCI. NO.	MANUFACTUR	Identification	LENGTH_FT
P856RA	Baxter	Explorer Monarch	3.00
P856RA	Corometrics	556.00	3.00
P856RA	Critikon	9700 series	3.00
P856RA	Draegar	Vitalert 1000	3.00
		Narkomed Cicero	
P856RA	3.00	BP806	3.00
P856RA	Entri Medical	Propaq 102 104 106	3.00
P856RA	Fukuda Denshi	HG302 Module	3.00
P856RA	Gambro	BP806	3.00
P856RA	Hewlett Packard	Miracle Adult system	3.00
		Merlin M1020A	
		Neonatal	
P856RA	Invivo	4500+2	3.00
P856RA	Ivy Medical	405A	3.00
P856RA	Marquette	Eagle 3100	3.00
P856RA	Medical Data	2800-50 OP11A	3.00
		OP30 Escort 100	
		Escort II	
P856RA	Mennen	800-060-730 Horizon	3.00
		XL	
P856RA	Nellcor	C-13-200M N-200	3.00
		DS100A	
P856RA	3.00	1281.00	3.00
P856RA	Propaq	Propaq	3.00
P856RA	S & W Vickers	4621643 4700708	3.00
		DSNT	
P856RA	Siemens	1281.00	3.00
P856RA	Spacelabs	015-0150-00 90489	
		90465 90466 90467	
		90651A-08	
	3.00		
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	3.00	N6000,N180E,N180	
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	Nellcor	N6000 N180E N180	0.00
P860RA	Baxter	600333-003,ASAT	
	11.00		
P859RA	Pacetek	520.00	0.00
P860RA	Baxter	600333-003 ASAT	
P861RA	0.00	3044,3300,3301,3302	
P860RA	Simed	S-100e	
	0.00		
P861RA	Biochem (BCI)	3044 3300 3301	3.00
		3302	
3044.00	Elmed	6006194,PSO1000	
	3.00		


P861RA	Bruker	Bruker	3.00
P861RA	Elmed	6006194 PSO1000	
	3.00		
P862RA	Criticare	511SD/934 1100	12.00
		1100-3 501 503 504	
		507 508 581 Poet	
		Plus	
P863RA	Datascope	0600-00-0026-02	
		Passport	
	12.00		
P864RA	Datascope	0600-00-0026-01	12.00
		Accustat	
P865RA	Sensor medics	766721 767367	
		Oxyshuttle Omnistat	
P867RA	Marquette	Tram system Modular	12.00
		Sys	
P867RA	Mennen	800-060-540,Horizon	
		XL	
	3.00		
P867RA	Marquette	Tram system Modular	12.00
		Sys	
P867RA	Mennen	800-060-540 Horizon	
		XL	
4700.00	12.00		
P867RA	Ohmeda	380-1000-042 6051-	
		0000-049 Oxystat	
		Modular CD Anaes	
		3700 4700	
,Oxystat Plus			
4700.00	12.00		
P870RA	Critikon	9084,8700 series	
	12.00		
P871RA	12.00	8604K,Onyx	0.00
P870RA	Critikon	9084 8700 series	
,CH-	12.00		
2S,SATLITE,ULTS0			
0,OSCAR II,M-			
ESTPR,CH-			
RS,ULTSVO,CARDI			
O CAP II			
P871RA	Nonin	8604K Onyx	0.00
P872RA	Datex	878579-1 SAS-F4	
		CG-2CS CH-S CH-	
		RS ULTSIO	
		ULTSVI AS/3 OSP-	
		200,	

CH-2S SATLITE ULTS00 OSCAR II M-ESTPR CH-RS ULTSVO CARDIO CAP II	12.00	369-083- 001,90419,90501,906 51A-03
P875RA	Novamatrix 0.00	8660 500 series
P875RA	Spacelabs 12.00	369-083-001 90419 90501 90651A-03
P876RA	0.00	015-0130-00,90651A- 08,90419- 01,90489,90467,9046 6,90465,90351-D/G
P876RA	Novamatrix 12.00	8774 8776 505 511 515 520 7000 Oxypleth
P876RA	Spacelabs	015-0130-00 90651A- 08 90419-01 90489 90467 90466 90465 90351-D/G
P877RA	12.00 Critikon	8997,Dinamap Oxytrak
P878RA	0.00 12.00	9383,Scout ES,Scout +4,4500Scout,Scout +,Scout+1,Scout+2,Sc out +3
P877RA	Critikon 0.00	8997 Dinamap Oxytrak
P878RA	Invivo	9383 Scout ES Scout +4 4500Scout Scout + Scout+1 Scout+2 Scout +3
P886RA	12.00 12.00 Kontron 12.00	78400 7840

**Description Chart To Aid The Identification & Selection
Of Viamed SpO₂ Probes**

Part Number	Main Compatible Manufacturers	Button Colour	Pad Colour	Cable Length		Connector Details				
				Feet (ft)	Metres (m)	Type	Colour Shell/Insert	Number Of Contact Positions	Gender	Pins/sockets Missing
P856RA	Nellcor	Blue	White	3	0.90	D	White	9	M	Pins 4 & 8
P857RA	Spacelabs	White	White	10	3.00	Hypertronics	Grey	12	M	Pins 5,6,7,9 & 10
P858RA	Nellcor	Blue	White	12	3.65	D	White	9	M	Pins 4 & 8
P859RA	Pace Tech	White	White	12	3.65	Switchcraft	Black	5	F	
P860RA	Baxter, Simed	White	White	12	3.65	Hypertronics	Black	9	M	Pins 2 & 5
P861RA	BCI	Orange	White	3	0.90	D	White	9	M	Pins 4 & 8
P862RA	Criticare	Purple	White	10	3.00	Lemo	Silver	5	M	
P863RA	Datascope	White	Black	3	0.90	D	White	9	M	Pin 8
P864RA	Datascope	White	Black	12	3.65	DIN	Black	8	M	
P865RA	Sensormedics	White	White	3	0.90	Compu-shield	Black	6	M	
P866RA	Baxter, Simed	White	White	12	3.65	Hypertronics	Black	9	M	Pins 2 & 5
P867RA	Ohmeda	Green	Black	12	3.65	Hypertronics	Black	9	M	Pins 3 & 6
P869RA	Criticare	Purple	White	3	0.90	D	White	9	M	Pins 3,7 & 8
P871RA	Nonin	White	White	3	0.90	D	White	9	M	Pins 4 & 8
P872RA	Datex	Tan	White	12	3.65	DPK1	Grey/Grey	10	F	Sockets 2,3 & 0
P873RA	Datex	Tan	White	3	0.90	D	White	9	M	Pins 2 & 8
P875RA	Novamatrix	Yellow	Black	10	3.00	Hypertronics	Black	7	M	Pin 7
P876RA	Novamatrix	Yellow	Black	10	3.00	Hypertronics	Black	7	M	Pin 7
P877RA	Critikon	White	White	3	0.90	D	White	9	M	Pins 7, 8 & 9
P878RA	Invivo	White	White	10	3.00	Lemo	Silver	7	M	
P886RA	Kontron	White	White	12	3.65	KPK1	Grey/Maroon	12	M	Pins 1,2,4,6,9 & 11
P887RA	Spacelabs	Yellow	Black	12	3.65	Hypertronics	Grey	12	M	Pins 5,6,7,8 & 12
P888RA	S & W	Green	Black	12	3.65	SPK1	Grey/Purple	10	F	Sockets 3 & 6
P889RA	Criticare	Purple	White	10	3.00	D	White	9	M	Pins 3,7 & 8
P890RA	Critikon	White	White	10	3.00	Thorkom	Grey	12	F	Sockets 5,8,9 & 10
P891RA	Critikon	White	White	10	3.00	Lemo	Silver	7	M	

14 June 1999

P856RA			P865RA			P878RA		
P856YA		T	P865YA		T	P878YA		T
P857RA			P867RA			P886RA		
P857YA		T	P867YA		T	P886YA		S
P858RA			P869RA			P887RA		
P858YA		T	P869YA		T	P887YA		T
P859RA			P871RA			P888RA		
P859YA		S	P871YA		T	P888YA		S
P860RA			P872RA			P889RA		
P860YA		T	P872YA		S	P889YA		T
P861RA			P873RA			P890RA		
P861YA		T	P873YA		T	P890YA		T
P862RA			P875RA			 Distributed By: Manufactured by: VIAMED Ltd 15 Station Rd, Crosshills, Kelghley, West Yorkshire, BD20 7DT, U.K. Tel : +44 (0)1535 634542, Fax : +44 (0)1535 635582 Web: www.viamed.co.uk, Email: info@viamed.co.uk Quality standards accreditations BS EN ISO 9002 & BS EN 46002 BS28343 BS28344		
P862YA		T	P875YA		T			
P863RA			P876RA					
P863YA		T	P876YA		T			
P864RA			P877RA					
P864YA		T	P877YA		T			

v: 'T' Top View, 'S' Side View, ● Pin/Socket Loaded, ○ Pin/Socket Unloaded Drawings not to scale.

Sp O₂ Finger Probe Details

Model	Compatible manufacturer	Length M (ft)	Connector type	Model	Notes Document Name VM3/COP/.....	Probe details (colours)			Labels	
						Pad support	Button	Pads	CE	Db
854RA	Medlab/Comdek	3 (10)	Redel M-7	854RA	VM3/COP/32.01	White	White	White	White	Orange
855RA	BCI	3 (10)	Redel M-7	855RA	VM3/COP/32.02	White	White	White	White	Or/Br
856RA	Nellcor	0.9 (3)	D-M-9	856RA	VM3/COP/32.03	Blue	Blue	White	White	Blue
857RA	Spacelabs	3 (10)	Hypertronics M-12	857RA	VM3/COP/32.04	White	White	White	White	Orange
858RA	Nellcor	3.65 (12)	D-M-9	858RA	VM3/COP/32.05	Blue	Blue	White	White	Blue
859RA	Pace Tech	3.65 (12)	Switchcraft F-5	859RA	VM3/COP/32.06	White	White	White	White	Orange
860RA	Simed Baxter	3.65 (12)	Hypertronics M-9	860RA	VM3/COP/32.07	White	White	White	White	Orange
861RA	BCI	0.9 (3)	D-M-9	861RA	VM3/COP/32.08	Orange	Orange	White	White	Or/Br
862RA	Criticare	3 (10)	Lemo M-5	862RA	VM3/COP/32.09	Purple	Purple	White	White	Purple
863RA	Datascope	0.9 (3)	D-M-9	863RA	VM3/COP/32.10	Red	Red	White	White	White
864RA	Datascope	3.65 (12)	DIN M-8	864RA	VM3/COP/32.11	Red	Red	White	White	White
865RA	Sensormedics	0.9 (3)	Compu-shield	865RA	VM3/COP/32.12	White	White	White	White	Orange
866RA	Simed Baxter	3.65 (12)	Hypertronics M-9	866RA	VM3/COP/32.13	White	White	White	White	Orange
867RA	Ohmeda	3.65 (12)	Hypertronics M-9	867RA	No diagram details	Green	Green	Black	White	Green
869RA	Criticare	0.9 (3)	D-M-9	869RA	VM3/COP/32.15	Purple	Purple	White	White	Purple
871RA	Nonin	0.9 (3)	D-M-9	871RA	VM3/COP/32.16	White	White	White	White	Orange
872RA	Datex	3.65 (12)	DPK1 F-10	872RA	VM3/COP/32.17	Tan	Tan	White	White	Brown
873RA	Datex	0.9 (3)	D-M-9	873RA	VM3/COP/32.18	Tan	Tan	White	White	Brown
874RA	Palco	0.9 (3)	Sub-min	874RA	VM3/COP/32.19	White	White	White	White	Orange
875RA	Novamatrix	3 (10)	Hypertronics M-7	875RA	VM3/COP/32.20	Yellow	Yellow	Black	White	Yellow
876RA	Novamatrix	3 (10)	Hypertronics M-7	876RA	VM3/COP/32.21	Yellow	Yellow	Black	White	Yellow
877RA	Critikon	0.9 (3)	D-M-9	877RA	VM3/COP/32.22	White	White	White	White	Orange
878RA	In vivo	3 (10)	Lemo M-7	878RA	VM3/COP/32.23	White	White	White	White	Orange
879RA	Artema	3.65 (12)	SPK1-A M-12	879RA	VM3/COP/32.24	White	White	White	White	Orange
885RA	Novamatrix	0.9 (3)	D-M-9	885RA	VM3/COP/32.25	Yellow	Yellow	Black	White	Yellow
886RA	Kontron	3.65 (12)	SPK1-A M-12	886RA	VM3/COP/32.26	White	White	White	White	Orange
887RA	Spacelabs	3 (10)	Hypertronics M-12	887RA	VM3/COP/32.27	White	White	Black	White	Orange
888RA	S&W	3.65 (12)	SPK1-O F-10	888RA	VM3/COP/32.28	White	White	Black	White	Orange
889RA	Criticare	3 (10)	D-M-9	889RA	VM3/COP/32.29	Purple	Purple	White	White	Purple
890RA	Critikon & Sensormedics	3 (10)	Thorkom F-12	890RA	VM3/COP/32.30	White	White	White	White	Orange
891RA	Critikon	3 (10)	Lemo M-7	891RA	VM3/COP/32.31	White	White	White	White	Orange
892RA	Marquette	3 (10)	Marquette M-11	892RA	VM3/COP/32.32	White	White	White	White	Orange
895RA	Hewlett Packard	3 (10)	HP M-12	895RA	No confirmation	White	White	White	White	Orange
896RA	S&W	3.65 (12)	SPK1-N F-10	896RA	VM3/COP/32.34	White	White	White	White	Orange

7 May 1998

Part number	Compatible manufacturer	Dot colour	Pad seat colour	Pad colour	Cable			Comments
					Colour & material	Length		
						Feet (ft)	Metres (m)	
P856RA	Nellcor	Blue	Blue	White	White PVC	3	0.90	
P857RA	Spacelabs	White	Blue	White	White PVC	10	3.00	
P858RA	Nellcor	Blue	Blue	White	White PVC	12	3.65	
P859RA	Pace Tech	White	Blue	White	White PVC	12	3.65	
P860RA	Baxter	White	Blue	White	White PVC	12	3.65	
P861RA	BCI	Orange	Orange	White	White PVC	3	0.90	
P862RA	Criticare	Purple	Purple	White	White PVC	10	3.00	
P863RA	Datascope	White	Green	Black	White PVC	3	0.90	
P864RA	Datascope	White	Green	Black	White PVC	12	3.65	
P865RA	Sensormedics	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P867RA	Ohmeda	Green	Green	Black	White PVC	12	3.65	
P869RA	Criticare	Purple	Purple	White	White PVC	3	0.90	
P870RA	Critikon	White	White	White	White PVC	3	0.90	Discontinued, replaced by P865RA.
P871RA	Nonin	White	Tan	White	White PVC	3	0.90	
P872RA	Datex	Tan	Tan	White	White PVC	12	3.65	
P873RA	Datex	Tan	Tan	White	White PVC	3	0.90	
P875RA	Novamatrix	Yellow	Yellow	Black	White PVC	10	3.00	
P876RA	Novamatrix	Yellow	Yellow	Black	White PVC	10	3.00	
P877RA	Critikon	White	White	White	White PVC	3	0.90	Non in stock until we receive new batch of Sensormedics LED's.
P878RA	Invivo	White	Blue	White	White PVC	10	3.00	
P886RA	Kontron	?	?	White	White PVC	12	3.65	Not released yet due to unavailability of connectors.
P887RA	Spacelabs	Yellow	Yellow	Black	White PVC	12	3.65	Unissued part number.

7 May 1998

Part number	Connector				Supplied by MCI	Derivation	Comments
	Type	Number of contacts	Gender	Comments			
P856RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P857RA	Hypertronics	12	M	Pins 5,6,7,9 & 10 missing.		P858	
P858RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P859RA	Switchcraft	5	F			P858	
P860RA	Hypertronics	9	M	Pins 2 & 5 missing?		P858	
P861RA	D	9	M	Pins 4 & 8 missing.	✓	Original	
P862RA	Lemo	5	M		✓	Original	
P863RA	D	9	M	Pin 8 missing.		P867	
P864RA	DIN	8	M			P867	
P865RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P867RA	Hypertronics	9	M	Pins 3 & 6 missing.	✓	Original	
P869RA	D	9	M			P862	
P870RA	Compu-shield	6	M			Kit	Sensormedics LED's & cable.
P871RA	D (Amp)	9	M			P873?	
P872RA	DPK1	10	F	Pins 2,3 & 0 missing.	✓	Original	No connector supplied
P873RA	D	9	M	Pins 2 & 8 missing.	✓	Original	
P875RA	Hypertronics	7	M	Pin 7 missing.		P876	Wire cross over
P876RA	Hypertronics	7	M	Pin 7 missing.	✓	Original	
P877RA	D	9	M	Pins 7, 8 & 9 missing		Kit	Sensormedics LED's & cable.
P878RA	Lemo	7	M			P858	Query cable length
P886RA					✓	Original	
P887RA	Hypertronics	12	M	Pins 5,6,7,8 & 12 missing.		P876	

7 May 1998

Part number	Connection		Extension cable part number
	Monitor	Extension cable	
P856RA	✓	✓	Nellcor 958E8 or P958E4. Spacelabs P926E10
P857RA	✓		
P858RA	✓	✓	
P859RA	✓		
P860RA	✓		
P861RA		✓	
P862RA	✓		
P863RA		✓	P963E10
P864RA	✓		
P865RA		✓	
P867RA	✓		
P869RA		✓	P968E8 or P969E8
P870RA		✓	
P871RA		✓	P971E3 or P971E10
P872RA	✓		
P873RA		✓	P973E10
P875RA	✓	✓	Spacelabs P925E2
P876RA	✓	✓	Spacelabs P925E2
P877RA		✓	
P878RA	✓		
P886RA	✓		
P887RA	✓		

MCI. NO.	MANUFACTUR	Identification	LENGTH_FT
P856RA	Baxter	Explorer Monarch	3.00
P856RA	Corometrics	556.00	3.00
P856RA	Critikon	9700 series	3.00
P856RA	Draegar	Vitalert 1000	3.00
		Narkomed Cicero	
P856RA		3.00 BP806	3.00
P856RA	Entri Medical	Propaq 102 104 106	3.00
P856RA	Fukuda Denshi	HG302 Module	3.00
P856RA	Gambro	BP806	3.00
P856RA	Hewlett Packard	Miracle Adult system	3.00
		Merlin M1020A	
		Neonatal	
P856RA	Invivo	4500+2	3.00
P856RA	Ivy Medical	405A	3.00
P856RA	Marquette	Eagle 3100	3.00
P856RA	Medical Data	2800-50 OP11A	3.00
		OP30 Escort 100	
		Escort II	
P856RA	Mennen	800-060-730 Horizon	3.00
		XL	
P856RA	Nellcor	C-13-200M N-200	3.00
		DS100A	
P856RA		3.00 1281.00	3.00
P856RA	Propaq	Propaq	3.00
P856RA	S & W Vickers	4621643 4700708	3.00
		DSNT	
P856RA	Siemens	1281.00	3.00
P856RA	Spacelabs	015-0150-00 90489	
		90465 90466 90467	
		90651A-08	
		3.00	
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA		3.00 N6000,N180E,N180	
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	Nellcor	N6000 N180E N180	0.00
P860RA	Baxter	600333-003,ASAT	
		11.00	

P859RA	Pacetek	520.00	0.00
P860RA	Baxter	600333-003 ASAT	
P861RA		0.00 3044,3300,3301,3302	
	P860RA Simed	S-100e	
		0.00	
P861RA	Biochem (BCI)	3044 3300 3301 3302	3.00
	3044.00 Elmed	6006194,PSO1000	
		3.00	
P861RA	Bruker	Bruker	3.00
P861RA	Elmed	6006194 PSO1000	
		3.00	
P862RA	Criticare	511SD/934 1100 1100-3 501 503 504 507 508 581 Poet Plus	12.00
P863RA	Datascope	0600-00-0026-02 Passport	
		12.00	
P864RA	Datascope	0600-00-0026-01 Accustat	12.00
P865RA	Sensor medics	766721 767367 Oxyshuttle Omnistat	
P867RA	Marquette	Tram system Modular Sys	12.00
P867RA	Mennen	800-060-540,Horizon XL	
		3.00	
P867RA	Marquette	Tram system Modular Sys	12.00
P867RA	Mennen	800-060-540 Horizon XL	
	4700.00	12.00	
P867RA	Ohmeda	380-1000-042 6051- 0000-049 Oxystat Modular CD Anaes 3700 4700	
,Oxystat Plus			
	4700.00	12.00	
P870RA	Critikon	9084,8700 series	
		12.00	
P871RA		12.00 8604K,Onyx	0.00

P870RA	Critikon	9084 8700 series	
,CH-2S,SATLITE,ULTS00,OSCAR II,M-ESTPR,CH-RS,ULTSVO,CARDIO CAP II		12.00	
P871RA	Nonin	8604K Onyx	0.00
P872RA	Datex	878579-1 SAS-F4 CG-2CS CH-S CH-RS ULTSIO ULTSVI AS/3 OSP-200,	
CH-2S SATLITE ULTS00 OSCAR II M-ESTPR CH-RS ULTSVO CARDIO CAP II		12.00 369-083-001,90419,90501,90651A-03	
P875RA	Novametrix	8660 500 series	0.00
P875RA	Spacelabs	369-083-001 90419 90501 90651A-03	12.00
P876RA		0.00 015-0130-00,90651A-08,90419-01,90489,90467,90466,90465,90351-D/G	
P876RA	Novametrix	8774 8776 505 511 515 520 7000 Oxypleth	12.00
P876RA	Spacelabs	015-0130-00 90651A-08 90419-01 90489 90467 90466 90465 90351-D/G	12.00
P877RA	Critikon	8997,Dinamap Oxytrak	0.00
P878RA		12.00 9383,Scout ES,Scout +4,4500Scout,Scout +,Scout+1,Scout+2,Scout +3	
P877RA	Critikon	8997 Dinamap Oxytrak	0.00

P878RA	Invivo	9383 Scout ES Scout
		+4 4500Scout Scout
		+ Scout+1 Scout+2 Scout +3
		12.00
		12.00
P886RA	Kontron	78400 7840
		12.00












VIAMED SpO₂ FINGER PROBES

Part number	*Compatible oximeter manufacturer	Cable length		Connector		
		Feet	Metres	Type	Number of contacts	Gender
P856RA	Nellcor	3	0.90	D	9	M
P857RA	Spacelabs	10	3.00	Hypertronics	12	M
P858RA	Nellcor	12	3.65	D	9	M
P859RA	Pace Tech	12	3.65	Switchcraft	5	F
P860RA	Simed, Baxter	12	3.65	Hypertronics	9	M
P861RA	BCI	3	0.90	D	9	M
P862RA	Criticare	10	3.00	Lemo	5	M
P863RA	Datascope	3	0.90	D	9	M
P864RA	Datascope	12	3.65	DIN	8	M
P865RA	Sensormedics	3	0.90	Compu-shield	6	M
P866RA	Simed, Baxter	12	3.65	Hypertronics	9	M
P867RA	Ohmeda	12	3.65	Hypertronics	9	M
P869RA	Criticare	3	0.90	D	9	M
P871RA	Nonin	3	0.90	D	9	M
P872RA	Datex	12	3.65	DPK1	10	F
P873RA	Datex	3	0.90	D	9	M
P875RA	Novametrix	10	3.00	Hypertronics	7	M
P876RA	Novametrix	10	3.00	Hypertronics	7	M
P877RA	Critikon	3	0.90	D	9	M
P878RA	Invivo	10	3.00	Lemo	7	M
P886RA	Kontron	12	3.65	DPK1	12	M
P887RA	Spacelabs	10	3.00	Hypertronics	12	M
P888RA	S & W	12	3.65	SPK1	10	F
P889RA	Criticare	10	3.00	D	9	M
P890RA	Critikon/Sensormedics	10	3.00	Thorkom	12	F
P891RA	Critikon	10	3.00	Lemo	7	M




























CE 0086

* Main compatible oximeter manufacturer for guidance purposes only. For full details of other compatible manufacturers, please refer to cross reference list. Cable lengths are approximate.







Finger Probe Chart

			Nihon / Kohden	852RA 0018520	10 ft 3m
			Odam Bruker	853RA 0018530	10 ft 3m
			Medlab/ Comdek	854RA 0018540	10 ft 3m
			BCI	855RA 0018550	10 ft 3m
			Nellcor	856RA 0018560	3 ft 0.9m
			Spacelabs	857RA 0018570	10 ft 3m
			Nellcor	858RA 0018580	12ft 3.65m
			Pacetek	859RA 0018590	12ft 3.65m
			Simed Baxter	860RA 0018600	12ft 3.65m
			BCI	861RA 0018610	3 ft 0.9m
			Criticare	862RA 0018620	10 ft 3m
			Datascope	863RA 0018630	3 ft 0.9m
			Datascope	864RA 0018640	12ft 3.65m
			Sensormedics	865RA 0018650	3 ft 0.9m
			Simed Baxter	866RA 0018660	12ft 3.65m
			Ohmeda	867RA 0018670	12ft 3.65m
			Criticare	868RA 0018680	3 ft 0.9m
			Criticare	869RA 0018690	3 ft 0.9m





















Finger Probe Chart

			Criticare	870RA 0018700	10 ft 3m
			Nonin	871RA 0018710	3 ft 0.9m
			Datex	872RA 0018720	12ft 3.65m
			Datex	873RA 0018730	3 ft 0.9m
			Palco	874RA 0018740	3 ft 0.9m
			Novametrix	875RA 0018750	10 ft 3m
			Novametrix	876RA 0018760	10 ft 3m
			Critikon	877RA 0018770	3 ft 0.9m
			Invivo	878RA 0018780	10 ft 3m
			Artema	879RA 0018790	12ft 3.65m
			Novametrix	885RA 0018850	3ft 0.9m
			Kontron	886RA 0018860	12ft 3.65m
			Spacelabs	887RA 0018870	10 ft 3m
			S&W	888RA 0018880	12ft 3.65m
			Criticare	889RA 0018890	10 ft 3m
			Critikon & Sensormedics	890RA 0018900	10 ft 3m
			Critikon	891RA 0018910	10 ft 3m































Finger Probe Chart

			Marquette	892RA 0018920	10 ft 3m
			Hewlett Packard	895RA 0018950	10 ft 3m
			S&W	896RA 0018960	12ft 3.65m







Finger Probe Chart

			Nihon / Kohden	852RA 0018520	10 ft 3m
			Odam Bruker	853RA 0018530	10 ft 3m
			Medlab/ Comdek	854RA 0018540	10 ft 3m
			BCI	855RA 0018550	10 ft 3m
			Nellcor	856RA 0018560	3 ft 0.9m
			Spacelabs	857RA 0018570	10 ft 3m
			Nellcor	858RA 0018580	12ft 3.65m
			Pacetech	859RA 0018590	12ft 3.65m
			Simed Baxter	860RA 0018600	12ft 3.65m
			BCI	861RA 0018610	3 ft 0.9m
			Criticare	862RA 0018620	10 ft 3m
			Datascope	863RA 0018630	3 ft 0.9m
			Datascope	864RA 0018640	12ft 3.65m
			Sensormedics	865RA 0018650	3 ft 0.9m
			Simed Baxter	866RA 0018660	12ft 3.65m
			Ohmeda	867RA 0018670	12ft 3.65m
			Criticare	868RA 0018680	3 ft 0.9m
			Criticare	869RA 0018690	3 ft 0.9m













Finger Probe Chart

			Criticare	870RA 0018700	10 ft 3m
			Nonin	871RA 0018710	3 ft 0.9m
			Datex	872RA 0018720	12ft 3.65m
			Datex	873RA 0018730	3 ft 0.9m
			Palco	874RA 0018740	3 ft 0.9m
			Novametrix	875RA 0018750	10 ft 3m
			Novametrix	876RA 0018760	10 ft 3m
			Critikon	877RA 0018770	3 ft 0.9m
			Invivo	878RA 0018780	10 ft 3m
			Artema	879RA 0018790	12ft 3.65m
			Novametrix	885RA 0018850	3ft 0.9m
			Kontron	886RA 0018860	12ft 3.65m
			Spacelabs	887RA 0018870	10 ft 3m
			S&W	888RA 0018880	12ft 3.65m
			Criticare	889RA 0018890	10 ft 3m
			Critikon & Sensormedics	890RA 0018900	10 ft 3m
			Critikon	891RA 0018910	10 ft 3m






























Finger Probe Chart

			Marquette	892RA 0018920	10 ft 3m
			Hewlett Packard	895RA 0018950	10 ft 3m
			S&W	896RA 0018960	12ft 3.65m







Finger Probe Chart

			Nihon / Kohden	852RA 0018520	10 ft 3m
			Odam Bruker	853RA 0018530	10 ft 3m
			Medlab/ Comdek	854RA 0018540	10 ft 3m
			BCI	855RA 0018550	10 ft 3m
			Nellcor	856RA 0018560	3 ft 0.9m
			Spacelabs	857RA 0018570	10 ft 3m
			Nellcor	858RA 0018580	12ft 3.65m
			Pacetek	859RA 0018590	12ft 3.65m
			Simed Baxter	860RA 0018600	12ft 3.65m
			BCI	861RA 0018610	3 ft 0.9m
			Criticare	862RA 0018620	10 ft 3m
			Datascope	863RA 0018630	3 ft 0.9m
			Datascope	864RA 0018640	12ft 3.65m
			Sensormedics	865RA 0018650	3 ft 0.9m
			Simed Baxter	866RA 0018660	12ft 3.65m
			Ohmeda	867RA 0018670	12ft 3.65m
			Criticare	868RA 0018680	3 ft 0.9m
			Criticare	869RA 0018690	3 ft 0.9m

Finger Probe Chart

			Criticare	870RA 0018700	10 ft 3m
			Nonin	871RA 0018710	3 ft 0.9m
			Datex	872RA 0018720	12ft 3.65m
			Datex	873RA 0018730	3 ft 0.9m
			Palco	874RA 0018740	3 ft 0.9m
			Novametrix	875RA 0018750	10 ft 3m
			Novametrix	876RA 0018760	10 ft 3m
			Critikon	877RA 0018770	3 ft 0.9m
			Invivo	878RA 0018780	10 ft 3m
			Artema	879RA 0018790	12ft 3.65m
			Novametrix	885RA 0018850	3ft 0.9m
			Kontron	886RA 0018860	12ft 3.65m
			Spacelabs	887RA 0018870	10 ft 3m
			S&W	888RA 0018880	12ft 3.65m
			Criticare	889RA 0018890	10 ft 3m
			Critikon & Sensormedics	890RA 0018900	10 ft 3m
			Critikon	891RA 0018910	10 ft 3m

Finger Probe Chart

			Marquette	892RA 0018920	10 ft 3m
			Hewlett Packard	895RA 0018950	10 ft 3m
			S&W	896RA 0018960	12ft 3.65m

MCI. NO.	MANUFACTUR	Identification	LENGTH_FT
P856RA	Baxter	Explorer Monarch	3.00
P856RA	Corometrics	556.00	3.00
P856RA	Critikon	9700 series	3.00
P856RA	Draegar	Vitalert 1000	3.00
		Narkomed Cicero	
P856RA	3.00	BP806	3.00
P856RA	Entri Medical	Propaq 102 104 106	3.00
P856RA	Fukuda Denshi	HG302 Module	3.00
P856RA	Gambro	BP806	3.00
P856RA	Hewlett Packard	Miracle Adult system	3.00
		Merlin M1020A	
		Neonatal	
P856RA	Invivo	4500+2	3.00
P856RA	Ivy Medical	405A	3.00
P856RA	Marquette	Eagle 3100	3.00
P856RA	Medical Data	2800-50 OP11A	3.00
		OP30 Escort 100	
		Escort II	
P856RA	Mennen	800-060-730 Horizon	3.00
		XL	
P856RA	Nellcor	C-13-200M N-200	3.00
		DS100A	
P856RA	3.00	1281.00	3.00
P856RA	Propaq	Propaq	3.00
P856RA	S & W Vickers	4621643 4700708	3.00
		DSNT	
P856RA	Siemens	1281.00	3.00
P856RA	Spacelabs	015-0150-00 90489	
		90465 90466 90467	
		90651A-08	
	3.00		
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	3.00	N6000,N180E,N180	
P857RA	Spacelabs	90651A-08	12.00
P858RA	Critikon	9700 series	11.00
P858RA	Nellcor	N6000 N180E N180	0.00
P860RA	Baxter	600333-003,ASAT	
	11.00		
P859RA	Pacetek	520.00	0.00
P860RA	Baxter	600333-003 ASAT	
P861RA	0.00	3044,3300,3301,3302	
P860RA	Simed	S-100e	
	0.00		
P861RA	Biochem (BCI)	3044 3300 3301	3.00
		3302	
3044.00	Elmed	6006194,PSO1000	
	3.00		

P861RA	Bruker	Bruker	3.00
P861RA	Elmed	6006194 PSO1000	
	3.00		
P862RA	Criticare	511SD/934 1100	12.00
		1100-3 501 503 504	
		507 508 581 Poet	
		Plus	
P863RA	Datascope	0600-00-0026-02	
		Passport	
	12.00		
P864RA	Datascope	0600-00-0026-01	12.00
		Accustat	
P865RA	Sensor medics	766721 767367	
		Oxyshuttle Omnistat	
P867RA	Marquette	Tram system Modular	12.00
		Sys	
P867RA	Mennen	800-060-540,Horizon	
		XL	
	3.00		
P867RA	Marquette	Tram system Modular	12.00
		Sys	
P867RA	Mennen	800-060-540 Horizon	
		XL	
4700.00	12.00		
P867RA	Ohmeda	380-1000-042 6051-	
		0000-049 Oxystat	
		Modular CD Anaes	
		3700 4700	
,Oxystat Plus			
4700.00	12.00		
P870RA	Critikon	9084,8700 series	
	12.00		
P871RA	12.00	8604K,Onyx	0.00
P870RA	Critikon	9084 8700 series	
,CH-	12.00		
2S,SATLITE,ULTS0			
0,OSCAR II,M-			
ESTPR,CH-			
RS,ULTSVO,CARDI			
O CAP II			
P871RA	Nonin	8604K Onyx	0.00
P872RA	Datex	878579-1 SAS-F4	
		CG-2CS CH-S CH-	
		RS ULTSIO	
		ULTSVI AS/3 OSP-	
		200,	

CH-2S SATLITE ULTS00 OSCAR II M-ESTPR CH-RS ULTSVO CARDIO CAP II	12.00	369-083- 001,90419,90501,906 51A-03
P875RA	Novamatrix 0.00	8660 500 series
P875RA	Spacelabs 12.00	369-083-001 90419 90501 90651A-03
P876RA	0.00	015-0130-00,90651A- 08,90419- 01,90489,90467,9046 6,90465,90351-D/G
P876RA	Novamatrix 12.00	8774 8776 505 511 515 520 7000 Oxypleth
P876RA	Spacelabs	015-0130-00 90651A- 08 90419-01 90489 90467 90466 90465 90351-D/G
P877RA	12.00 Critikon	8997,Dinamap Oxytrak
P878RA	0.00 12.00	9383,Scout ES,Scout +4,4500Scout,Scout +,Scout+1,Scout+2,Sc out +3
P877RA	Critikon 0.00	8997 Dinamap Oxytrak
P878RA	Invivo	9383 Scout ES Scout +4 4500Scout Scout + Scout+1 Scout+2 Scout +3
P886RA	12.00 12.00 Kontron 12.00	78400 7840

Distributed by

Replacement SpO₂ Probe

This compatible replacement adult oxygen saturation and pulse rate monitoring SpO₂ probe, is for use when continuous non-invasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing greater than 40kg.

For use with product compatible oximeter (See table on reverse)

Instructions for Use

An index finger is the preferred location for the sensor.

Orientate the sensor so that the fingernail is directly under the sensor. Plug the sensor into the appropriate oximeter and verify proper operation in accordance with the original manufacturers oximeter manual.

NB: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff or intravascular infusion line.

Reusable sensors may be

used on the same site for a maximum of 4 hours, provided the site is inspected routinely to ensure skin integrity and correct positioning. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

NB: If the patient wears fingernail polish or artificial fingernails choose a different probe site, or remove the fingernail polish or artificial fingernails.

Troubleshooting

The probe may not reliably track pulse and SpO₂ if the finger is incorrectly positioned within the sensor, or the sensor site is too thick, thin, or deeply pigmented to permit appropriate light transmission. If any of these situations occur, reposition the sensor.

Avoid any site that may have excessive ambient light, excessive motion, or electrical interference

Not for use within MRI

scanning suites.

Cleaning

Sensor may be surface cleaned with 70% isopropyl alcohol.

Do not immerse the sensor or connector in any liquid solution.

Do not autoclave this probe.

The SpO₂ Finger probe will not be harmed by exposure to an ethylene oxide mixture at 49° to 54°C (120° to 130°F)

Note

Refer to the instrument operator's manual for complete instructions for use of the probe with the monitor.

Accuracy

Compatible to Original manufacturers specification

WARRANTY

This product is warranted for a period of 12 months from the date of purchase against defects in workmanship and materials by Viamed whose sole obligation under this

warranty is to repair or replace the product.

To return a product under warranty the local distributor should be contacted for a return authorisation.

Internal components are manufactured to our proprietary design specifications. Viamed does not have a licensing agreement with any Original Oximeter manufacturer

Manufactured By MCI USA.

Viamed Finger Probe Re-Cycle programme

This year will see the introduction of EEC legislation regarding re-cycling of waste and as an medical engineering based company we understand the frustration of broken pulse oximeter probes collecting dust whilst you earnestly wish you could repair them.

This is now possible as Viamed have brought onto line a finger probe re-cycling service and we are able to repair most probes. Moulded connectors obviously cause problems and may prevent us effecting a clean profession repair. We also have a learning curve and therefore with less popular versions of finger probes we may exceed our target of a maximum 7 day turn around.

The aim is to return to the user a finger probe that looks new and has not just been glued together or spliced and covered with heatshrink tubing. Our object is to re-cycle as much as possible of the existing finger probe instead of consigning it to the waste bin. Where required new parts will be fitted. These probes are used in critical care areas and it is our intention not only to repair them but to return them clean and professionally finished so they look the part.

Our warranty covers the parts and labour used to repair your probe. i.e. if we replace the clip and the clip breaks in 6 months we repair it free. However if we repair the clip then later the instrument plug develops a fault this is treated as a separate repair.

Viamed is ISO9002/EN46002 certified and all repairs will be carried out in accordance with these procedures. Before despatch Every re-cycled finger probe is tested electronically and in many cases on the original manufacturers equipment. As the number of different types of probe we repair increases so does the number of O.E.M' SpO2 instruments we acquire. Full reports and tests carried out are documented and supplied with each probe.

The proof is always easy to establish by trying us out at no risk. If we cannot repair to your satisfaction there is no charge. If we are not satisfied, under our BS EN ISO 9002/EN46002 registration, with the finished product or its accuracy we will not claim to have repaired it.

Our repair Policy is :-

- We will endeavour to provide quality repairs
- All repairs are made using, when required, new replacement parts
- Prompt service and fast turn around
- Six months warranty on all repairs
- There is a fixed cost regardless of type of cable or extent of repair
- We can usually only repair Finger Probes at present
- We cannot always repair OEM cables with special moulded plugs although we try.
- If we cannot repair a probe a service replacement will be offered at 20% off our list price.
- Most Finger Probes can be re-cycled
- The probes should be decontaminated before they are returned to us in accordance with the official recommendations.

The fixed cost of the repair/ re-cycle service is £80.00. + VAT & Carriage

John S Lamb.

MANUFACTURER	COMMENTS	MODEL1	EPIC_NO	EPIC_NEW_	Tele No
Baxter	Finger Sensor	ASAT	E333	E112-12	P860
Biochem (BCF)	Finger Sensor	3301	E344	E103-06	P861
Criticare	Finger Sensor	1100/1100-3	E511*	E112-05	P862
Criticare	Adapter Cable	Poet Plus	E5118L	E790-05	P869
Criticare	Extension cable	503 Spot	E5118D	E710-05	P868
Critikon	Finger Sensor		E7667	E112-10	P870
Critikon	Finger Sensor	8997	E8997	E103-10	P877
Datascope	Ext.Cable 3ft	Passport	E012	E710-07	P855
Datex	Finger Sensor	CG-2CS	E8785	E112-09	P872
Draeger	Pre Amp Cable	Narkomed	EN300	E610-01	P879
Invivo	Finger Sensor	4500	E9000	E112-14	P878
J & J Medical	Finger Sensor		E180	E111-01	P858
Kontron	Finger sensor	7840	E8010	E112-13	P886
Nellcor	Finger Sensor	All Oximeters	E100A	E103-01	P856
Nellcor	Extension cable	N180E	E004	E704-01	P853
Nellcor	Extension cable	N6000	E008	E708-01	P854
Nonin	Finger Sensor	8500/8600/Oximeter	E8000*	E103-08	P871
Nonin	Extension Cable 3ft		E8500	E703-08	P873
Nonin	Extension Cable 10ft		E8604	E710-08	P874
Novamatrix	Finger sensor	500	E8700	E192-03	P875
Novamatrix	Finger Sensor		E8700/5	E112-03	P876
Ohmeda	Finger Sensor	Oxystat	E3000	E112-02	P867
Pacetek	Finger Sensor	520	E201	E112-11	P859
Spacelabs	Extension cable	90651A	E646	E702-04	P865
Spacelabs	Extension cable	90465	E706	E702-04	P866
Spacelabs	No adapter direct to oximeter	90467	E100SL	E112-04	P857
Nellcor	Flexi-site			E203-01	P915
Spacelabs	Flexi-site			E212-04	P916
Nellcor	Flexi-site			E211-01	P917
Pacetek	Flexi-site			E212-11	P918
Simed	Flexi-site			E212-12	P919
Biochem (BCF)	Flexi-site			E203-06	P920
Criticare	Flexi-site			E212-05	P921
Ohmeda	Flexi-site			E212-02	P922
Critikon	Flexi-site			E212-10	P923
Nonin	Flexi-site Adult			E203-08	P924
Kontron	Flexi-site			E212-13	P925
Novamatrix	Flexi-site			E212-03	P926
Datex	Flexi-site 1Mtr			E212-09	P927
Critikon	Flexi-site			E203-10	P928
Invivo	Flexi-site			E212-14	P929
S & W Vickers	SSI one piece	Poet	EF511	E212-05	P862
Datascope	Finger Sensor 12Ft	Passport	E600		P863
Datascope	Finger Probe 3Ft	Accustat	E601		P864
					P911