Service Repair Sheet SRS68380

Contact Name

Company/ Hospital Name

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

Mark Evans

Royal United Hospitals Bath

Medical Equipment Management

01225824140

01225 428331

00000350

mark.evans7@nhs.net

21/Feb/2023

Robert Connor

Viamed

Warranty

VIAMEDcle	ar
Goods In Only Decontamination certificate provided by customer	Ε
Cleaned by Viamed, if no declaration certificate from customer Signed:	
Goods Out Only Cleaned by Viamed beforeturning to customer	ore

Notes 02/Feb/2023 Steve Hardaker

02/Feb/2023 Steve Hardaker

Returning a Max-550E s/n FM82299148 sold 7/4/21 with a MaxO2 ME. Still under warranty.

Intermittent output, when moved or tilted it suddenly drops very low. He has tried it on anther MaxO2 ME and the fault m with the sensor.

21/Feb/2023 Robert Connor

Received 1 x MAX-550E s/n FM82299148, with flow diverter

Ready For quote

Repair Complete Signed

CGNOON 21.2.23

SRN	Equipment	Stock Ref	Serial Number	Warranty
SRN35678	Oxygen Sensor	0110452	FM82299148	Yes

0140013 × 1 FOC 01d S/N, SRS, SRN

New S/N, SRS, SRN

UPS × 1 FOC.

Royal United Hospitals Bath NHS

NHS Foundation Trust

Medical Equipment Management Service

Dispatch note (external service)

Our reference

266639

Supplier

VIAMED LTD

Order date

09/02/2023

Address

15 STATION ROAD, CROSS HILLS

Our Reference

WO266639

KEIGHLEY, W YORKSHIRE BD20 7DT

Contact person

MARK EVANS

Your reference Phone

SRS68380

Phone

Email

01225 824154

01535 634542

Fax

Fax

01535 635582

Mobile phone

mark.evans7@nhs.net

Package no Type of service

WARRANTY WORK

MEMS no.

49104

Serial no

FM72099004

Device type **Brand**

OXYGEN MONITOR

MAXTEC

MAX O2 ME

Model Owner

0707

NICU; RUH

Location

0707

NICU; RUH

Message

Hi

This is just the sensor MAX550E SN: FM82299148 which gives varying readings of air. It can be stable at 21% then fall as low as 6% or even 0% if moved. Have tested with two oxygen monitors with the same results.

Please replace under warranty.

Many thanks

Mark Evans

Return date

Signature

The am

Before commencing an in-depth repair, could you please advise us by telephone of the likely cost if over £ by asking for the contact person above. This letter serves as a certificate that the device has been decontminated to standard local protocols. PLEASE RETURN ITEM TO DELIVERY ADDRESS BELOW

The goods are to be returned to the address below

Invoice address

Finance Dept, Malvern House -

Building 6

Royal United Hospital

Combe Park **BATH** BA1 3NG

Delivery address

MEMS - Department A8 ROYAL UNITED HOSPITAL

COMBE PARK **BATH** BA1 3NG

> This item has been decontaminated

Si	an	at	ur	e
-	э.,	~-	ч.	~

Signed Date





Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdor Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

Service Repair Sheet 68380

Contact Name

Company/ Hospital Name

Department **Direct Phone** General Phone

Opera Account

Email

Date Received Booked in By

Mark Evans

Royal United Hospitals Bath

Medical Equipment Management Service

01225824140 01225 428331

00000350

mark.evans7@nhs.net

21/Feb/2023 Robert Connor

Repair SRN35678 Ref

S/N 0110452 FM82299148 **Equipment Type**

Oxygen Sensor

Under Warranty

Time: 0 Hour(s)

Parts Replaced

0110452 Qty: 1

Replaced under warranty with S/N HH82199042