


versastream luer lock male selection guide



Is the patient intubated?

Yes

Long-Term
≤ 72 hours

Airway Adapter



Adult/Paediatric
4420873

Infant
4420874



Airway adapter
connector sizes:
Adult/Paediatric 15 mm O.D.
to 15 mm I.D. / 22 mm O.D.
Infant 15 mm O.D.
to 4 mm I.D. / 15 mm I.D.



Luer Lock to Luer Lock
Sampling Line
for use with 4420531
4420876



Airway Adapter with Luer
Lock Female Connection
e.g. Adult/Paediatric
4420531

No

Long-Term
≤ 72 hours

Does the patient require oxygen?

Yes

Nasal



Adult
4420860

Paediatric
4420861

Nasal/Oral



Adult
4420869

Paediatric
4420870



For use with patients
that exceed 10 kg



For use with patients
that are less than 10 kg

No

Nasal



Adult
4420857

Paediatric
4420858

Nasal/Oral



Adult
4420867

Paediatric
4420868



Part Number: 4490117
Date: 03/23