

A6.4 Test Results P867RA

Oximeter Monitor: Ohmeda 3740Monitor S/N: FMQU02027SpO₂ TestingFinger Probe Model: P867RABioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		98	94	80	70	60	98	94	80	70	60
9731-4	7H00978	98	94	80	70	60	0	0	0	0	0
9731-4	7H00979	98	94	80	70	60	0	0	0	0	0
9731-4	7H00980	98	94	80	70	60	0	0	0	0	0
9731-4	7H00981	98	94	80	70	60	0	0	0	0	0
9731-4	7H00982	98	94	80	70	60	0	0	0	0	0
9731-4	7H00983	98	94	80	70	60	0	0	0	0	0
9731-4	7H00984	98	94	80	70	60	0	0	0	0	0
9731-4	7H00985	98	94	80	70	60	0	0	0	0	0
9731-4	7H00986	98	94	80	70	60	0	0	0	0	0
9731-4	7H00987	98	94	80	70	60	0	0	0	0	0
9731-4	7H00988	98	94	80	70	60	0	0	0	0	0
9731-4	7H00989	98	94	80	70	60	0	0	0	0	0
9731-4	7H00990	98	94	80	70	60	0	0	0	0	0
9731-4	7H00991	98	94	80	70	60	0	0	0	0	0
9731-4	7H00992	98	94	80	70	60	0	0	0	0	0
9731-4	7H00993	98	94	80	70	60	0	0	0	0	0
9731-4	7H00994	98	94	80	70	60	0	0	0	0	0
9731-4	7H00995	98	94	80	70	60	0	0	0	0	0
9731-4	7H00996	98	94	80	70	60	0	0	0	0	0
9731-4	7H00997	98	94	80	70	60	0	0	0	0	0
9731-4	7H00998	98	94	80	70	60	0	0	0	0	0
9731-4	7H00999	98	94	80	70	60	0	0	0	0	0
9731-4	7H01000	98	94	80	70	60	0	0	0	0	0
9731-4	7H01001	98	94	80	70	60	0	0	0	0	0
9731-4	7H01002	99	94	80	70	60	1	0	0	0	0
9731-4	7H01003	97	94	80	70	60	-1	0	0	0	0
9731-4	7H01004	97	94	80	70	60	-1	0	0	0	0
9731-4	7H01005	98	94	80	70	60	0	0	0	0	0
9731-4	7H01006	98	94	80	70	60	0	0	0	0	0
9731-4	7H01007	98	94	80	70	60	0	0	0	0	0
9731-4	7H01008	98	94	80	70	60	0	0	0	0	0
9731-4	7H01009	98	94	80	70	60	0	0	0	0	0
9731-4	7H01010	98	94	80	70	60	0	0	0	0	0
9731-4	7H01011	99	94	80	70	60	1	0	0	0	0
9731-4	7H01012	98	94	80	70	60	0	0	0	0	0
9731-4	7H01013	98	94	80	70	60	0	0	0	0	0
9731-4	7H01014	98	94	80	70	60	0	0	0	0	0
9731-4	7H01015	98	94	80	70	60	0	0	0	0	0
9731-4	7H01016	99	94	80	70	60	1	0	0	0	0
9731-4	7H01017	98	94	80	70	60	0	0	0	0	0
9731-4	7H01018	98	94	80	70	60	0	0	0	0	0
9731-4	7H01019	99	94	80	70	60	1	0	0	0	0
9731-4	7H01020	98	94	80	70	60	0	0	0	0	0

Oximeter Monitor: Ohmeda 3740
Monitor S/N: FMQU02027

Finger Probe Model: P867RA
BioTek Simulator S/N: 125672

9731-4	7H01021	98	94	80	70	60	0	0	0	0	0
9731-4	7H01022	99	94	80	70	60	1	0	0	0	0
9731-4	7H01023	99	94	80	70	60	1	0	0	0	0
9731-4	7H01024	98	94	80	70	60	0	0	0	0	0
9731-4	7H01025	98	94	80	70	60	0	0	0	0	0
9731-4	7H01026	98	94	80	70	60	0	0	0	0	0
9731-4	7H01027	99	94	80	70	60	1	0	0	0	0
9731-2	7H01028	98	94	80	70	60	0	0	0	0	0
9731-2	7H01029	99	94	80	70	60	1	0	0	0	0
9731-2	7H01030	99	94	80	70	60	1	0	0	0	0
9731-2	7H01031	99	94	80	70	60	1	0	0	0	0
9731-2	7H01032	98	94	80	70	60	0	0	0	0	0
9731-2	7H01033	98	94	80	70	60	0	0	0	0	0
9731-2	7H01034	99	94	80	70	60	1	0	0	0	0
9731-2	7H01035	99	94	80	70	60	1	0	0	0	0
9731-2	7H01036	98	94	80	70	60	0	0	0	0	0
9731-2	7H01037	97	94	80	70	60	-1	0	0	0	0
9731-2	7H01038	99	94	80	70	60	1	0	0	0	0
9731-2	7H01039	98	94	80	70	60	0	0	0	0	0
9731-2	7H01040	98	94	80	70	60	0	0	0	0	0
9731-2	7H01041	98	94	80	70	60	0	0	0	0	0
9731-2	7H01042	98	94	80	70	60	0	0	0	0	0
9731-2	7H01043	99	94	80	70	60	1	0	0	0	0
9731-2	7H01044	99	94	80	70	60	1	0	0	0	0
9731-2	7H01045	99	94	80	70	60	1	0	0	0	0
9731-2	7H01046	98	94	80	70	60	0	0	0	0	0
9731-2	7H01047	99	94	80	70	60	1	0	0	0	0
9731-2	7H01048	98	94	80	70	60	0	0	0	0	0
9731-2	7H01049	99	94	80	70	60	1	0	0	0	0
9731-2	7H01050	98	94	80	70	60	0	0	0	0	0
9731-2	7H01051	98	94	80	70	60	0	0	0	0	0
9731-2	7H01052	98	94	80	70	60	0	0	0	0	0
9731-5	7H01053	98	94	80	70	60	0	0	0	0	0
9731-5	7H01054	98	94	80	70	60	0	0	0	0	0
9731-5	7H01055	98	94	80	70	60	0	0	0	0	0
9731-5	7H01056	98	94	80	70	60	0	0	0	0	0
9731-5	7H01057	98	94	80	70	60	0	0	0	0	0
9731-5	7H01058	98	94	80	70	60	0	0	0	0	0
9731-5	7H01059	98	94	80	70	60	0	0	0	0	0
9731-5	7H01060	98	94	80	70	60	0	0	0	0	0

A6.4 P867RA Test Results
Oximeter Monitor: Ohmeda 3740
Monitor S/N: FMQU02027

Pulse Rate Testing
Finger Probe Model: P867RA
BioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		250	150	95	80	50	250	150	95	80	50
9731-4	7H00978	250	150	95	80	50	0	0	0	0	0
9731-4	7H00979	250	150	95	80	50	0	0	0	0	0
9731-4	7H00980	250	150	95	80	50	0	0	0	0	0
9731-4	7H00981	250	150	95	80	50	0	0	0	0	0
9731-4	7H00982	250	150	95	80	50	0	0	0	0	0
9731-4	7H00983	250	150	95	80	50	0	0	0	0	0
9731-4	7H00984	250	150	95	80	50	0	0	0	0	0
9731-4	7H00985	250	150	95	80	50	0	0	0	0	0
9731-4	7H00986	250	150	95	80	50	0	0	0	0	0
9731-4	7H00987	250	150	95	80	50	0	0	0	0	0
9731-4	7H00988	250	150	95	80	50	0	0	0	0	0
9731-4	7H00989	250	150	95	80	50	0	0	0	0	0
9731-4	7H00990	252	150	95	80	50	2	0	0	0	0
9731-4	7H00991	250	150	95	80	50	0	0	0	0	0
9731-4	7H00992	250	150	95	80	50	0	0	0	0	0
9731-4	7H00993	250	150	95	80	50	0	0	0	0	0
9731-4	7H00994	250	150	95	80	50	0	0	0	0	0
9731-4	7H00995	250	150	95	80	50	0	0	0	0	0
9731-4	7H00996	250	150	95	80	50	0	0	0	0	0
9731-4	7H00997	250	150	95	80	50	0	0	0	0	0
9731-4	7H00998	250	150	95	80	50	0	0	0	0	0
9731-4	7H00999	250	150	95	80	50	0	0	0	0	0
9731-4	7H01000	250	150	95	80	50	0	0	0	0	0
9731-4	7H01001	250	150	95	80	50	0	0	0	0	0
9731-4	7H01002	250	150	95	80	50	0	0	0	0	0
9731-4	7H01003	250	150	95	80	50	0	0	0	0	0
9731-4	7H01004	250	150	95	80	50	0	0	0	0	0
9731-4	7H01005	251	150	95	80	50	1	0	0	0	0
9731-4	7H01006	250	150	95	80	50	0	0	0	0	0
9731-4	7H01007	250	150	95	80	50	0	0	0	0	0
9731-4	7H01008	250	150	95	80	50	0	0	0	0	0
9731-4	7H01009	250	150	95	80	50	0	0	0	0	0
9731-4	7H01010	250	150	95	80	50	0	0	0	0	0
9731-4	7H01011	250	150	95	80	50	0	0	0	0	0
9731-4	7H01012	250	150	95	80	50	0	0	0	0	0
9731-4	7H01013	250	150	95	80	50	0	0	0	0	0
9731-4	7H01014	250	150	95	80	50	0	0	0	0	0
9731-4	7H01015	250	150	95	80	50	0	0	0	0	0
9731-4	7H01016	252	150	95	80	50	2	0	0	0	0
9731-4	7H01017	252	150	95	80	50	2	0	0	0	0
9731-4	7H01018	250	150	95	80	50	0	0	0	0	0
9731-4	7H01019	250	150	95	80	50	0	0	0	0	0
9731-4	7H01020	250	150	95	80	50	0	0	0	0	0

Oximeter Monitor: Ohmeda 3740
Monitor S/N: FMQU02027

Finger Probe Model: P867RA
BioTek Simulator S/N: 125672

9731-4	7H01021	250	150	95	80	50	0	0	0	0	0
9731-4	7H01022	250	150	95	80	50	0	0	0	0	0
9731-4	7H01023	250	150	95	80	50	0	0	0	0	0
9731-4	7H01024	250	150	95	80	50	0	0	0	0	0
9731-4	7H01025	250	150	95	80	50	0	0	0	0	0
9731-4	7H01026	250	150	95	80	50	0	0	0	0	0
9731-4	7H01027	250	150	95	80	50	0	0	0	0	0
9731-2	7H01028	250	150	95	80	50	0	0	0	0	0
9731-2	7H01029	250	150	95	80	50	0	0	0	0	0
9731-2	7H01030	250	150	95	80	50	0	0	0	0	0
9731-2	7H01031	250	150	95	80	50	0	0	0	0	0
9731-2	7H01032	250	150	95	80	50	0	0	0	0	0
9731-2	7H01033	250	150	95	80	50	0	0	0	0	0
9731-2	7H01034	250	150	95	80	50	0	0	0	0	0
9731-2	7H01035	250	150	95	80	50	0	0	0	0	0
9731-2	7H01036	250	150	95	80	50	0	0	0	0	0
9731-2	7H01037	250	150	95	80	50	0	0	0	0	0
9731-2	7H01038	250	150	95	80	50	0	0	0	0	0
9731-2	7H01039	250	150	95	80	50	0	0	0	0	0
9731-2	7H01040	250	150	95	80	50	0	0	0	0	0
9731-2	7H01041	250	150	95	80	50	0	0	0	0	0
9731-2	7H01042	250	150	95	80	50	0	0	0	0	0
9731-2	7H01043	250	150	95	80	50	0	0	0	0	0
9731-2	7H01044	250	150	95	80	50	0	0	0	0	0
9731-2	7H01045	250	150	95	80	50	0	0	0	0	0
9731-2	7H01046	250	150	95	80	50	0	0	0	0	0
9731-2	7H01047	250	150	95	80	50	0	0	0	0	0
9731-2	7H01048	250	150	95	80	50	0	0	0	0	0
9731-2	7H01049	250	150	95	80	50	0	0	0	0	0
9731-2	7H01050	250	150	95	80	50	0	0	0	0	0
9731-2	7H01051	250	150	95	80	50	0	0	0	0	0
9731-2	7H01052	250	150	95	80	50	0	0	0	0	0
9731-5	7H01053	252	150	95	80	50	2	0	0	0	0
9731-5	7H01054	251	150	95	80	50	1	0	0	0	0
9731-5	7H01055	250	150	95	80	50	0	0	0	0	0
9731-5	7H01056	250	150	95	80	50	0	0	0	0	0
9731-5	7H01057	250	150	95	80	50	0	0	0	0	0
9731-5	7H01058	251	150	95	80	50	1	0	0	0	0
9731-5	7H01059	252	150	95	80	50	2	0	0	0	0
9731-5	7H01060	250	150	95	80	50	0	0	0	0	0

Oximeter Monitor: Ohmeda 3740
Monitor S/N: FMOU02027

Finger Probe Model: P867RA
BioTek Simulator S/N: 125672

9731-5	7H01061	98	94	80	70	60	0	0	0	0	0
9731-5	7H01062	98	94	80	70	60	0	0	0	0	0
9731-5	7H01063	98	94	80	70	60	0	0	0	0	0
9731-5	7H01064	98	94	80	70	60	0	0	0	0	0
9731-5	7H01065	99	94	80	70	60	1	0	0	0	0
9731-5	7H01066	98	94	80	70	60	0	0	0	0	0
9731-5	7H01067	98	94	80	70	60	0	0	0	0	0
9731-5	7H01068	98	94	80	70	60	0	0	0	0	0
9731-5	7H01069	98	94	80	70	60	0	0	0	0	0
9731-5	7H01070	98	94	80	70	60	0	0	0	0	0
9731-5	7H01071	98	94	80	70	60	0	0	0	0	0
9731-5	7H01072	98	94	80	70	60	0	0	0	0	0
9731-5	7H01073	98	94	80	70	60	0	0	0	0	0
9731-5	7H01074	98	94	80	70	60	0	0	0	0	0
9731-5	7H01075	98	94	80	70	60	0	0	0	0	0
9731-5	7H01076	98	94	80	70	60	0	0	0	0	0

Tester: [Signature]

Date: 9/26/97

Q.A.: [Signature]

Date: 9/26/97

Oximeter Monitor: Ohmeda 3740
Monitor S/N: FMQU02027

Finger Probe Model: P867RA
BioTek Simulator S/N: 125672

9731-5	7H01061	252	150	95	80	50	2	0	0	0	0
9731-5	7H01062	251	150	95	80	50	1	0	0	0	0
9731-5	7H01063	250	150	95	80	50	0	0	0	0	0
9731-5	7H01064	251	150	95	80	50	1	0	0	0	0
9731-5	7H01065	250	150	95	80	50	0	0	0	0	0
9731-5	7H01066	252	150	95	80	50	2	0	0	0	0
9731-5	7H01067	250	150	95	80	50	0	0	0	0	0
9731-5	7H01068	250	150	95	80	50	0	0	0	0	0
9731-5	7H01069	250	150	95	80	50	0	0	0	0	0
9731-5	7H01070	252	150	95	80	50	2	0	0	0	0
9731-5	7H01071	250	150	95	80	50	0	0	0	0	0
9731-5	7H01072	250	150	95	80	50	0	0	0	0	0
9731-5	7H01073	250	150	95	80	50	0	0	0	0	0
9731-5	7H01074	251	150	95	80	50	1	0	0	0	0
9731-5	7H01075	250	150	95	80	50	0	0	0	0	0
9731-5	7H01076	251	150	95	80	50	1	0	0	0	0

Tester: [Signature]
Q.A.: [Signature]

Date: 9/24/97
Date: 9/26/97

A6.5 P876RA Test Results
Oximeter Monitor: Novamatrix 515A
Monitor S/N: 69-1010A

SpO₂ Testing
Finger Probe Model: P876RA
BioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		98	94	80	70	60	98	94	80	70	60
9731-1	7I01380	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01382	98	94	80	70	60	0	0	0	0	0
9731-1	7I01385	98	94	80	70	60	0	0	0	0	0
9731-1	7I01386	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01387	98	94	80	70	60	0	0	0	0	0
9731-1	7I01388	98	94	80	70	60	0	0	0	0	0
9731-1	7I01389	98	94	80	70	60	0	0	0	0	0
9731-1	7I01390	98	94	80	70	60	0	0	0	0	0
9731-1	7I01391	98	94	80	70	60	0	0	0	0	0
9731-1	7I01392	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01393	98	94	80	70	60	0	0	0	0	0
9731-1	7I01394	98	94	80	70	60	0	0	0	0	0
9731-1	7I01395	98	94	80	70	60	0	0	0	0	0
9731-1	7I01396	98	94	80	70	60	0	0	0	0	0
9731-1	7I01397	98	94	80	70	60	0	0	0	0	0
9731-1	7I01398	98	94	80	70	60	0	0	0	0	0
9731-1	7I01399	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01400	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01401	98	94	80	70	60	0	0	0	0	0
9731-1	7I01402	98	94	80	70	60	0	0	0	0	0
9731-1	7I01403	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01404	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01405	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01406	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01407	98	94	80	70	60	0	0	0	0	0
9731-1	7I01408	98	94	80	70	60	0	0	0	0	0
9731-1	7I01409	98	94	80	70	60	0	0	0	0	0
9731-1	7I01410	98	94	80	70	60	0	0	0	0	0
9731-1	7I01411	98	94	80	70	60	0	0	0	0	0
9731-1	7I01412	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01413	98	94	80	70	60	0	0	0	0	0
9731-1	7I01414	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01415	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01416	98	94	80	70	60	0	0	0	0	0
9731-1	7I01417	98	94	80	70	60	0	0	0	0	0
9731-1	7I01418	98	94	80	70	60	0	0	0	0	0
9731-1	7I01419	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01420	98	94	80	70	60	0	0	0	0	0
9731-1	7I01421	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01422	98	94	80	70	60	0	0	0	0	0
9731-1	7I01423	98	94	80	70	60	0	0	0	0	0
9731-1	7I01424	98	94	80	70	60	0	0	0	0	0
9731-1	7I01425	98	94	80	70	60	0	0	0	0	0
9731-1	7I01426	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01427	98	94	80	70	59	0	0	0	0	-1
9731-1	7I01428	98	94	80	70	60	0	0	0	0	0
9731-1	7I01429	98	94	80	70	60	0	0	0	0	0
9731-1	7I01430	98	94	80	70	60	0	0	0	0	0
9731-1	7I01431	98	94	80	70	60	0	0	0	0	0

Tester: *[Signature]*

Date: 9/26/97

Q.A.: *[Signature]*

Date: 9/26/97

A6.5 P876RA Test Results

Oximeter Monitor: Novamatrix 515A
Monitor S/N: 69-1010A

Pulse Rate Testing

Finger Probe Model: P876RA
BioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		250	150	95	80	50	250	150	95	80	50
9731-1	7I01380	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01382	200	160	100	80	60	0	0	0	0	0
9731-1	7I01385	200	160	100	80	60	0	0	0	0	0
9731-1	7I01386	200	160	100	80	60	0	0	0	0	0
9731-1	7I01387	200	160	100	80	60	0	0	0	0	0
9731-1	7I01388	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01389	200	160	100	80	60	0	0	0	0	0
9731-1	7I01390	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01391	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01392	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01393	200	160	100	80	60	0	0	0	0	0
9731-1	7I01394	200	160	100	80	60	0	0	0	0	0
9731-1	7I01395	200	160	100	80	60	0	0	0	0	0
9731-1	7I01396	200	160	100	80	60	0	0	0	0	0
9731-1	7I01397	200	160	100	80	60	0	0	0	0	0
9731-1	7I01398	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01399	200	160	100	80	60	0	0	0	0	0
9731-1	7I01400	200	160	100	80	60	0	0	0	0	0
9731-1	7I01401	200	160	100	80	60	0	0	0	0	0
9731-1	7I01402	200	160	100	80	60	0	0	0	0	0
9731-1	7I01403	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01404	200	160	100	80	60	0	0	0	0	0
9731-1	7I01405	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01406	200	160	100	80	60	0	0	0	0	0
9731-1	7I01407	200	160	100	80	60	0	0	0	0	0
9731-1	7I01408	200	160	100	80	60	0	0	0	0	0
9731-1	7I01409	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01410	200	160	100	80	60	0	0	0	0	0
9731-1	7I01411	200	160	100	80	60	0	0	0	0	0
9731-1	7I01412	200	160	100	80	60	0	0	0	0	0
9731-1	7I01413	200	160	100	80	60	0	0	0	0	0
9731-1	7I01414	200	160	100	80	60	0	0	0	0	0
9731-1	7I01415	200	160	100	80	60	0	0	0	0	0
9731-1	7I01416	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01417	200	160	100	80	60	0	0	0	0	0
9731-1	7I01418	200	160	100	80	60	0	0	0	0	0
9731-1	7I01419	200	160	100	80	60	0	0	0	0	0
9731-1	7I01420	200	160	100	80	60	0	0	0	0	0
9731-1	7I01421	200	160	100	80	60	0	0	0	0	0
9731-1	7I01422	200	160	100	80	60	0	0	0	0	0
9731-1	7I01423	200	160	100	80	60	0	0	0	0	0
9731-1	7I01424	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01425	200	160	100	80	60	0	0	0	0	0
9731-1	7I01426	200	160	100	80	60	0	0	0	0	0
9731-1	7I01427	200	160	100	80	60	0	0	0	0	0
9731-1	7I01428	200	160	100	80	60	0	0	0	0	0
9731-1	7I01429	200	160	100	80	60	0	0	0	0	0
9731-1	7I01430	199	160	100	80	60	-1	0	0	0	0
9731-1	7I01431	200	160	100	80	60	0	0	0	0	0

Tester: [Signature]
Q.A.: [Signature]

Date: 9/26/97
Date: 9/26/97