






Spectrometer	Prema SPM 9001
Manufacturer of probe	Nelcor
Type of probe	Finger
Part number	DS-100A
Serial number	C07895
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nelcor and compatible monitors
Tested on	Nelcor N-200
Red peak wavelength (nm)	658
Infra red peak wavelength (nm)	879.8
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998



Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	C08921
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength (nm)	661.3
Infra red peak wavelength (nm)	876.5
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998




Spectrometer: Prema SPM 9001

Manufacturer of probe: Nellcor
Part number: DS-100A
Serial number:
Cable length: 3 ft
Connector: 9 pin D male - pins 4 & 8 missing
Connects to: Extension cable


For use on: Nellcor
Tested on: Nellcor N-200

Red peak wavelength: 661
Infra red peak wavelength: 883


LED configuration: 2 lead




Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00176
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength	661.3
Infra red peak wavelength	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998




Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00236
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength (nm)	658
Infra red peak wavelength (nm)	870.1
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998




Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00326
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength (nm)	658
Infra red peak wavelength (nm)	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998




Spectrometer	Prema SPM 9001
Manufacturer of probe	Nelcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00376
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nelcor and compatible monitors
Tested on	Nelcor N-200
Red peak wavelength (nm)	664.6
Infra red peak wavelength (nm)	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998




Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00387
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength	661.3
Infra red peak wavelength	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998



Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	W00435
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength	661.3
Infra red peak wavelength	886.2
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998



Spectrometer	Prema SPM 9001
Manufacturer of probe	Neillcor
Type of probe	Finger
Part number	DS-100A
Serial number	C08091
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Neillcor and compatible monitors
Tested on	Neillcor N-200
Red peak wavelength	661.3
Infra red peak wavelength	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998



Spectrometer	Prema SPM 9001
Manufacturer of probe	Nellcor
Type of probe	Finger
Part number	DS-100A
Serial number	C08091
Cable length	3 ft
Connector	D type 9 pin male
Connects to	Monitor or extension cable
For use on	Nellcor and compatible monitors
Tested on	Nellcor N-200
Red peak wavelength	661.3
Infra red peak wavelength	883
Emitter configuration	2 lead (parallel, reverse)
Carried out by (name printed)	S Nixon
Signed	
Date	31 July 1998

NEUCOL DS100A

16/1/01

SW

