

A6.1 Test Results P856RA
Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

SpO₂ Testing
Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		98	94	80	70	60	98	94	80	70	60
9732-1	7H00001	98	94	80	70	60	0	0	0	0	0
9732-1	7H00002	98	94	80	70	60	0	0	0	0	0
9732-1	7H00003	98	94	80	70	60	0	0	0	0	0
9732-1	7H00004	98	94	80	70	60	0	0	0	0	0
9732-1	7H00005	98	94	80	70	60	0	0	0	0	0
9732-1	7H00006	98	94	80	70	60	0	0	0	0	0
9732-1	7H00007	98	94	80	70	60	0	0	0	0	0
9732-1	7H00008	98	94	80	69	60	0	0	0	-1	0
9732-1	7H00009	98	94	80	70	60	0	0	0	0	0
9732-1	7H00010	98	94	80	70	60	0	0	0	0	0
9732-1	7H00011	98	94	80	70	60	0	0	0	0	0
9732-1	7H00012	98	94	80	70	60	0	0	0	0	0
9732-1	7H00013	98	94	80	70	60	0	0	0	0	0
9732-1	7H00014	98	94	80	70	60	0	0	0	0	0
9732-1	7H00015	98	94	80	70	60	0	0	0	0	0
9732-1	7H00016	98	94	80	70	60	0	0	0	0	0
9732-1	7H00017	98	94	80	70	60	0	0	0	0	0
9732-1	7H00018	98	94	80	70	60	0	0	0	0	0
9732-1	7H00019	98	94	80	70	60	0	0	0	0	0
9732-1	7H00020	98	94	80	70	60	0	0	0	0	0
9732-1	7H00021	98	94	80	70	60	0	0	0	0	0
9732-1	7H00022	98	94	80	70	60	0	0	0	0	0
9732-1	7H00023	98	94	80	70	60	0	0	0	0	0
9732-1	7H00024	98	94	80	70	60	0	0	0	0	0
9732-1	7H00025	98	94	80	70	60	0	0	0	0	0
9732-1	7H00026	98	94	80	70	60	0	0	0	0	0
9732-1	7H00027	98	94	80	70	60	0	0	0	0	0
9732-1	7H00028	98	94	80	70	60	0	0	0	0	0
9732-1	7H00029	98	94	80	70	60	0	0	0	0	0
9732-1	7H00030	98	94	80	70	60	0	0	0	0	0
9732-1	7H00031	98	94	80	70	60	0	0	0	0	0
9732-1	7H00032	98	94	80	70	60	0	0	0	0	0
9732-1	7H00033	98	94	80	70	60	0	0	0	0	0
9732-1	7H00034	98	94	80	70	60	0	0	0	0	0
9732-1	7H00035	98	94	80	70	60	0	0	0	0	0
9732-1	7H00036	98	94	80	70	60	0	0	0	0	0
9732-1	7H00037	98	94	80	70	60	0	0	0	0	0
9732-1	7H00038	98	94	80	70	60	0	0	0	0	0
9732-1	7H00039	98	94	80	70	60	0	0	0	0	0
9732-1	7H00040	98	94	80	70	60	0	0	0	0	0
9732-1	7H00041	98	94	80	70	60	0	0	0	0	0
9732-1	7H00042	98	94	80	70	60	0	0	0	0	0
9732-1	7H00043	98	94	80	70	60	0	0	0	0	0

Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

9732-1	7H00044	98	94	80	70	60	0	0	0	0	0
9732-1	7H00045	98	94	80	70	60	0	0	0	0	0
9732-1	7H00046	98	94	80	70	60	0	0	0	0	0
9732-1	7H00047	98	94	80	70	60	0	0	0	0	0
9732-1	7H00048	98	94	80	70	60	0	0	0	0	0
9732-1	7H00049	98	94	80	70	60	0	0	0	0	0
9732-1	7H00050	98	94	80	70	60	0	0	0	0	0
9732-2	7H00051	98	94	80	70	60	0	0	0	0	0
9732-2	7H00052	98	94	80	70	60	0	0	0	0	0
9732-2	7H00053	98	94	80	70	60	0	0	0	0	0
9732-2	7H00054	98	94	80	70	60	0	0	0	0	0
9732-2	7H00055	98	94	80	70	60	0	0	0	0	0
9732-2	7H00056	98	94	80	70	60	0	0	0	0	0
9732-2	7H00057	98	94	80	70	60	0	0	0	0	0
9732-2	7H00058	98	94	80	70	60	0	0	0	0	0
9732-2	7H00059	98	94	80	70	60	0	0	0	0	0
9732-2	7H00060	98	94	80	70	60	0	0	0	0	0
9732-2	7H00061	98	94	80	70	60	0	0	0	0	0
9732-2	7H00062	98	94	80	70	60	0	0	0	0	0
9732-2	7H00063	98	94	80	70	60	0	0	0	0	0
9732-2	7H00064	98	94	80	70	60	0	0	0	0	0
9732-2	7H00065	98	94	80	70	60	0	0	0	0	0
9732-2	7H00066	98	94	80	70	60	0	0	0	0	0
9732-2	7H00067	98	94	80	70	60	0	0	0	0	0
9732-2	7H00068	98	94	80	70	60	0	0	0	0	0
9732-2	7H00069	98	94	80	70	60	0	0	0	0	0
9732-2	7H00070	98	94	80	70	60	0	0	0	0	0
9732-2	7H00071	98	94	80	70	60	0	0	0	0	0
9732-2	7H00072	98	94	80	70	60	0	0	0	0	0
9732-2	7H00073	98	94	80	70	60	0	0	0	0	0
9732-2	7H00074	98	94	80	70	60	0	0	0	0	0
9732-2	7H00075	98	94	80	70	60	0	0	0	0	0
9732-2	7H00076	98	94	80	70	60	0	0	0	0	0
9732-2	7H00077	98	94	80	70	60	0	0	0	0	0
9732-2	7H00078	98	94	80	70	60	0	0	0	0	0
9732-2	7H00079	98	94	80	70	60	0	0	0	0	0
9732-2	7H00080	98	94	80	70	60	0	0	0	0	0
9732-2	7H00081	98	94	80	70	60	0	0	0	0	0
9732-2	7H00082	98	94	80	70	60	0	0	0	0	0
9732-2	7H00083	98	94	80	70	60	0	0	0	0	0
9732-2	7H00084	98	94	80	70	60	0	0	0	0	0
9732-2	7H00085	98	94	80	70	60	0	0	0	0	0
9732-2	7H00086	98	94	80	70	60	0	0	0	0	0

Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

9732-2	7H00087	98	94	80	70	60	0	0	0	0	0
9732-2	7H00088	98	94	80	70	60	0	0	0	0	0
9732-2	7H00089	98	94	80	70	60	0	0	0	0	0
9732-2	7H00090	98	94	80	70	60	0	0	0	0	0
9732-2	7H00091	98	94	80	70	60	0	0	0	0	0
9732-2	7H00092	98	94	80	70	60	0	0	0	0	0
9732-2	7H00093	98	94	80	70	60	0	0	0	0	0
9732-2	7H00094	98	94	80	70	60	0	0	0	0	0
9732-2	7H00095	98	94	80	70	60	0	0	0	0	0
9732-2	7H00096	98	94	80	70	60	0	0	0	0	0
9732-2	7H00097	98	94	80	70	60	0	0	0	0	0
9732-2	7H00098	98	94	80	70	60	0	0	0	0	0
9732-2	7H00099	98	94	80	70	60	0	0	0	0	0

Tester: [Signature]

Date: 9/23/97
10/10/97

Q.A.: [Signature]

Date: 9/23/97

A6.1 Test Results P856RA

Pulse Rate Testing

Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

Lot Number	Serial Number	Test Results					Bias Error				
		250	150	95	80	50	250	150	95	80	50
9732-1	7H00001	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00002	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00003	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00004	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00005	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00006	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00007	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00008	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00009	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00010	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00011	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00012	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00013	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00014	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00015	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00016	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00017	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00018	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00019	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00020	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00021	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00022	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00023	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00024	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00025	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00026	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00027	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00028	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00029	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00030	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00031	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00032	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00033	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00034	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00035	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00036	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00037	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00038	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00039	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00040	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00041	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00042	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00043	249	150	95	80	50	-1	0	0	0	0

Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

9732-1	7H00044	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00045	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00046	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00047	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00048	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00049	249	150	95	80	50	-1	0	0	0	0
9732-1	7H00050	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00051	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00052	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00053	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00054	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00055	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00056	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00057	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00058	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00059	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00060	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00061	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00062	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00063	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00064	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00065	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00066	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00067	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00068	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00069	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00070	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00071	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00072	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00073	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00074	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00075	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00076	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00077	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00078	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00079	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00080	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00081	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00082	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00083	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00084	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00085	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00086	249	149	95	80	50	-1	-1	0	0	0

Pulse Oximeter Monitor: Nellcor N200
Monitor S/N: 20619220

Finger Probe Model: P856RA
BioTek Simulator S/N: 125672

9732-2	7H00087	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00088	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00089	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00090	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00091	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00092	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00093	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00094	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00095	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00096	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00097	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00098	249	150	95	80	50	-1	0	0	0	0
9732-2	7H00099	249	150	95	80	50	-1	0	0	0	0

Tester: [Signature]

Date: 9/23/97

Q.A.: [Signature]

Date: 9/23/97