

Training Feedback Form

Training Course Completed:					
Date:	22-5-23	Time/Length:	30mins	Trainer:	S. Harinder
Content	Yes	No	Unsure		
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Relevance	Yes	No	Unsure		
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Trainer	Yes	No	Unsure		
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Comments					
Do you require any further training in this area?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
If so, what would you like this training to cover?					
Further comments:					
Name:	JANINE GILL				
Date:	22-5-23				

*Please delete as applicable

Competency Assessment Questions – Envitec Oxygen Monitoring

Please refer to any system resources that you have access to in order to locate the information. The training materials will be made available on the system.

If you are unable to find the information, please make notes at the end of this document detailing where you struggled.

1) What is the difference between an oxygen analyser and an oxygen monitor?

Oxygen analyser - Not for continuous monitoring
Oxygen monitor - For Continuous monitoring

2) What is the stated accuracy of both devices?

Less than 1%

3) What is the purpose of the Perfusion Adapter Kit?

This kit helps connect chambers to a wide variety of ID tubing etc.

4) What is the part number of a replacement sensor for the OxiQuant S?

0110660

5) What is the stated response time for a MySign O for a 90% change in value?

13 Seconds

6) What is the part number of the MySign O kit, complete with power supply?

011277

Notes or Comments:

Name:

JANINE GILL

Date:

22-5-23

