Training Feedback Form

Training Course Completed: TOF 3D Technical Training					
Date: 27/2/23	Time/Length: 45 mins	Trainer:	Steve Hardaker		
Content		Yes	No	Unsure	
Was the course content presented in a logical					
manner?					
Was the course content and material complete and					
comprehensive?					
Will this information be useful to you in your job		/			
role?					
Relevance		Yes	No	Unsure	
Do you feel you now have a better understanding of		/			
the product/procedure/training area*?					
Did the course challenge y	our thinking and				
understanding of the product/procedure/training					
area*?					
Do you feel the training is beneficial to your team?					
Trainer		Yes	No	Unsure	
Did the trainer communicate and explain the		103	NO	Offsure	
material clearly?					
Did you feel the instructor was knowledgeable in the			-		
area covered?					
Did the trainer enopurage discussions and					
questions?					
Comments					
Do you require any further training in this area?					
If so, what would you like the			-		
•					
•					
Further comments:					
•					
	1				
Name: M. GLOCK	7				
Date: 1-3-2	5		*Dia-		
			"Please delet	e as applicable	

TOF 3D Technical Training - Competency Assessment Questions

Please refer to any system resources that you have access to in order to locate the information. The training materials will be made available on the system and will appear as required reading. If you are unable to find the information, please make notes at the end of this document detailing where you struggled.

1) What type of monitoring method does the TOF 3D use: subjective or quantitative?

1 It was quantitative but can use subjective.
2) How does the TOF 3D measure twitch response?
Accelero myography
3) Where is the most common place on the patient to measure twitch response? The thunk
4) What is the warranty period? $2 + month 5$
5) Does the TOF 3D or accessories contain any Latex?
No
6) Where can you find cleaning guidelines? Stock pase - technical documents 7) What is the name for the stimulation mode that delivers 4 short stimuli over a 2s period? TOF - train of four
8) Which accessories are single-patient use: hand adapter / thumb adapter / eye adapter? thumb adapter + eye adapter
Notes or Comments:
•
Name: