

Training Feedback Form

Training Course Completed: 400 A C Temperature Probe			
Date:	Time/Length:	Trainer:	
Content	Yes	No	Unsure
Was the course content presented in a logical manner?	/		
Was the course content and material complete and comprehensive?	/		
Will this information be useful to you in your job role?			/
Relevance	Yes	No	Unsure
Do you feel you now have a better understanding of the product/procedure/training area*?	/		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	/		
Do you feel the training is beneficial to your team?	/		
Trainer	Yes	No	Unsure
Did the trainer communicate and explain the material clearly?	/		
Did you feel the instructor was knowledgeable in the area covered?	/		
Did the trainer encourage discussions and questions?	/		
Comments			
Do you require any further training in this area?		/	
If so, what would you like this training to cover?			
Further comments:			
Name: Michael Lamb			
Date:			

*Please delete as applicable

Competency Assessment Questions – 400-AC Temperature Probes

Please refer to any system resources that you have access to in order to locate the information. The training materials will be made available on the system and will appear as required reading. If you are unable to find the information, please make notes at the end of this document detailing where you struggled.

1) List the 3 patient sizes of 400-AC probes available.

neonatal
paediatric
adult

2) What are the 2 main types of probes available? (in terms of placement / application)

Skin
Oropharyngeal/Rectal

3) Which placement of temperature monitoring is considered to be the most accurate indication of core body temperature?

Rectal

4) What is the warranty period?

12 months

5) Do 400-AC probes contain any Latex?

no

6) Where can you find cleaning guidelines?

attached to the product pack.
in instructions for use.

7) What is another name for a steam sterilizer?

autoclave

Notes or Comments:

Name:

Michael Lamb

Date:

