

### Competency Assessment Questions – Phototherapy Masks

Please refer to any system resources that you have access to in order to locate the information. The training materials will be made available on the system.

If you are unable to find the information, please make notes at the end of this document detailing where you struggled.

1) What is the colour of light typically emitted by phototherapy lights?

Blue

2) List the 3 sizes in the EyeMax2 range.

- Regular
- Premie
- Micro

3) How many phototherapy masks are supplied in a pack?

20

4) Can phototherapy masks be re-used on more than one patient?

No

5) Where on the patient do you measure the occipital front circumference?

over the most prominent part on the back of the head (occiput) and just above the eyebrows (supraorbital ridges).

6) Which part number NeoMask should be used on a patient with an occipital front circumference of 30 cm?

1114016 - Medium

Notes or Comments:

Name: Sophie Lines Date: 27/3/2023