

Training Feedback Form

Training Course Completed: Temperature Probes (Technical).			
Date: 6/2/23	Time/Length: 20 mins	Trainer: Steve H	
Content	Yes	No	Unsure
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Relevance	Yes	No	Unsure
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Trainer	Yes	No	Unsure
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Comments			
Do you require any further training in this area?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If so, what would you like this training to cover? N/A			
Further comments: N/A			
Name: Sophie Lines			
Date: 6-2-2023			

Competency Assessment Questions – 400-AC Temperature Probes

Please refer to any system resources that you have access to in order to locate the information. The training materials will be made available on the system and will appear as required reading. If you are unable to find the information, please make notes at the end of this document detailing where you struggled.

- 1) List the 3 patient sizes of 400-AC probes available. Adult
Paediatric
Neonatal
- 2) What are the 2 main types of probes available? (in terms of placement / application)
Skin/surface
oesophageal/rectal
- 3) Which placement of temperature monitoring is considered to be the most accurate indication of core body temperature?
oesophageal/rectal
- 4) What is the warranty period? 12 months from date of invoice
- 5) Do 400-AC probes contain any Latex? NO
- 6) Where can you find cleaning guidelines? In the Instructions for use.
- 7) What is another name for a steam sterilizer? Auto-clave

Notes or Comments:

Name: Sophie Lines Date: 6/2/23

