





































‘Y’ Probe Chart

856 Nellcor	D-type 9-pin			
857 Spaelabs	Hypertronic s 12-pin			
858 Nellcor	D-type 9-pin			
859 Pacotech	Switchcraft 5-pin			
860 Simed Baxter	Hypertronic s 9-pin			
861 BCI	D-type 9-pin			
864 Datascope	Din 8-pin			
866 Simed Baxter	Hypertronic s 9-pin			
867 Ohmeda	Hypertronic s 9-pin			
869 Criticare	D-type 9-pin			
871 Nonin	D-type 9-pin			
873 Datex	D-type 9-pin			
875 Novametri x	Hypertronic s 7-pin			
876 Novametri x	Hypertronic s 7-pin			
878 Invivo	Lemo 7-Pin			
886 Kontron	DPK1 12-Pin			

‘Y’ Probe Chart

				
887 Spaelabs	Hypertronic s 12-pin			
888 S&W	SPK1			

SpO₂ 'Y' Probe Details

Model	Compatible manufacturer	Length M (ft)	Connector type	Model	Notes Document Name VM3/COP/.....	Colours			
						Probes	Button	CE	Db
856YA	Nellcor	0.9 (3)	D-M-9	856YA	VM3/COP/35.03	White	Blue	White	Blue
857YA	Spacelabs	3 (10)	Hypertronics M-12	857YA	VM3/COP/35.04	White	White	White	Orange
858YA	Nellcor	3.65 (12)	D-M-9	858YA	VM3/COP/35.05	White	Blue	White	Blue
859YA	Pace Tech	3.65 (12)	Switchcraft F-5	859YA	VM3/COP/35.06	White	White	White	Orange
860YA	Simed Baxter	3.65 (12)	Hypertronics M-9	860YA	VM3/COP/35.07	White	White	White	Orange
861YA	BCI	0.9 (3)	D-M-9	861YA	VM3/COP/35.08	White	Orange	White	Or/Br
862YA	Criticare	3 (10)	Lemo M-5	862YA	Do Not Make	-	-	-	-
863YA	Datascope	0.9(3)	D-M-9	863YA	Do Not Make	-	-	-	-
864YA	Datascope	3.65 (12)	DIN M-8	864YA	VM3/COP/35.11	Black	Red	White	White
865YA	Sensormedics	0.9(3)	Compushield	865YA	Do Not Make	-	-	-	-
866YA	Simed Baxter	3.65 (12)	Hypertronics M-9	866YA	VM3/COP/35.13	White	White	White	Orange
867YA	Ohmeda	3.65 (12)	Hypertronics M-9	867YA	VM3/COP/35.14	Black	Green	White	Green
869YA	Criticare	0.9 (3)	D-M-9	869YA	VM3/COP/35.15	White	Purple	White	Purple
871YA	Nonin	0.9 (3)	D-M-9	871YA	VM3/COP/35.16	White	White	White	Orange
872YA	Datex	3.65(12)	DPK1	872YA	Do Not Make	-	-	-	-
873YA	Datex	0.9 (3)	D-M-9	873YA	VM3/COP/35.18	White	Tan	White	Brown
875YA	Novametrix	3 (10)	Hypertronics M-7	875YA	VM3/COP/35.20	Black	Yellow	White	Yellow
876YA	Novametrix	3 (10)	Hypertronics M-7	876YA	VM3/COP/35.21	Black	Yellow	White	Yellow
877YA	Critikon	0.9(3)	D-M-9	877YA	Do Not Make	-	-	-	-
878YA	Invivo	3 (10)	Lemo M-7	878YA	VM3/COP/35.23	White	White	White	Orange
886YA	Kontron	3.65 (12)	SPK1-A M-12	886YA	VM3/COP/35.26	White	White	White	Orange
887YA	Spacelabs	3 (10)	Hypertronics M-12	887YA	VM3/COP/35.27	Black	White	White	Orange
888YA	S&W	3.65 (12)	SPK1-O F-10	888YA	VM3/COP/35.28	White	White	White	Orange
889YA	Criticare	3(10)	D-M-9	889YA	Do Not Make	-	-	-	-
890YA	Critikon & Sensormedics	3(10)	Thorkom F-12	890YA	Do Not Make	-	-	-	-
891YA	Critikon	3(10)	Lemo M-7	891YA	Do Not Make	-	-	-	-
Next Available Number					VM/COP/35.29				

Drawing No	Brief Description	Part No	Drawn	Date	Rev	No	Date	Change
SPY-856	Nellcor - D type M9	0018562	K.J.R.	04.11.03				
SPY-857	Spacelabs - Hype M12	0018572	K.J.R.	04.11.03				
SPY-858	Nellcor - D type M9	0018582	K.J.R.	04.11.03				
SPY-859	Pacetek - Switchcraft F5	0018592	K.J.R.	04.11.03				
SPY-860	Simed-Baxter - Hype M9	0018602	K.J.R.	04.11.03				
SPY-861	BCI - D type M9	0018612	K.J.R.	04.11.03				
SPY-862	Criticare - Lemo M9	0018622	K.J.R.	04.11.03				
SPY-863	Datascope - D type M9	0018632	K.J.R.	04.11.03				
SPY-864	Datascope - DIN M8	0018642	K.J.R.	04.11.03				
SPY-865	S/Medics - Compushield M6	0018652	K.J.R.	04.11.03				
SPY-866	Simed-Baxter - Hype M9	0018662	K.J.R.	04.11.03				
SPY-867	Ohmeda - Hype M9	0018672	K.J.R.	04.11.03				
SPY-869	Criticare - D type M9	0018692	K.J.R.	04.11.03				
SPY-871	Nonin - D type M9	0018712	K.J.R.	04.11.03				
SPY-872	Datex - DPK1 - F10	0018722	K.J.R.	04.11.03				
SPY-873	Datex - D type M9	0018732	K.J.R.	04.11.03				
SPY-875	Novamatrix - Hype M7	0018752	K.J.R.	04.11.03				
SPY-876	Novamatrix - Hype M7	0018762	K.J.R.	04.11.03				
SPY-877	Critikon - D type M9	0018772	K.J.R.	04.11.03				
SPY-878	Invivo - Lemo M7	0018782	K.J.R.	04.11.03				
SPY-886	Kontron - SPK1 A - M12	0018582	K.J.R.	04.11.03				
SPY-887	Spacelabs - Hype M12	0018872	K.J.R.	04.11.03				
SPY-888	S & W - SPK1 O - F12	0018882	K.J.R.	04.11.03				
SPY-889	Criticare - D type M9	0018892	K.J.R.	04.11.03				
SPY-890	Critikon & S/med - Thorkon F12	0018902	K.J.R.	04.11.03				
SPY-891	Critikon - Lemo M7	0018912	K.J.R.	04.11.03				

SpO₂ 'Y' Probes Index
Drawings Confirmation

Model	Compatible manufacturer	Length M (ft)	Drawing No.	Notes	Date	Model
856YA	Nellcor	0.9 (3)	SPY-856			856YA
857YA	Spacelabs	3 (10)	SPY-857			857YA
858YA	Nellcor	3.65 (12)	SPY-858			858YA
859YA	Pace Tech	3.65 (12)	SPY-859			859YA
860YA	Simed Baxter	3.65 (12)	SPY-860			860YA
861YA	BCI	0.9 (3)	SPY-861			861YA
862YA	Criticare			Don't Make		862YA
863YA	Datascope			Don't Make		863YA
864YA	Datascope	3.65 (12)	SPY-864			864YA
865YA	Sensormedics			Don't Make		865YA
866YA	Simed Baxter	3.65 (12)	SPY-866			866YA
867YA	Ohmeda	3.65 (12)	SPY-867	Yet To Confirm		867YA
869YA	Criticare	0.9 (3)	SPY-869			869YA
871YA	Nonin	0.9 (3)	SPY-871			871YA
872YA	Datex			Don't Make		872YA
873YA	Datex	0.9 (3)	SPY-873			873YA
875YA	Novamatrix	3 (10)	SPY-875			875YA
876YA	Novamatrix	3 (10)	SPY-876			876YA
877YA	Critikon			Don't Make		877YA
878YA	Invivo	3 (10)	SPY-878			878YA
886YA	Kontron	3.65 (12)	SPY-886			886YA
887YA	Spacelabs	3 (10)	SPY-887			887YA
888YA	S&W	3.65 (12)	SPY-888			888YA
889YA	Criticare			Don't Make		889YA
890YA	Critikon & Sensormedics			Don't Make		890YA
891YA	Critikon			Don't Make		891YA

All the above documents have been checked and can be confirmed as being correct as of the date specified.

Approved By:

Name:

Job Title:

Dept:

Signature:

Date:/...../.....